碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301 或 800-8283301
订货 e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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	Apoptosis	Necroptosis	Pyroptosis
Trigger	Death receptor activation	Death receptor activation	Activation of inflammasomes

		when caspase-8 is inhibited	by pathogens
Inflammatory response	No	Yes	Yes
Key proteins	Caspase family	RIPK1、RIPK3、MLKL	Caspases 1、4、5 and 11、Gasdermin D
Morphology	Cell rounding、blebbing、formation of apoptotic bodies	Organ swelling、membrane rupture、release of cell contents	Plasma membrane rupture、release of cell contents、maintained mitochondrial integrity
Programmed cell death	Yes	Yes	Yes

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- Chen X,Hao K,Yu X,Huang A,Zhu B,Wang GX,Ling F.Magnolol protects Ctenopharyngodon idella kidney cells from apoptosis induced by grass carp reovirus.FISH SHELLFISH IMMUN . 2018 Mar;74:426-435
- Lei Q,Tan J,Yi S,Wu N,Wang Y,Wu H.Mitochondric acid 5 activates the MAPK-ERK-yap signaling pathways to protect mouse microglial BV-2 cells against TNF α -induced apoptosis via increased Bnip3-related mitophagy.Cell Mol Biol Lett . 2018 Apr 5;23:14
- Wang H,Zhao X,Ni C,Dai Y,Guo Y.Zearalenone regulates endometrial stromal cell apoptosis and migration via the promotion of mitochondrial fission by activation of the JNK/Drp1 pathway.Mol Med Rep . 2018 Jun;17(6):7797-7806
- Wang Y,Xu H,Lu Z,Yu X,Lv C,Tian Y,Sui D.Pseudo-Ginsenoside Rh2 induces A549 cells apoptosis via the Ras/Raf/ERK/p53 pathway.Exp Ther Med . 2018 Jun;15(6):4916-4924
- Zhu G,Zhang W,Liu Y,Wang S.miR-371b-5p inhibits endothelial cell apoptosis in monocrotaline-induced pulmonary arterial hypertension via PTEN/PI3K/Akt signaling pathways.Mol Med Rep . 2018 Dec;18(6):5489-5501
- Chen C,Wang S,Liu P.Deferoxamine Enhanced Mitochondrial Iron Accumulation and Promoted Cell Migration in Triple-Negative MDA-MB-231 Breast Cancer Cells Via a ROS-Dependent Mechanism.Int J Mol Sci. 2019 Oct 8;20(19). pii: E4952

Version 2021.09.01