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- Zhang C, Wu R, Zhu H, Hu YZ, Jiang H, Lin NM, He QJ, Yang B. Enhanced anti-tumor activity by the combination of TRAIL/Apo-2L and combretastatin A-4 against human colon cancer cells via induction of apoptosis in vitro and in vivo. *Cancer Lett.* 2011 Mar 1;302(1):11-9.
- Xu X, Gao X, Jin L, Bhadury PS, Yuan K, Hu D, Song B, Yang S. Antiproliferation and cell apoptosis inducing bioactivities of constituents from *Dyosma versipellis* PC3 and Bcap-37 cell lines. *Cell Div.* 2011 Jun 15;6(1):14.
- Liu M, Yang S, Jin L, Hu D, Wu Z, Yang S. Chemical Constituents of the Ethyl Acetate Extract of *Belamcanda chinensis* (L.) DC Roots and Their Antitumor Activities. *Molecules.* 2012 May 24;17(5):6156-69.
- Peng YR, Ding YF, Wei YJ, Shu B, Li YB, Liu XD. Caudatin-2,6-dideoxy-3-O-methyl-β-D-cymaropyranoside 1 induced apoptosis through caspase 3-dependent pathway in human hepatoma cell line SMMC7721. *Phytother Res.* 2011 May;25(5):631-7.
- Liu MC, Yang SJ, Jin LH, Hu DY, Xue W, Song BA, Yang S. Synthesis and cytotoxicity of novel ursolic acid derivatives containing an acyl piperazine moiety. *Eur J Med Chem.* 2012 Dec;58:128-35.
- Liu L, Wang P, Liang C, He D, Yu Y, Liu X. Distinct effects of Nampt inhibition on mild and severe models of lipopolysaccharide-induced myocardial impairment. *Int Immunopharmacol.* 2013 Oct;17(2):342-9.
- Zhu Z, Li K, Xu D, Liu Y, Tang H, Xie Q, Xie L, Liu J, Wang H, Gong Y, Hu Z, Zheng J. Erratum to: ZFX regulates glioma cell proliferation and survival in vitro and in vivo. *J Neurooncol.* 2013 Jul;113(3):533.
- Qi R, Yang F, Huang J, Peng H, Liu Y, Liu Z. Supplementation with conjugated linoleic acid decreases pig back fat deposition by inducing adipocyte apoptosis. *BMC Vet Res.* 2014 Jun 26;10:141.
- Wang N, Wei H, Yin D, Lu Y, Zhang Y, Jiang D, Jiang Y, Zhang S. Cyclin D1b overexpression inhibits cell proliferation and induces cell apoptosis in cervical cancer cells in vitro and in vivo. *Int J Clin Exp Pathol.* 2014 Jun 15;7(7):4016-23.
- You H, Li T, Zhang J, Lei Q, Tao X, Xie P, Lu W. Reduction in ischemic cerebral infarction is mediated through golgi phosphoprotein 3 and Akt/mTOR signaling following salvianolate administration. *Curr Neurovasc Res.* 2014 May;11(2):107-13.
- Wu Z, Kong X, Zhang T, Ye J, Fang Z, Yang X. Pseudoephedrine/ephedrine shows potent anti-inflammatory activity against TNF-α-mediated acute liver failure induced by lipopolysaccharide/D-galactosamine. *Eur J Pharmacol.* 2014 Feb 5;724:112-21.
- Hao R, Hu X, Wu C, Li N. Hypoxia-induced miR-15a promotes mesenchymal ablation and adaptation to hypoxia during lung development in chicken. *PLoS One.* 2014 Jun 2;9(6):e98868.
- Shi Y, Song Q, Yu S, Hu D, Zhuang X. Microvascular invasion in hepatocellular carcinoma overexpression promotes cell proliferation and inhibits cell apoptosis of hepatocellular carcinoma via inhibiting miR-199a expression. *Oncotargets Ther.* 2015 Aug 27;8:2303-10.
- Ke J, Yao YL, Zheng J, Wang P, Liu YH, Ma J, Li Z, Liu XB, Li ZQ, Wang ZH, Xue YX. Knockdown of long non-coding RNA HOTAIR inhibits malignant biological behaviors of human glioma cells via modulation of miR-326. *Oncotarget.* 2015 Sep 8;6(26):21934-49.

15. Lu C, Jiang Q, Hu M, Tan C, Yu H, Hua Z. Preliminary biological evaluation of ¹⁸F-FBEM-Cys-Annex V a novel apoptosis imaging agent. *Molecules*. 2015 Mar 17;20(3):4902-14.
16. Xu P, Meng Q, Sun H, Yin Q, Yu H, Zhang Z, Cao M, Zhang Y, Li Y. Shrapnel nanoparticles loading docetaxel inhibit metastasis and growth of breast cancer. *Biomaterials*. 2015 Sep;64:10-20.
17. Li Z, Yin H, Hao S, Wang L, Gao J, Tan X, Yang Z. miR-200 family promotes podocyte differentiation through repression of RSAD2. *Sci Rep*. 2016 Jun 2;6:27105.
18. Han R, Wang R, Zhao Q, Han Y, Zong S, Miao S, Song W, Wang L. Trim69 regulates zebrafish brain development by ap-1 pathway. *Sci Rep*. 2016 Apr 6;6:24034.
19. Meng C, Song L, Wang J, Li D, Liu Y, Cui X. Propofol induces proliferation partially via down regulation of p53 protein and promotes migration via activation of the Nrf2 pathway in human breast cancer cell line MDA-MB-231. *Oncol Rep*. 2017 Feb;37(2):841-848.
20. Liu W, Gao Y, Chang N. Nurr1 overexpression exerts neuro protective and anti-inflammatory roles via down-regulating CCL2 expression in both in vivo and in vitro Parkinson's disease models. *Biochem Biophys Res Commun*. 2017 Jan 22;482(4):1312-9.
21. Hu SY, Zhang Y, Zhu PJ, Zhou H, Chen YD. Liraglutide directly protects cardiomyocytes against reperfusion injury possibly via modulation of intracellular calcium homeostasis. *J Geriatr Cardiol*. 2017 Jan;14(1):57-66.
22. Meng C, Song L, Wang J, Li D, Liu Y, Cui X. Propofol induces proliferation partially via downregulation of p53 protein and promotes migration via activation of the Nrf2 pathway in human breast cancer cell line MDA-MB-231. *Oncol Rep*. 2017 Feb;37(2):841-848.
23. Yan Y, Wang L, He J, Liu P, Lv X, Zhang Y, Xu X, Zhang L, Zhang Y. Synergy with interferon-lambda 3 and sorafenib suppresses hepatocellular carcinoma proliferation. *Biomed Pharmacother*. 2017 Apr;88:395-402.
24. Geng YD, Zhang C, Lei JL, Yu P, Xia YZ, Zhang H, Yang L, Kong LY. Walsurionoid B induces mitochondrial and lysosomal dysfunction leading to apoptotic rather than autophagic cell death via ROS/p53 signaling pathways in liver cancer. *Biochem Pharmacol*. 2017 Oct 15;142:71-86.
25. Lu C, Jiang Q, Hu M, Tan C, Yu H, Hua Z. Preliminary biological evaluation of ¹⁸F-AIF-NOTA-MAL-Cys-Annexin V as a novel apoptosis imaging agent. *Oncotarget*. 2017 Apr 10;8(31):51086-51095.
26. Liu W, Gao Y, Chang N. Nurr1 overexpression exerts neuroprotective and anti-inflammatory roles via down-regulating CCL2 expression in both in vivo and in vitro Parkinson's disease models. *Biochem Biophys Res Commun*. 2017 Jan 22;482(4):1312-1319.
27. Wang J, Wang WN, Xu SB, Wu H, Dai B, Jian DD, Yang M, Wu YT, Feng Q, Zhu JH, Zhang L, Zhang L. MicroRNA-214-3p: A link between autophagy and endothelial cell dysfunction in atherosclerosis. *Acta Physiol (Oxf)*. 2017 Sep 9.
28. Zhao K, Li G, Yao Y, Zhou Y, Li Z, Guo Q, Lu N. Activation of phospholipase C- γ 1 and translocation of phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase contribute to GL-V9-induced apoptosis in human gastric cancer cells. *Exp Cell Res*. 2017 Jul 1;356(1):8-19.
29. Pan T, Zhang M, Zhang F, Yan G, Ru Y, Wang Q, Zhang Y, Wei X, Xu X, Shen L, Zhang J, Wu K, Yao L, Li X. NDRG2 overexpression suppresses hepatoma cells survival during metabolic stress through disturbing the activation of fatty acid oxidation. *Biochem Biophys Res Commun*. 2017 Feb 5;483(2):860-866.
30. Song J, Li Y, Song J, Hou F, Liu B, Li A. Mangiferin protects mitochondrial function by preserving mitochondrial hexokinase-II in vessel endothelial cells. *Biochim Biophys Acta*. 2017 Jul;1863(7):1829-1839.
31. Li R, Ma J, Wu Y, Nangle M, Zou S, Li Y, Yin J, Zhao Y, Xu H, Zhang H, Li X, Ye QS, Wang J, Xiao J. Dual Delivery of NGF and bFGF Coacervate Ameliorates Diabetic Peripheral Neuropathy via Inhibiting Schwann Cells Apoptosis. *Int J Biol Sci*. 2017 May 16;13(5):640-651.
32. Liu ZY, Hu S, Zhong QW, Tian CN, Ma HM, Yu JJ. N-Methyl-D-Aspartate Receptor-Driven Calcium Influx Potentiates the Adverse Effects of Myocardial Ischemia-Reperfusion Injury Ex Vivo. *J Cardiovasc Pharmacol*. 2017 Nov;70(5):329-338.
33. Dai H, Zhang X, Yang Z, Li J, Zheng J. Effects of Baicalin on Blood Pressure and Left Ventricular Remodeling in Rats with Renovascular Hypertension. *Med Sci Monit*. 2017 Jun 16;23:2939-2948.
34. Chen D, Tao X, Zhou L, Sun F, Sun M, Fang X. Spaghetti, a homolog of human RPAP3 (RNA polymerase II-associated protein 3), determines the fate of germline stem cells in *Drosophila* ovary. *Cell Biol Int*. 2017 Nov 7.
35. Zhang X, Wang X, Liu T, Mo M, Ao L, Liu J, Cao J, Cui Z. ZnSO₄ rescued vimentin from collapse in DBP-exposed Sertoli cells by attenuating ER stress and apoptosis. *Toxicol In Vitro*. 2018 Apr;48:195-204.
36. Zheng X, Bi W, Yang G, Zhao J, Wang J, Li X, Zhou X. Hyperglycemia Induced by Chronic Restraint Stress in Mice Is Associated With Nucleus Tractus Solitarius Injury and Not Just the Direct Effect of Glucocorticoids. *Front Neurosci*. 2018 Dec 19;12:983.
37. Huang Y, Wu S, Zhang Y, Wang L, Guo Y. Antitumor effect of triptolide in T-cell lymphoblastic lymphoma by inhibiting cell viability, invasion, and epithelial-mesenchymal transition via regulating the PI3K/AKT/mTOR pathway. *Oncol Targets Ther*. 2018 Feb 12;11:769-779.
38. Jiang C, Yang Y, Yang Y, Guo L, Huang J, Liu X, Wu C, Zou J. Long Noncoding RNA (lncRNA) HOTAIR Affects Tumorigenesis and Metastasis of Non-Small Cell Lung Cancer by Upregulating miR-613. *Oncol Res*. 2018 Jun 11;26(5):725-734.
39. Feng Q, He P, Wang Y. MicroRNA-223-3p regulates cell chemo-sensitivity by targeting FOXO3 in prostatic cancer. *Gene*. 2018 Jun 5;658:152-158.
40. Luo J, Zhu H, Jiang H, Cui Y, Wang M, Ni X, Ma C. The effects of aberrant expression of lncRNA DGCR5/miR-873-5p/TUSC3 in lung cancer cell progression. *Cancer Med*. 2018 May 23.
41. Liu J, Zhou S, Zhang Y, Li X, Qian X, Tao W, Jin L, Zhao J. Bax inhibitor-1 suppresses early brain injury following experimental subarachnoid hemorrhage in rats. *Int J Mol Med*. 2018 Nov;42(5):2891-2902.
42. Liu Z, Li H, Hong C, Chen M, Yue T, Chen C, Wang Z, You Q, Li C, Weng Q, Xie H, Hu R. ALS-Associated E478G Mutation in Human OPTN (Optineurin) Promotes Inflammation and Induces Neuronal Cell Death. *Front Immunol*. 2018 Nov 14;9:2647.
43. Tang P, Dang H, Huang J, Xu T, Yuan P, Hu J, Sheng JF. NADPH oxidase NOX4 is a glycolytic regulator through mROS-HIF1 α axis in thyroid carcinomas. *Sci Rep*. 2018 Oct 26;8(1):15897.
44. Lou Z, Wang AP, Duan XM, Hu GH, Zuo ML, Yang ZB. Role of ALK5/SMAD2/3 signaling in the regulation of NOX expression in cerebral ischemia/reperfusion injury. *Exp Ther Med*. 2018 Sep;16(3):1671-1678.
45. Zhang H, Deng Z, Zhang D, Li H, Zhang L, Niu J, Zuo W, Fu R, Fan L, Ye JH, She J. High expression of leucine-rich repeat-containing 8A is indicative of a worse outcome of colon cancer patients by enhancing cancer cell growth and metastasis. *Oncol Rep*. 2018 Sep;40(3):1275-1286.
46. Zhang G, Wang Z, Chen W, Cao Y, Wu J, Qiang G, Ji A, Wu J, Jiang C. Dual effects of gossypol on human hepatocellular carcinoma via endoplasmic reticulum stress and autophagy. *Int J Biochem Cell Biol*. 2019 Aug;113:48-57.
47. Zheng Y, Zhang P, Zhang C, Zeng W. Surgery-induced cryptorchidism induces apoptosis and autophagy of spermatogenic cells in mice. *Zygote*. 2019 Apr;27(2):101-110.
48. Yuan Y, Wang YY, Liu X, Luo B, Zhang L, Zheng F, Li XY, Guo LY, Wang L, Jiang M, Pan YM, Yan YW, Yang JY, Chen SY, Wang JN, Tang JM. KPC1 alleviates hypoxia/reoxygenation-induced apoptosis in rat cardiomyocyte cells through BAX degradation. *J Cell Physiol*. 2019 Dec;234(12):22921-22934.
49. Xiao Q, Zeng FL, Tang GY, Lei CY, Zou XX, Liu XL, Peng BL, Qin S, Li HX. Expression of galectin-3 and apoptosis in placental villi from patients with missed abortion during early pregnancy. *Exp Ther Med*. 2019 Apr;17(4):2623-2631.
50. Li W, Qinghai S, Kai L, Xue M, Lili N, Jihua R, Zhengxiang L, Xiaoling L, Di G, Qi Y, Mengyun D, Jianfeng F. Oral administration of Ginkgolide B alleviates hypoxia-induced neuronal damage in rat hippocampus by inhibiting oxidative stress and apoptosis. *Iran J Basic Med Sci*. 2019 Feb;22(2):140-145.
51. Li J, Tan YJ, Wang MZ, Sun Y, Li GY, Wang QL, Yao JC, Yue J, Liu Z, Zhang GM, Ren YS. Loganetin protects against rhabdomyolysis-induced acute kidney injury by modulating the toll-like receptor 4 signalling pathway. *Br J Pharmacol*. 2019 Apr;176(8):1106-1121.
52. Ma W, Shen H, Li Q, Song H, Guo Y, Li F, Zhou X, Guo X, Shi J, Cui Q, Xing J, Deng J, Yu Y, Liu W, Zhao H. MARVELD1 attenuates arsenite-induced apoptosis in liver cancer cells by inhibiting reactive oxygen species production. *Ann Transl Med*. 2019 May;7(9):200.
53. Nie C, Qin X, Li X, Tian B, Zhao Y, Jin Y, Li Y, Wang Q, Zeng D, Hong

- A,Chen X.CACNA2D3 Enhances the Chemosensitivity of Esophageal Squamous Cell Carcinoma to Cisplatin via Inducing Ca²⁺-Mediated Apoptosis and Suppressing PI3K/Akt Pathways.Front Oncol. 2019 Apr 2;9:185.
54. Yang Y,Chen A,Ma J,Wu A,Xu F.Effects of radioactive 125I on apoptosis of HGC-27 gastric cancer cells.Oncol Lett. 2019 Nov;18(5):4916-4922.
 55. Zhang J,Fan J,Li Y,Liang S,Huo S,Wang X,Zuo Y,Cui D,Li W,Zhong Z,Zhong F.Porcine Parvovirus Infection Causes Pig Placenta Tissue Damage Involving Nonstructural Protein 1 (NS1)-Induced Intrinsic ROS/Mitochondria-Mediated Apoptosis.Viruses. 2019 Apr 25;11(4). pii: E389.
 56. Ngekre M X K,Jiang J,Enayatullah H,Ennab W,Mustafa S,Rodeni S,Wei Q,Shi F.Sweet taste receptor agonists alter ovarian functions and ovarian cycles in aged mice.Reprod Biol. 2019 Sep;19(3):230-236.
 57. Liu M,Ai W,Sun L,Fang F,Wang X,Chen S,Wang H.Triclosan-induced liver injury in zebrafish (Danio rerio) via regulating MAPK/p53 signaling pathway.Comp Biochem Physiol C Toxicol Pharmacol. 2019 Aug;222:108-117.
 58. Liang M,Jia J,Chen L,Wei B,Guan Q,Ding Z,Yu J,Pang R,He G.LncRNA MCM3AP-AS1 promotes proliferation and invasion through regulating miR-211-5p/SPARC axis in papillary thyroid cancer.Endocrine. 2019 Aug;65(2):318-326
 59. Liang Q,Li F,Li Y,Liu Y,Lan M,Wu S,Wu X,Ji Y,Zhang R,Yin L..Self-assisted membrane-penetrating helical polypeptides mediate anti-inflammatory RNAi against myocardial ischemic reperfusion (IR) injury.Biomater Sci. 2019 Aug 20;7(9):3717-3728.
 60. Wang H,Yu D,Fang J,Cao C,Liu Z,Ren J,Qu X.Renal-Clearable Porphyrinic Metal-Organic Framework Nanodots for Enhanced Photodynamic Therapy.ACS Nano. 2019 Aug 27;13(8):9206-9217.
 61. Wang A,Zhao C,Gao Y,Duan G,Yang Y,Fan B,Wang X,Wang K.LEF1-AS1 contributes to proliferation and invasion through regulating miR-544a/ FOXP1 axis in lung cancer.Invest New Drugs. 2019 Dec;37(6):1127-1134.
 62. Tang KK,Liu XY,Wang ZY,Qu KC,Fan RF.Trehalose alleviates cadmium-induced brain damage by ameliorating oxidative stress, autophagy inhibition, and apoptosis.Metallomics. 2019 Dec 11;11(12):2043-2051.
 63. Tang Y,Cao G,Zhao G,Wang C,Qin Q.LncRNA differentiation antagonizing non-protein coding RNA promotes proliferation and invasion through regulating miR-135a/NLRP37 axis in pancreatic cancer.Invest New Drugs. 2019 Jul 3
 64. Zheng Y,Liu W,Chen Y,Li C,Jiang H,Wang X.Conjugating gold nanoclusters and antimicrobial peptides: From aggregation-induced emission to antibacterial synergy.J Colloid Interface Sci. 2019 Jun 15;546:1-10.
 65. Zheng S,Li M,Miao K,Xu H.SNHG1 contributes to proliferation and invasion by regulating miR-382 in breast cancer.Cancer Manag Res. 2019 Jun 18;11:5589-5598.
 66. An W,Yu Y,Zhang Y,Zhang Z,Yu Y,Zhao X.Exogenous IL-19 attenuates acute ischaemic injury and improves survival in male mice with myocardial infarction.Br J Pharmacol. 2019 Mar;176(5):699-710.
 67. Huang S,Qi P,Zhang T,Li F,He X.The HIF-1 α /miR-224-3p/ATG5 axis affects cell mobility and chemosensitivity by regulating hypoxia-induced protective autophagy in glioblastoma and astrocytoma.Oncol Rep. 2019 Mar;41(3):1759-1768.
 68. Wu G,Zhu Q,Zeng J,Gu X,Miao Y,Xu W,Lv T,Song Y.Extracellular mitochondrial DNA promote NLRP3 inflammasome activation and induce acute lung injury through TLR9 and NF- κ B.J Thorac Dis. 2019 Nov;11(11):4816-4828.
 69. Yang W,Hao Y,Hou W,Fang X,Fang P,Jiang T,Ma C,Zhang QJieduan-Niwan Formula Reduces Liver Apoptosis in a Rat Model of Acute-on-Chronic Liver Failure by Regulating the E2F1-Mediated Intrinsic Apoptosis Pathway.Evid Based Complement Alternat Med. 2019 Nov 11;2019:8108503.
 70. Bai Y,Zhang G,Cheng R,Yang R,Chu H.CASC15 contributes to proliferation and invasion through regulating miR-766-5p/ KLK12 axis in lung cancer.Cell Cycle. 2019 Sep;18(18):2323-2331.
 71. Liu J,Zhang Y,Zeng Q,Zeng H,Liu X,Wu P,Xie H,He L,Long Z,Lu X,Xiao M,Zhu Y,Bo H,Cao K.Delivery of RIPK4 small interfering RNA for bladder cancer therapy using natural halloysite nanotubes.Sci Adv. 2019 Sep 25;5(9):eaaw6499.
 72. Gao F,Wu J,Gao H,Hu X,Liu L,Midgley AC,Liu Q,Sun Z,Liu Y,Ding D,Wang Y,Kong D,Huang X.Hypoxia-tropic nanozymes as oxygen generators for tumor-favoring theranostics.Biomaterials. 2020 Feb;230:119635.
 73. Zhan H,Zhao H,Muhammad N,Li T,Liu Y,Wang J.Lytic peptide-grafted beta-cyclodextrin polymer based nano-scaled drug delivery system with enhanced camptothecin anti-cancer efficacy.Nanotechnology. 2020 Feb 7;31(7):075101.

Version 2021.09.01