碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

NADP⁺/NADPH检测试剂盒(WST-8法)

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$$[\text{NADP}^+] = [\text{NADP}_{\text{total}}] - [\text{NADPH}]$$

$$[\text{NADP}^+]/[\text{NADPH}] = ([\text{NADP}_{\text{total}}] - [\text{NADPH}])/[\text{NADPH}]$$

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- Ma Y,Qi Y,Wang L,Zheng Z,Zhang Y,Zheng J. SIRT5-mediated SDHA desuccinylation promotes clear cell renal cell carcinoma tumorigenesis. FREE RADICAL BIO MED. 2019 Apr;134:458-467.
- Xie GH,Dai HJ,Liu F,Zhang YP,Zhu L,Nie JJ,Wu JH. A Dual Role of ATM in Ischemic Preconditioning and Ischemic Injury. Cell Mol Neurobiol. 2019 Dec 16;
- Jialei Sun,Chenhao Zhou,Yue Zhao,Xiaofei Zhang,Wanyong Chen,Qiang Zhou,Bo Hu,Dongmei Gao,Lisa Raatz,Zhefang Wang,Peter J Nelson,Yuchao Jiang,Ning Ren,Christiane J Bruns,Haijun Zhou. Quiescin sulfhydryl oxidase 1 promotes sorafenib-induced ferroptosis in hepatocellular carcinoma by driving EGFR endosomal trafficking and inhibiting NRF2 activation. Redox Biol. 2021 May;41:101942.
- Fenglan Qin,Huige Zhou,Jiayang Li,Jing Liu,Yaling Wang,Ru Bai,Shihui Liu,Manman Ma,Tao Liu,Fene Gao,Peiyao Du,Xiaoquan Lu,Chunying Chen. Hypoxia and pH co-triggered oxidative stress amplifier for tumor therapy. Eur J Pharmacol. 2021 Aug 15;905:174187.
- Yi Wen,Hansen Chen,Lu Zhang,Meiling Wu,Feng Zhang,Dan Yang,Jiangang Shen,Jianping Chen. Glycyrrhetic acid induces oxidative/nitrative stress and drives ferroptosis through activating NADPH oxidases and iNOS, and depriving glutathione in triple-negative breast cancer cells. Free Radic Biol Med. 2021 Sep;173:41-51.
- Mingming Sun,Hao Sheng,Tingfeng Wu,Jiaqi Song,Huanran Sun,Yingzhi Wang,Jiyan Wang,Zhen Li,Huifang Zhao,Junzhen Tan,Yanping Li,Guo Chen,Qingrong Huang,Yuan Zhang,Bei Lan,Shuangping Liu,Changliang Shan,Shuai Zhang. PIKE-A promotes glioblastoma growth by driving PPP flux through increasing G6PD expression mediated by phosphorylation of STAT3. Biochem Pharmacol. 2021 Oct;192:114736.
- Yan Wang,Haojie Jin,Yafang Wang,Yang Yao,Cuixia Yang,Jihong Meng,Xiaomu Tan,Yu Nie,Lixiang Xue,Baohui Xu,Heng Zhao,Feng Wang. Sult2b1 deficiency exacerbates ischemic stroke by promoting pro-inflammatory macrophage polarization in mice. Theranostics. 2021 Nov 1;11(20):10074-10090.
- Jian Luo,Qiang He,Jin-Zhi Xu,Chen Xu,Yin-Ze Han,Hai-Long Gao,Xian-Zhi Meng,Guo-Qing Pan,Tian Li,Ze-Yang Zhou. Microsporidia infection upregulates host energy metabolism but maintains ATP homeostasis. J Invertebr Pathol. 2021 Nov;186:107596.
- Ting Liu,Qianqian Ma,Wenjie Li,Yan Hu,Jun Yang,Qi Yao. Ubiquilin 1 suppresses the cancer stem cell-like traits of non-small cell lung cancer cells by regulating reactive oxygen species homeostasis. Bioengineered. 2021 Dec;12(1):7143-7155.