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	Apoptosis	Necroptosis	Pyroptosis
Trigger	Death receptor activation	Death receptor activation when caspase-8 is	Activation of inflammasomes by

		inhibited	pathogens
Inflammatory response	No	Yes	Yes
Key proteins	Caspase family	RIPK1、RIPK3、MLKL	Caspases 1、4、5 and 11、Gasdermin D
Morphology	Cell rounding、blebbing、formation of apoptotic bodies	Organ swelling、membrane rupture、release of cell contents	Plasma membrane rupture、release of cell contents、maintained mitochondrial integrity
Programmed cell death	Yes	Yes	Yes

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1. Chen X,Hao K,Yu X,Huang A,Zhu B,Wang GX,Ling F. Magnolol protects Ctenopharyngodon idella kidney cells from apoptosis induced by grass carp reovirus. FISH SHELLFISH IMMUN. 2018 Mar;74:426-435.
2. Lei Q,Tan J,Yi S,Wu N,Wang Y,Wu H. Mitochondrial acid 5 activates the MAPK-ERK-yap signaling pathways to protect mouse microglial BV-2 cells against TNF $\alpha$ -induced apoptosis via increased Bnip3-related mitophagy. Cell Mol Biol Lett. 2018 Apr 5;23:14.
3. Wang H,Zhao X,Ni C,Dai Y,Guo Y. Zearalenone regulates endometrial stromal cell apoptosis and migration via the promotion of mitochondrial fission by activation of the JNK/Drp1 pathway. Mol Med Rep. 2018 Jun;17(6):7797-7806.
4. Wang Y,Xu H,Lu Z,Yu X,Lv C,Tian Y,Sui D. Pseudo-Ginsenoside Rh2 induces A549 cells apoptosis via the Ras/Raf/ERK/p53 pathway. Exp Ther Med. 2018 Jun;15(6):4916-4924.
5. Zhu G,Zhang W,Liu Y,Wang S. miR-371b-5p inhibits endothelial cell apoptosis in monocrotaline-induced pulmonary arterial hypertension via PTEN/PI3K/Akt signaling pathways. Mol Med Rep. 2018 Dec;18(6):5489-5501.
6. Chen C,Wang S,Liu P. Deferoxamine Enhanced Mitochondrial Iron Accumulation and Promoted Cell Migration in Triple-Negative MDA-MB-231 Breast Cancer Cells Via a ROS-Dependent Mechanism. Int J Mol Sci. 2019 Oct 8;20(19). pii: E4952.
7. Wen Chen, Shuangfeng Chen, Chenghua Yan, Yaguang Zhang, Ronghua Zhang, Min Chen, Shufen Zhong, Weiguo Fan, Songling Zhu, Danyan Zhang, Xiao Lu, Jia Zhang, Yuying Huang, Lin Zhu, Xuezhen Li, Dawei Lv, Yadong Fu, Houkun Iv, Zhiyang Ling, Liyan Ma, Hai Jiang, Gang Long, Jinfang Zhu, Dong Wu, Bin Wu, Bing Sun. Allergen protease-activated stress granule assembly and gasdermin D fragmentation control interleukin-33 secretion. Nat Immunol. 2022 Jul;23(7):1021-1030.
8. Feng He, Yuanjun Ma, Shi Li, Haozhe Ren, Qian Liu, Xiaohua Chen, Hui Miao, Tao Ye, Qian Lu, Zuge Yang, Tianle Li, Xin Tong, Hongxu Yang, Mian Zhang, Helin Wang, Yazhou Wang, Shibin Yu. Necroptotic TNF $\alpha$ -Syndecan 4-TNF $\alpha$  Vicious Cycle as a Therapeutic Target for Preventing Temporomandibular Joint Osteoarthritis. J Bone Miner Res. 2022 May;37(5):1044-1055.

Version 2024.03.12