碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Zhang C, Wu R, Zhu H, Hu YZ, Jiang H, Lin NM, He QJ, Yang B. Enhanced anti-tumor activity by the combination of TRAIL/Apo-2L and combretastatin A-4 against human colon cancer cells via induction of apoptosis in vitro and in vivo. *Cancer Lett.* 2011 Mar 1;302(1):11-9.
2. Xu X, Gao X, Jin L, Bhadury PS, Yuan K, Hu D, Song B, Yang S. Antiproliferation and cell apoptosis inducing bioactivities of constituents from *Dysosma versipellisin* PC3 and Bcap-37 cell lines. *Cell Div.* 2011 Jun 15;6(1):14.
3. Liu M, Yang S, Jin L, Hu D, Wu Z, Yang S. Chemical Constituents of the Ethyl Acetate Extract of *Belamcanda chinensis* (L.) DC Roots and Their Antitumor Activities. *Molecules.* 2012 May 24;17(5):6156-69.
4. Peng YR, Ding YF, Wei YJ, Shu B, Li YB, Liu XD. Caudatin-2,6-dideoxy-3-O-methyl-β-D-cymaropyranoside 1 induced apoptosis through caspase 3-dependent pathway in human hepatoma cell line SMMC7721. *Phytother Res.* 2011 May;25(5):631-7.
5. Liu MC, Yang SJ, Jin LH, Hu DY, Xue W, Song BA, Yang S. Synthesis and cytotoxicity of novel ursolic acid derivatives containing an acyl piperazine moiety. *Eur J Med Chem.* 2012 Dec;58:128-35.
6. Liu L, Wang P, Liang C, He D, Yu Y, Liu X. Distinct effects of Nampt inhibition on mild and severe models of lipopolysaccharide-induced myocardial impairment. *Int Immunopharmacol.* 2013 Oct; 17(2):342-9.
7. Zhu Z, Li K, Xu D, Liu Y, Tang H, Xie Q, Xie L, Liu J, Wang H, Gong Y, Hu Z, Zheng J. Erratum to: ZFX regulates glioma cell proliferation and survival in vitro and in vivo. *J Neurooncol.* 2013 Jul;113(3):533.
8. Qi R, Yang F, Huang J, Peng H, Liu Y, Liu Z. Supplementation with conjugated linoleic acid decreases pig back fat deposition by inducing adipocyte apoptosis. *BMC Vet Res.* 2014 Jun 26;10:141.
9. Wang N, Wei H, Yin D, Lu Y, Zhang Y, Jiang D, Jiang Y, Zhang S. Cyclin D1b overexpression inhibits cell proliferation and induces cell apoptosis in cervical cancer cells in vitro and in vivo. *Int J Clin Exp Pathol.* 2014 Jun 15;7(7):4016-23.
10. You H, Li T, Zhang J, Lei Q, Tao X, Xie P, Lu W. Reduction in ischemic cerebral infarction is mediated through golgi phosphoprotein 3 and Akt/mTOR signaling following salvianolate administration. *Curr Neurovasc Res.* 2014 May;11(2):107-13.
11. Wu Z, Kong X, Zhang T, Ye J, Fang Z, Yang X. Pseudoephedrine/ephedrine shows potent anti-inflammatory activity against TNF-α-mediated acute liver failure induced by lipopolysaccharide/D-galactosamine. *Eur J Pharmacol.* 2014 Feb 5; 724: 112-21.
12. Hao R, Hu X, Wu C, Li N. Hypoxia-induced miR-15a promotes mesenchymal ablation and adaptation to hypoxia during lung development in chicken. *PLoS One.* 2014 Jun 2;9(6):e98868.
13. Shi Y, Song Q, Yu S, Hu D, Zhuang X. Microvascular invasion in hepatocellular carcinoma overexpression promotes cell proliferation and inhibits cell apoptosis of hepatocellular carcinoma via inhibiting miR-199a expression. *Oncotargets Ther.* 2015 Aug 27;8:2303-10.
14. Ke J, Yao YL, Zheng J, Wang P, Liu YH, Ma J, Li Z, Liu XB, Li ZQ, Wang ZH, Xue YX. Knockdown of long non-coding RNA HOTAIR inhibits malignant biological behaviors of human glioma cells via modulation of miR-326. *Oncotarget.* 2015 Sep 8;6(26):21934-49.
15. Lu C, Jiang Q, Hu M, Tan C, Yu H, Hua Z. Preliminary biological evaluation of ¹⁸F-FBEM-Cys-Annex V a novel apoptosis imaging agent. *Molecules.* 2015 Mar 17;20(3):4902-14.
16. Xu P, Meng Q, Sun H, Yin Q, Yu H, Zhang Z, Cao M, Zhang Y, Li Y.

- Shrapnel nanoparticles loading docetaxel inhibit metastasis and growth of breast cancer. *Biomaterials*. 2015 Sep;64:10-20.
17. Li Z, Yin H, Hao S, Wang L, Gao J, Tan X, Yang Z. miR-200 family promotes podocyte differentiation through repression of RSAD2. *Sci Rep*. 2016 Jun 2;6:27105.
 18. Han R, Wang R, Zhao Q, Han Y, Zong S, Miao S, Song W, Wang L. Trim69 regulates zebrafish brain development by ap-1 pathway. *Sci Rep*. 2016 Apr 6;6:24034.
 19. Meng C, Song L, Wang J, Li D, Liu Y, Cui X. Propofol induces proliferation partially via down regulation of p53 protein and promotes migration via activation of the Nrf2 pathway in human breast cancer cell line MDA-MB-231. *Oncol Rep*. 2017 Feb;37(2):841-848.
 20. Liu W, Gao Y, Chang N. Nurr1 overexpression exerts neuroprotective and anti-inflammatory roles via down-regulating CCL2 expression in both *in vivo* and *in vitro* Parkinson's disease models. *Biochem Biophys Res Commun*. 2017 Jan 22;482(4):1312-1319.
 21. Hu SY, Zhang Y, Zhu PJ, Zhou H, Chen YD. Liraglutide directly protects cardiomyocytes against reperfusion injury possibly via modulation of intracellular calcium homeostasis. *J Geriatr Cardiol*. 2017 Jan;14(1):57-66.
 22. Meng C, Song L, Wang J, Li D, Liu Y, Cui X. Propofol induces proliferation partially via downregulation of p53 protein and promotes migration via activation of the Nrf2 pathway in human breast cancer cell line MDA-MB-231. *Oncol Rep*. 2017 Feb;37(2):841-848.
 23. Yan Y, Wang L, He J, Liu P, Lv X, Zhang Y, Xu X, Zhang L, Zhang Y. Synergy with interferon-lambda 3 and sorafenib suppresses hepatocellular carcinomaproiferation. *Biomed Pharmacother*. 2017 Apr;88:395-402.
 24. Geng YD, Zhang C, Lei JL, Yu P, Xia YZ, Zhang H, Yang L, Kong LY. Walsuronoid B induces mitochondrial and lysosomal dysfunction leading to apoptotic ratherthan autophagic cell death via ROS/p53 signaling pathways in liver cancer. *Biochem Pharmacol*. 2017 Oct 15;142:71-86.
 25. Lu C, Jiang Q, Hu M, Tan C, Yu H, Hua Z. Preliminary biological evaluation of 18F-AIF-NOTA-MAL-Cys-Annexin V as a novel apoptosisimaging agent. *Oncotarget*. 2017 Apr 10;8(31):51086-51095.
 26. Liu W, Gao Y, Chang N. Nurr1 overexpression exerts neuroprotective and anti-inflammatory roles via down-regulatingCCL2 expression in both *in vivo* and *in vitro* Parkinson's disease models. *Biochem Biophys Res Commun*. 2017 Jan 22;482(4):1312-1319.
 27. Wang J, Wang WN, Xu SB, Wu H, Dai B, Jian DD, Yang M, Wu YT, Feng Q, Zhu JH, Zhang L, Zhang L. MicroRNA-214-3p:A link between autophagy and endothelial cell dysfunction in atherosclerosis. *Acta Physiol (Oxf)*. 2017 Sep 9.
 28. Zhao K, Li G, Yao Y, Zhou Y, Li Z, Guo Q, Lu N. Activation of phospholipase C-γ 1 and translocation of phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase contribute to GL-V9-induced apoptosis in human gastric cancer cells. *Exp Cell Res*. 2017 Jul 1;356(1):8-19.
 29. Pan T, Zhang M, Zhang F, Yan G, Ru Y, Wang Q, Zhang Y, Wei X, Xu X, Shen L, Zhang J, Wu K, Yao L, Li X. NDRG2 overexpression suppresses hepatoma cells survival during metabolic stress through disturbing the activation of fatty acid oxidation. *Biochem Biophys Res Commun*. 2017 Feb 5;483(2):860-866.
 30. Song J, Li Y, Song J, Hou F, Liu B, Li A. Mangiferin protects mitochondrial function by preserving mitochondrial hexokinase-II in vesselendothelial cells. *Biochim Biophys Acta*. 2017 Jul;1863(7):1829-1839.
 31. Li R, Ma J, Wu Y, Nangle M, Zou S, Li Y, Yin J, Zhao Y, Xu H, Zhang H, Li X, Ye QS, Wang J, Xiao J. Dual Delivery of NGF and bFGF Coacervater Ameliorates Diabetic Peripheral Neuropathy viaInhibiting Schwann Cells Apoptosis. *Int J Biol Sci*. 2017 May 16;13(5):640-651.
 32. Liu ZY, Hu S, Zhong QW, Tian CN, Ma HM, Yu JJ. N-Methyl-D-Aspartate Receptor-Driven Calcium Influx Potentiates the Adverse Effects of Myocardial Ischemia-Reperfusion Injury Ex Vivo. *J Cardiovasc Pharmacol*. 2017 Nov;70(5):329-338.
 33. Dai H, Zhang X, Yang Z, Li J, Zheng J. Effects of Baicalin on Blood Pressure and Left Ventricular Remodeling in Rats with RenovascularHypertension. *Med Sci Monit*. 2017 Jun 16;23:2939-2948.
 34. Chen D, Tao X, Zhou L, Sun F, Sun M, Fang X. Spaghetti,a homolog of human RPAP3 (RNA polymerase II-associated protein 3), determines the fate of germline stem cells in *Drosophila* ovary. *Cell Biol Int*. 2017 Nov 7.
 35. Zhang X, Wang X, Liu T, Mo M, Ao L, Liu J, Cao J, Cui Z. ZnSO4 rescued vimentin from collapse in DBP-exposed Sertoli cells by attenuating ER stress and apoptosis. *Toxicol In Vitro*. 2018 Apr;48:195-204.
 36. Zheng X, Bi W, Yang G, Zhao J, Wang J, Li X, Zhou X. Hyperglycemia Induced by Chronic Restraint Stress in Mice Is Associated With Nucleus Tractus Solitarius Injury and Not Just the Direct Effect of Glucocorticoids. *Front Neurosci*. 2018 Dec 19;12:983.
 37. Huang Y, Wu S, Zhang Y, Wang L, Guo Y. Antitumor effect of triptolide in T-cell lymphoblastic lymphoma by inhibiting cell viability, invasion, and epithelial-mesenchymal transition via regulating the PI3K/AKT/mTOR pathway. *Onco Targets Ther*. 2018 Feb 12;11:769-779.
 38. Jiang C, Yang Y, Yang Y, Guo L, Huang J, Liu X, Wu C, Zou J. Long Noncoding RNA (lncRNA) HOTAIR Affects Tumorigenesis and Metastasis of Non-Small Cell Lung Cancer by Upregulating miR-613. *Oncol Res*. 2018 Jun 11;26(5):725-734.
 39. Feng Q, He P, Wang Y. MicroRNA-223-3p regulates cell chemo-sensitivity by targeting FOXO3 in prostatic cancer. *Gene*. 2018 Jun 5;658:152-158.
 40. Luo J, Zhu H, Jiang H, Cui Y, Wang M, Ni X, Ma C. The effects of aberrant expression of lncRNA DGCR5/miR-873-5p/TUSC3 in lung cancer cell progression. *Cancer Med*. 2018 May 23.
 41. Liu J, Zhou S, Zhang Y, Li X, Qian X, Tao W, Jin L, Zhao J. Bax inhibitor-1 suppresses early brain injury following experimental subarachnoid hemorrhage in rats. *Int J Mol Med*. 2018 Nov;42(5):2891-2902.
 42. Liu Z, Li H, Hong C, Chen M, Yue T, Chen C, Wang Z, You Q, Li C, Weng Q, Xie H, Hu R. ALS-Associated E478G Mutation in Human OPTN (Optineurin) Promotes Inflammation and Induces Neuronal Cell Death. *Front Immunol*. 2018 Nov 14;9:2647.
 43. Tang P, Dang H, Huang J, Xu T, Yuan P, Hu J, Sheng JF. NADPH oxidase NOX4 is a glycolytic regulator through mROS-HIF1 α axis in thyroid carcinomas. *Sci Rep*. 2018 Oct 26;8(1):15897.
 44. Lou Z, Wang AP, Duan XM, Hu GH, Zuo ML, Yang ZB. Role of ALK5/SMAD2/3 signaling in the regulation of NOX expression in cerebral ischemia/reperfusion injury. *Exp Ther Med*. 2018 Sep;16(3):1671-1678.
 45. Zhang H, Deng Z, Zhang D, Li H, Zhang L, Niu J, Zuo W, Fu R, Fan L, Ye JH, She J. High expression of leucine-rich repeat-containing 8A is indicative of a worse outcome of colon cancer patients by enhancing cancer cell growth and metastasis. *Oncol Rep*. 2018 Sep;40(3):1275-1286.
 46. Zhang G, Wang Z, Chen W, Cao Y, Wu J, Qiang G, Ji A, Wu J, Jiang C. Dual effects of gossypol on human hepatocellular carcinoma via endoplasmic reticulum stress and autophagy. *Int J Biochem Cell Biol*. 2019 Aug;113:48-57.
 47. Zheng Y, Zhang P, Zhang C, Zeng W. Surgery-induced cryptorchidism induces apoptosis and autophagy of spermatogenic cells in mice. *Zygote*. 2019 Apr;27(2):101-110.
 48. Yuan Y, Wang YY, Liu X, Luo B, Zhang L, Zheng F, Li XY, Guo LY, Wang L, Jiang M, Pan YM, Yan YW, Yang JY, Chen SY, Wang JN, Tang JM. KPC1 alleviates hypoxia/reoxygenation-induced apoptosis in rat cardiomyocyte cells though BAX degradation. *J Cell Physiol*. 2019 Dec;234(12):22921-22934.
 49. Xiao Q, Zeng FL, Tang GY, Lei CY, Zou XX, Liu XL, Peng BL, Qin S, Li HX. Expression of galectin-3 and apoptosis in placental villi from patients with missed abortion during early pregnancy. *Exp Ther Med*. 2019 Apr;17(4):2623-2631.
 50. Li W, Qinghai S, Kai L, Xue M, Lili N, Jihua R, Zhengxiang L, Xiaoling L, Di G, Qi Y, Mengyun D, Jianfeng F. Oral administration of Ginkgolide B alleviates hypoxia-induced neuronal damage in rat hippocampus by inhibiting oxidative stress and apoptosis. *Iran J Basic Med Sci*. 2019 Feb;22(2):140-145.
 51. Li J, Tan YJ, Wang MZ, Sun Y, Li GY, Wang QL, Yao JC, Yue J, Liu Z, Zhang GM, Ren YS. Loganatin protects against rhabdomyolysis-induced acute kidney injury by modulating the toll-like receptor 4 signalling pathway. *Br J Pharmacol*. 2019 Apr;176(8):1106-1121.
 52. Ma W, Shen H, Li Q, Song H, Guo Y, Li F, Zhou X, Guo X, Shi J, Cui Q, Xing J, Deng J, Yu Y, Liu W, Zhao Y, HMARVELD1 attenuates arsenic trioxide-induced apoptosis in liver cancer cells by inhibiting reactive oxygen species production. *Ann Transl Med*. 2019 May;7(9):200.
 53. Nie C, Qin X, Li X, Tian B, Zhao Y, Jin Y, Li Y, Wang Q, Zeng D, Hong A, Chen X. CACNA2D3 Enhances the Chemosensitivity of Esophageal Squamous Cell Carcinoma to Cisplatin via Inducing Ca2+-Mediated Apoptosis and Suppressing PI3K/Akt Pathways. *Front Oncol*. 2019 Apr 2;9:185.

54. Yang Y,Chen A, Ma J,Wu A,Xu F.Effects of radioactive ^{125}I on apoptosis of HGC-27 gastric cancer cells.Oncol Lett. 2019 Nov;18(5):4916-4922.
55. Zhang J,Fan J,Li Y,Liang S,Huo S,Wang X,Zuo Y,Cui D,Li W,Zhong Z,Zhong F.Porcine Parvovirus Infection Causes Pig Placenta Tissue Damage Involving Nonstructural Protein 1 (NS1)-Induced Intrinsic ROS/Mitochondria-Mediated Apoptosis.Viruses. 2019 Apr 25;11(4). pii: E389.
56. Ngkure M X K,Jiang J,Enayatullah H,Ennab W,Mustafa S,Roden S,Wei Q,Shi F.Sweet taste receptor agonists alter ovarian functions and ovarian cycles in aged mice.Reprod Biol. 2019 Sep;19(3):230-236.
57. Liu M,Ai W,Sun L,Fang F,Wang X,Chen S,Wang H.Triclosan-induced liver injury in zebrafish (*Danio rerio*) via regulating MAPK/p53 signaling pathway.Comp Biochem Physiol C Toxicol Pharmacol. 2019 Aug;222:108-117.
58. Liang M,Jia J,Chen L,Wei B,Guan Q,Ding Z,Yu J,Pang R,He G.LncRNA MCM3AP-AS1 promotes proliferation and invasion through regulating miR-211-5p/SPARC axis in papillary thyroid cancer.Endocrine. 2019 Aug;65(2):318-326
59. Liang Q,Li F,Li Y,Liu Y,Lan M,Wu S,Wu X,Ji Y,Zhang R,Yin L..Self-assisted membrane-penetrating helical polypeptides mediate anti-inflammatory RNAi against myocardial ischemic reperfusion (IR) injury.Biomater Sci. 2019 Aug 20;7(9):3717-3728.
60. Wang H,Yu D,Fang J,Cao C,Liu Z,Ren J,Qu X.Renal-Clearable Porphyrinic Metal-Organic Framework Nanodots for Enhanced Photodynamic Therapy.ACS Nano. 2019 Aug 27;13(8):9206-9217.
61. Wang A,Zhao C,Gao Y,Duan G,Yang Y,Fan B,Wang X,Wang K.LEF1-AS1 contributes to proliferation and invasion through regulating miR-544a/ FOXP1 axis in lung cancer.Invest New Drugs. 2019 Dec;37(6):1127-1134.
62. Tang KK,Liu XY,Wang ZY,Qu KC,Fan RF.Trehalose alleviates cadmium-induced brain damage by ameliorating oxidative stress, autophagy inhibition, and apoptosis.Metallomics. 2019 Dec 11;11(12):2043-2051.
63. Tang Y,Cao G,Zhao G,Wang C,Qin Q.LncRNA differentiation antagonizing non-protein coding RNA promotes proliferation and invasion through regulating miR-135a/NLRP37 axis in pancreatic cancer.Invest New Drugs. 2019 Jul 3.
64. Zheng Y,Liu W,Chen Y,Li C,Jiang H,Wang X.Conjugating gold nanoclusters and antimicrobial peptides: From aggregation-induced emission to antibacterial synergy.J Colloid Interface Sci. 2019 Jun 15;546:1-10.
65. Zheng S,Li M,Miao K,Xu H.SNHG1 contributes to proliferation and invasion by regulating miR-382 in breast cancer.Cancer Manag Res. 2019 Jun 18;11:5589-5598.
66. An W,Yu Y,Zhang Y,Zhang Z,Yu Y,Zhao X.Exogenous IL-19 attenuates acute ischaemic injury and improves survival in male mice with myocardial infarction.Br J Pharmacol. 2019 Mar;176(5):699-710.
67. Huang S,Qi P,Zhang T,Li F,He X.The HIF-1 α /miR-224-3p/ATG5 axis affects cell mobility and chemosensitivity by regulating hypoxia-induced protective autophagy in glioblastoma and astrocytoma.Oncol Rep. 2019 Mar;41(3):1759-1768.
68. Wu G,Zhu Q,Zeng J,Gu X,Miao Y,Xu W,Lv T,Song Y.Extracellular mitochondrial DNA promote NLRP3 inflammasome activation and induce acute lung injury through TLR9 and NF- κ B.J Thorac Dis. 2019 Nov;11(11):4816-4828.
69. Yang W,Hao Y,Hou W,Fang X,Fang P,Jiang T, Ma C,Zhang QJieduan-Niwan Formula Reduces Liver Apoptosis in a Rat Model of Acute-on-Chronic Liver Failure by Regulating the E2F1-Mediated Intrinsic Apoptosis Pathway.Evid Based Complement Alternat Med. 2019 Nov 11;2019:8108503.
70. Bai Y,Zhang G,Cheng R,Yang R,Chu H.CASC15 contributes to proliferation and invasion through regulating miR-766-5p/ KLK12 axis in lung cancer.Cell Cycle. 2019 Sep;18(18):2323-2331.
71. Liu J,Zhang Y,Zeng Q,Zeng H,Liu X,Wu P,Xie H,He L,Long Z,Lu X,Xiao M,Zhu Y,Bo H,Cao K.Delivery of RIPK4 small interfering RNA for bladder cancer therapy using natural halloysite nanotubes.Sci Adv. 2019 Sep 25;5(9):eaaw6499.
72. Gao F,Wu J,Gao H,Hu X,Liu L,Midgley AC,Liu Q,Sun Z,Liu Y,Ding D,Wang Y,Kong D,Huang X.Hypoxia-tropic nanozymes as oxygen generators for tumor-favoring theranostics.Biomaterials. 2020 Feb;230:119635.
73. Zhan H,Zhao H,Muhammad N,Li T,Liu Y,Wang J.Lytic peptide-grafted beta-cyclodextrin polymer based nano-scaled drug delivery system with enhanced camptothecin anti-cancer efficacy.Nanotechnology. 2020 Feb 7;31(7):075101.
74. Jie Liu,Wang Wang,Lei Wang,Xian-Mei Qi,Yu-Hui Sha,Ting Yang.3-Bromopyruvate alleviates the development of monocrotaline-induced rat pulmonary arterial hypertension by decreasing aerobic glycolysis, inducing apoptosis, and suppressing inflammation CHINESE MED J-PEKING. 2020 Jan 5;133(1):49-60.;doi: 10.1097/CM9.0000000000000577
75. Hailin Tan,Guiming Zhang,Xuecheng Yang,Tao Jing,Daqing Shen,Xinsheng Wang,Peimine inhibits the growth and motility of prostate cancer cells and induces apoptosis by disruption of intracellular calcium homeostasis through Ca^{2+} /CaMKII/JNK pathway J Cell Biochem. 2020 Jan;121(1):81-92.;doi: 10.1002/jcb.28870
76. Gao F,Wu J,Gao H,Hu X,Liu L,Midgley AC,Liu Q,Sun Z,Liu Y,Ding D,Wang Y,Kong D,Huang X.Hypoxia-tropic nanozymes as oxygen generators for tumor-favoring theranostics.Biomaterials. 2020 Feb;230:119635
77. Zhan H,Zhao H,Muhammad N,Li T,Liu Y,Wang J.Lytic peptide-grafted beta-cyclodextrin polymer based nano-scaled drug delivery system with enhanced camptothecin anti-cancer efficacy.Nanotechnology. 2020 Feb 7;31(7):075101
78. Rui Wu,Xi Yang,Qian Zhou,Wenyan Yu,Minwei Li,Jianer Wo,Weiguang Shan,Hong Zhao,Yan Chen,Zhajun Zhan.Aurovertin B exerts potent antitumor activity against triple-negative breast cancer in vivo and in vitro via regulating ATP synthase activity and DUSP1 expression Pharmazie. 2020 Jun 1;75(6):261-265.;doi: 10.1691/ph.2020.0380
79. Yu-Shan Ren,Hong-Hua Li,Jing-Chun Yao,Yu-Jun Tan,Li-Hong Pan,Tao Peng,Li-Li Zhao,Gui-Min Zhang,Jiang Yue,Xue-Mei Hu,Zhong Liu,Jie Li.Application quantitative proteomics approach to identify differentially expressed proteins associated with cardiac protection mediated by cycloastragenol in acute myocardial infarction rats J Proteomics. 2020 Jun 30;222:103691.;doi: 10.1016/j.jprot.2020.103691
80. Xiaoping Guo,Hongkun Lin,Jingjing Liu,Dongxia Wang,Dan Li,Chunjie Jiang,Yuhan Tang,Jun Wang,Tingrui Zhang,Yanyan Li,Ping Yao.1,25-Dihydroxyvitamin D attenuates diabetic cardiac autophagy and damage by vitamin D receptor-mediated suppression of FoxO1 translocation J Nutr Biochem. 2020 Jun;80:108380.;doi: 10.1016/j.jnutbio.2020.108380
81. Xia Zhao,Shuai Li,Uma Gaur,Wenhua Zheng.Artemisinin Improved Neuronal Functions in Alzheimer's Disease Animal Model 3xtg Mice and Neuronal Cells via Stimulating the ERK/CREB Signaling Pathway Aging Dis. 2020 Jul 23;11(4):801-819.;doi: 10.14336/AD.2019.0813
82. Yi-Jun Zhou,Jian-Min Chen,Kiran Sapkota,Jiang-Yi Long,Yuan-Jun Liao,Jun-Jun Jiang,Bing-Yu Liang,Jin-Bin Wei,Yan Zhou.Panax notoginseng saponins attenuate CCL2-induced cognitive deficits in rats via anti-inflammation and anti-apoptosis effects that involve suppressing over-activation of NMDA receptors Biomed Pharmacother. 2020 Jul;127:110139.;doi: 10.1016/j.biopha.2020.110139
83. Yan Xiao,Wanying Chen,Zehao Zhong,Liang Ding,Hua Bai,Hao Chen,Hongru Zhang,Yihuang Gu,Shengfeng Lu.Electroacupuncture preconditioning attenuates myocardial ischemia-reperfusion injury by inhibiting mitophagy mediated by the mTORC1-ULK1-FUND1 pathway Biomed Pharmacother. 2020 Jul;127:110148.;doi: 10.1016/j.biopha.2020.110148
84. Guanwen Huang,Jiwen Bao,Xinghua Shao,Wenyan Zhou,Bei Wu,Zhaohui Ni,Ling Wang.Inhibiting pannexin-1 alleviates sepsis-induced acute kidney injury via decreasing NLRP3 inflammasome activation and cell apoptosis Life Sci. 2020 Aug 1;254:117791.;doi: 10.1016/j.lfs.2020.117791
85. Yajun Song,Wengang Hu,Ya Xiao,Yang Li,Xu Wang,Weifeng He,Jingming Hou,Yun Liu,Guangping Liang,Chibing Huang.Keratinocyte growth factor ameliorates mycophenolate mofetil-induced intestinal barrier disruption in mice Mol Immunol. 2020 Aug;124:61-69.;doi: 10.1016/j.molimm.2020.04.012
86. Ying Luo,Yue Tian,Chunrong Zhao.Taurine attenuates liver autophagy and injury of offspring in gestational diabetic mellitus rats Life Sci. 2020 Sep 15;257:117889.;doi: 10.1016/j.lfs.2020.117889
87. Zhibing Qiu,Wen Chen,Yafeng Liu,Ben Jiang,Li Yin,Xin Chen.LncRNA AC061961.2 overexpression inhibited endoplasmic reticulum stress induced apoptosis in dilated cardiomyopathy rats and cardiomyocytes via activating wnt/ β -catenin pathway J RECEPT SIG TRANSD. 2020 Oct 22;1-10.;doi: 10.1080/10799893.2020.1828915
88. Wenchong Sun,Jian Wang,Dasheng Cai,Ling Pei.Neuroprotection of the Developing Brain by Dexmedetomidine Is Mediated by Attenuating Single Propofol-induced Hippocampal Apoptosis and Synaptic Plasticity Deficits Exp Neurobiol. 2020 Oct 31;29(5):356-375.;doi:

89. Wenjing Wang,Wenxiao Ding,Hanpeng Huang,Yina Zhu,Ning Ding,Ganzhu Feng,Xilong Zhang.The role of mitophagy in the mechanism of genioglossal dysfunction caused by chronic intermittent hypoxia and the protective effect of adiponectin *Sleep Breath.* 2020 Oct 8.;doi: 10.1007/s11325-020-02211-0
90. Jun Zhang,Li Liu,Furong Li,Zongqian Wang,Jinghong Zhao.Treatment with catalpol protects against cisplatin-induced renal injury through Nrf2 and NF- κ B signaling pathways *Exp Ther Med.* 2020 Oct;20(4):3025-3032.;doi: 10.3892/etm.2020.9077
91. Bin Huai,Jiyu Ding,Atractylenolide III attenuates bleomycin-induced experimental pulmonary fibrosis and oxidative stress in rat model via Nrf2/NQO1/HO-1 pathway activation *IMMUNOPHARM IMMUNOT.* 2020 Oct;42(5):436-444.;doi: 10.1080/08923973.2020.1806871
92. Xiaoli Li,Wenhui Chen,Jinfang Feng,Bo Zhao.The effects of HIF-1 α overexpression on renal injury, immune disorders and mitochondrial apoptotic pathways in renal ischemia/reperfusion rats *Transl Androl Urol.* 2020 Oct;9(5):2157-2165.;doi: 10.21037/tau-20-918
93. Wenjing Liu,Yang Liu,Zhaohai Pan,Xin Zhang,Yao Qin,Xiaojie Chen,Minjing Li,Xiaoyu Chen,Qiusheng Zheng,Xiaona Liu,Defang Li.Systematic Analysis of rRNA-Derived Small RNAs Discloses New Therapeutic Targets of Caloric Restriction in Myocardial Ischemic Rats *Front Cell Dev Biol.* 2020 Nov 3;8:568116.;doi: 10.3389/fcell.2020.568116
94. Zhezhe Sun,Ke Wu,Lei Gu,Lijie Huang,Qichuan Zhuge,Su Yang,Zhenzhong Wang.IGF-1R stimulation alters microglial polarization via TLR4/NF- κ B pathway after cerebral hemorrhage in mice *Brain Res Bull.* 2020 Nov;164:221-234.;doi: 10.1016/j.brainresbull.2020.08.026
95. Dejun Xu,Lingbin Liu,Yongju Zhao,Li Yang,Jianyong Cheng,Rongmao Hua,Zelin Zhang,Qingwang Li.Melatonin protects mouse testes from palmitic acid-induced lipotoxicity by attenuating oxidative stress and DNA damage in a SIRT1-dependent manner *J Pineal Res.* 2020 Nov;69(4):e12690.;doi: 10.1111/jpi.12690
96. Zhenrong Liu,Yan Huang,Fan Zhang,Huirong Tang,Youjuan Wang.Autophagy-related gene 7 deficiency caused by miR-154-5p overexpression suppresses the cell viability and tumorigenesis of retinoblastoma by increasing cell apoptosis *Ann Transl Med.* 2020 Nov;8(21):1451.;doi: 10.21037/atm-20-6009
97. Hua Xue,Manxiang Li.Protective effect of pterostilbene on sepsis-induced acute lung injury in a rat model via the JAK2/STAT3 pathway *Ann Transl Med.* 2020 Nov;8(21):1452.;doi: 10.21037/atm-20-5814
98. Zhilong Zhang,Yuting Zhao,Yanwen Zhang,Ruqian Zhao,Bin He.Paternal systemic inflammation induces offspring programming of growth and liver regeneration in association with Igf2 upregulation *Mol Cell Endocrinol.* 2020 Dec 1;518:111001.;doi: 10.1016/j.mce.2020.111001
99. Yixiong Wang,Baozhu Xu,Jianying Zhou,Xianyan Wu.Propofol activates AMPK to inhibit the growth of HepG2 cells in vitro and hepatocarcinogenesis in xenograft mouse tumor models by inducing autophagy *J Gastrointest Oncol.* 2020 Dec;11(6):1322-1332.;doi: 10.21037/jgo-20-472
- 100.Qiushi Xu,Yunhui Fan,Juan J Loor,Yusheng Liang,Hongming Lv,Xudong Sun,Hongdou Jia,Chuang Xu,Aloin protects mice from diet-induced non-alcoholic steatohepatitis via activation of Nrf2/HO-1 signaling *Food Funct.* 2021 Jan 21;12(2):696-705.;doi: 10.1039/d0fo02684k
- 101.Yuan Tian,Yue-Hai Xiao,Tian Geng,Chao Sun,Jiang Gu,Kai-Fa Tang,Bei Liu,Yi-Meng Liu,Fa Sun.Clusterin suppresses spermatogenic cell apoptosis to alleviate diabetes-induced testicular damage by inhibiting autophagy via the PI3K/AKT/mTOR axis *Biol Cell.* 2021 Jan;113(1):14-27.;doi: 10.1111/boc.202000030
- 102.Zhi Yang,Weigang Wu,Pengcheng Ou,Minna Wu,Furong Zeng,Boping Zhou,Shipin Wu.Mir-122-5p knockdown protects against APAP-mediated liver injury through up-regulating NDRG3 *Mol Cell Biochem.* 2021 Feb;476(2):1257-1267.;doi: 10.1007/s11010-020-03988-0
- 103.Zhe Xu,Xiaoguang Liu,Yanqing Li,Hongliang Gao,Tao He,Chunlei Zhang,Wei Hao,Xu Teng.Shuxuetong injection simultaneously ameliorates dexamethasone-driven vascular calcification and osteoporosis *Exp Ther Med.* 2021 Mar;21(3):197.;doi: 10.3892/etm.2021.9630
- 104.Xia Xu,Hua Feng,Chaochao Dai,Weida Lu,Jun Zhang,Xiaosun Guo,Qihui Yin,Jianli Wang,Xiaopei Cui,Fan Jiang.Therapeutic efficacy of the novel selective RNA polymerase I inhibitor CX-5461 on pulmonary arterial hypertension and associated vascular remodelling *BRIT J PHARMACOL.* 2021 Apr;178(7):1605-1619.;doi: 10.1111/bph.15385
- 105.Yamin Wang,Fengrong Liu,Peng Liu.23-Hydroxytormentic acid reduces cerebral ischemia/reperfusion damage in rats through anti-apoptotic, antioxidant, and anti-inflammatory mechanisms *N-S ARCH PHARMACOL.* 2021 May;394(5):1045-1054.;doi: 10.1007/s00210-020-02038-2