碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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- Li J, Xia X, Ke Y, Nie H, Smith MA, Zhu X. Trichosanthin induced apoptosis in HL-60 cells via mitochondrial and endoplasmic reticulum stress signaling pathways. *Biochim Biophys Acta.* 2007 Aug; 1770(8):1169-80.
- Chen R, Jin R, Wu L, Ye X, Yang Y, Luo K, Wang W, Wu D, Ye X, Huang L, Huang T, Xiao G. Reticulon 3 attenuates the clearance of cytosolic prion aggregates via inhibiting autophagy. *Autophagy.* 2011 Feb; 7(2):205-16.
- Chen X, Ran ZH, Tong JL, Nie F, Zhu MM, Xu XT, Xiao SD. RNA interference (RNAi) of Ufd1 protein can sensitize a hydroxycamptothecin-resistant coloncancer cell line SW1116/HCPT to hydroxycamptothecin. *J Dig Dis.* 2011 Apr; 12(2):110-6.
- Ning B, Bai M, Shen W. Reduced glutathione protects human hepatocytes from palmitate-mediated injury by suppressing endoplasmic reticulum stress response. *Hepatogastroenterology.* 2011 Sep-Oct; 58(110-111):1670-9.
- Xiao G, Zhou J, Lu X, Huang R, Zhang H. Excessive UDPG resulting from the mutation of UAP1 causes programmed cell death by triggering reactive oxygen species accumulation and caspase-like activity in rice. *New Phytol.* 2018 Jan; 217(1):332-343.
- Xiao G, Zhou J, Lu X, Huang R, Zhang H. Excessive UDPG resulting from the mutation of UAP1 causes programmed cell death by triggering reactive oxygen species accumulation and caspase-like activity in rice. *New Phytol.* 2018 Jan; 217(1):332-343.

Version 2021.09.01