



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

组织线粒体分离试剂盒

产品编号	产品名称	包装
C3606	组织线粒体分离试剂盒	50-100次

产品简介:

- 组织线粒体分离试剂盒(Tissue Mitochondria Isolation Kit)是用于快速便捷分离动物组织线粒体的试剂盒。
- 本试剂盒在分离线粒体的同时可以获得去除线粒体的细胞浆蛋白，可用于研究细胞色素c等线粒体蛋白向胞浆的释放。
- 使用本试剂盒分离获得的线粒体纯度较高，并且绝大部分分离获得的线粒体都含有完整的内膜和外膜，并具有线粒体的生理功能。因此本试剂盒分离得到的线粒体可以用于线粒体的生理功能等方面的研究。例如可以使用碧云天的C2006线粒体膜电位检测试剂盒(JC-1)测定分离得到的线粒体的膜电位。
- 本试剂盒分离得到的线粒体也可以被试剂盒中的线粒体裂解液或其它适当裂解液裂解后用于SDS-PAGE、Western、双向电泳等蛋白分析。
- 提供了两种不同的线粒体分离试剂，适用于不同的组织样品。
- 本试剂盒提供了用于组织消化的胰酶消化液，使对组织样品的处理更加便捷。
- 如果每个组织样品的重量为50-100mg，本试剂盒可以处理50-100个样品。

包装清单:

产品编号	产品名称	包装
C3606-1	线粒体分离试剂A	60ml
C3606-2	线粒体分离试剂B	60ml
C3606-3	胰酶消化液	50ml
C3606-4	线粒体储存液	3ml
C3606-5	线粒体裂解液	15ml
—	说明书	1份

保存条件:

-20°C保存，一年有效。其中胰酶消化液可以4°C保存。

注意事项:

- 试剂盒中的试剂对于不同的实验目的不必全部使用。
- 如果用于制备蛋白样品，需自备PMSF，否则不必使用PMSF。PMSF(ST506)可以向碧云天订购。PMSF一定要在线粒体分离试剂或线粒体裂解液加入到样品中前2-3分钟内加入，以免PMSF在水溶液中很快失效。
- 使用说明中的操作步骤按照组织重量为50-100mg进行说明，如果组织块较大，可以按比例加大各溶液的使用量。
- 分离线粒体的所有步骤均需在冰上或4°C进行，所用溶液需冰浴或4°C预冷。
- 通常在分离线粒体时前后两次离心速度选取600g和11,000g，如果希望纯度更高，但对线粒体的得率要求不高，前后两次离心速度可以采用1000g和3500g。
- 使用本试剂盒中的线粒体储存液稀释或储存的线粒体样品应及时使用，以免线粒体膜电位受影响。如果不能及时使用，建议在-80°C保存。冻存后的线粒体样品不推荐用于膜电位的检测，但可以用线粒体蛋白或核酸的相关检测。
- 如果出现本试剂盒中线粒体储存液不够用的情况，可以单独订购线粒体储存液(C3609)。
- Bradford法蛋白浓度测定试剂盒(P0006)、Bradford蛋白浓度测定试剂盒(去垢剂兼容型)(P0006C)和BCA蛋白浓度测定试剂盒(P0009/P0010/P0010S/P0011/P0012/P0012S)可以向碧云天订购。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明:

1. **准备溶液：**室温融解试剂盒中的各种溶液，溶解后立即置于冰上并混匀。第一次使用时，把1.5ml PMSF(溶剂)加入到PMSF(晶体)中，溶解并混匀，即得到1.5ml 100mM PMSF。配制好的100mM PMSF溶液-20°C保存。如果最终目的是制备线粒体蛋白，根据样品数量，取适量相应的线粒体分离试剂备用，在线粒体分离试剂加入到组织样品中前数分钟内加入PMSF，使PMSF的最终浓度为1mM。取适量线粒体裂解液备用，在线粒体裂解液加入到线粒体样品中前数分钟内加入PMSF，使PMSF的最终浓度为1mM。

2. 从软组织(例如肝脏和脑)中分离线粒体:

- a. 使用新鲜的动物组织，通常要求动物处死后不超过一小时，并且保存在冰上的组织。请勿使用经过冷冻保存的组织。
- b. 剪取一小块组织，重量约为50-100mg。在1.5ml离心管内对剪取的组织进行称重。用PBS洗涤组织一次。
- c. 把组织放在一个置于冰上的离心管或培养皿中，用剪刀或刀片把组织剪切成非常细小的组织碎片。
- d. 加入10体积预冷的线粒体分离试剂A或临用前添加了PMSF的线粒体分离试剂A，在冰浴上进行匀浆，匀浆10次左右。
注：如果组织的重量为80mg，可以大致认为组织的体积接近80微升，此时需加入800微升线粒体分离试剂A。
- e. 把匀浆在600g, 4°C离心5分钟。
注：如需获得纯度更高的线粒体，可以将此步骤的离心速度改为1000g，其它离心条件不变；获得更高纯度线粒体的缺点是相同数量细胞的线粒体抽提得率会下降。
- f. 小心把上清转移到另一离心管中，在11,000g, 4°C离心10分钟。
注：如需获得纯度更高的线粒体，可以将此步骤的离心速度改为3500g，其它离心条件不变；获得更高纯度线粒体的缺点是相同数量细胞的线粒体抽提得率会下降。
- g. 小心去除上清。沉淀即为分离得到的线粒体。
注：如果希望获得去除线粒体的细胞浆蛋白，在本步骤需收集上清，并且在收集上清时注意勿触及沉淀。随后把收集的上清12,000g, 4°C离心10分钟。取上清即为去除线粒体的细胞浆蛋白。该细胞浆蛋白如果需要测定蛋白浓度，必须使用PBS或生理盐水至少稀释2倍，然后才能用BCA法或Bradford法进行检测(须特别注意，标准品也需要配制在含有相应比例线粒体分离试剂的PBS或生理盐水中)。

3. 从硬组织(例如心肌和骨骼肌)中分离线粒体:

- a. 对于心肌组织采用线粒体分离试剂A，对于骨骼肌组织采用线粒体分离试剂B，对于其它类似组织可以优先考虑使用线粒体分离试剂A，如遇到效果不佳的情况可以试用线粒体分离试剂B。
- b. 使用新鲜的动物组织，通常要求动物处死后不超过一小时，并且保存在冰上的组织。请勿使用经过冷冻保存的组织。
- c. 剪取一小块组织，重量约为50-100mg。在1.5ml离心管内对剪取的组织进行称重。用PBS洗涤组织一次。
- d. 把组织放在一个置于冰上的离心管或培养皿中，用剪刀或刀片把组织剪切成非常细小的组织碎片。
- e. 加入10体积预冷的PBS，冰浴3分钟。
注：如果组织的重量为80mg，可以大致认为组织的体积接近80微升，此时需加入800微升PBS。
- f. 600g离心10-20秒，沉淀组织样品，弃上清。
- g. 再加入8体积预冷的胰酶消化液，冰浴20分钟。
注：如果组织的重量为80mg，此时需加入约640微升相应的胰酶消化液。
- h. 600g离心10-20秒，沉淀组织样品，弃上清。
- i. 加入2体积相应线粒体分离试剂，重悬组织，用于洗去残余的胰酶。
注：如果组织的重量为80mg，此时需加入约160微升线粒体分离试剂。
- j. 600g离心10-20秒，沉淀组织样品，弃上清。
- k. 加入8体积预冷的相应线粒体分离试剂或临用前添加了PMSF的线粒体分离试剂，在冰浴上进行匀浆，匀浆20-30次。
注：如果组织的重量为80mg，此时需加入约640微升相应的线粒体分离试剂。
- l. 把匀浆在600g, 4°C离心5分钟。
注：如需获得纯度更高的线粒体，可以将此步骤的离心速度改为1000g，其它离心条件不变；获得更高纯度线粒体的缺点是相同数量细胞的线粒体抽提得率会下降。
- m. 小心把上清转移到另一离心管中，在11,000g, 4°C离心10分钟。
注：如需获得纯度更高的线粒体，可以将此步骤的离心速度改为3500g，其它离心条件不变；获得更高纯度线粒体的缺点是相同数量细胞的线粒体抽提得率会下降。
- n. 小心去除上清。沉淀即为分离得到的线粒体。
注：如果希望获得去除线粒体的细胞浆蛋白，在本步骤需收集上清，并且在收集上清时注意勿触及沉淀。随后把收集的上清12,000g, 4°C离心10分钟。取上清即为去除线粒体的细胞浆蛋白。该细胞浆蛋白如果需要测定蛋白浓度，请参考步骤2g。

4. 线粒体的使用:

- a. 如果用于完整线粒体的功能或酶活性研究，可以加入适量线粒体储存液，重悬线粒体。通常每100mg组织获得的线粒体可以用40微升相应的线粒体储存液重悬，预期分离获得的线粒体样品的蛋白浓度为10-20mg/ml。该蛋白样品可以直接使用Bradford法测定蛋白浓度或离心沉淀后裂解再用BCA法测定蛋白浓度。
注：线粒体储存液中含有的DTT会干扰BCA法测定蛋白浓度。
- b. 如果用于线粒体的蛋白分析，初始重量为50-100mg的组织样品分离得到的线粒体样品中可以加入150-200微升临用前添加了PMSF的线粒体裂解液裂解线粒体。裂解后的线粒体可以用于PAGE、Western、IP以及线粒体中的一些酶活性的测定等。裂解后的蛋白样品可以使用BCA法或去垢剂兼容型Bradford法(P0006C)测定蛋白浓度。
- c. 如果用于双向电泳，请使用适当的用于双向电泳的裂解液处理线粒体。

使用本产品的文献:

1. Shen J, Huang C, Jiang L, Gao F, Wang Z, Zhang Y, Bai J, Zhou H, Chen Q. Enhancement of cisplatin induced apoptosis by suberoylanilide hydroxamic acid in human oral squamous cell carcinoma cell lines. *Biochem Pharmacol*. 2007 Jun 15;73(12):1901-9.
2. Zhang WC, Peng YJ, He WQ, Lv N, Chen C, Zhi G, Chen HQ, Zhu MS. Identification and functional characterization of an aggregation domain in long myosin light chain kinase. *FEBS J*. 2008 May;275(10):2489-500.
3. Chen K, Zhang Q, Wang J, Liu F, Mi M, Xu H, Chen F, Zeng K. Taurine protects transformed rat retinal ganglion cells from hypoxia-induced apoptosis by preventing mitochondrial dysfunction. *Brain Res*. 2009 Jul 7;1279:131-8.
4. Xu WN, Liu WB, Liu ZP. Trichlorfon-induced apoptosis in hepatocyte primary cultures of *Carassius auratus gibelio*. *Chemosphere*. 2009 Nov;77(7):895-901.

5. Zhang Z, Wang S, Qiu H, Duan C, Ding K, Wang Z. Waltonitone induces human hepatocellular carcinoma cells apoptosis in vitro and in vivo. *Cancer Lett.* 2009 Dec 28;286(2):223-31.
6. Yang Q, Gong ZJ, Zhou Y, Yuan JQ, Cheng J, Tian L, Li S, Lin XD, Xu R, Zhu ZR, Mao C. Role of *Drosophila* alkaline ceramidase (Dacer) in *Drosophila* development and longevity. *Cell Mol Life Sci.* 2010 May;67(9):1477-90.
7. Zhong Y, Hu YJ, Yang Y, Peng W, Sun Y, Chen B, Huang X, Kong WJ. Contribution of common deletion to total deletion burden in mitochondrial DNA from inner ear of d-galactose-induced aging rats. *Mutat Res.* 2011 Jul 1;712(1-2):11-9.
8. Zhong Y, Hu YJ, Chen B, Peng W, Sun Y, Yang Y, Zhao XY, Fan GR, Huang X, Kong WJ. Mitochondrial transcription factor A overexpression and base excision repair deficiency in the inner ear of rats with D-galactose-induced aging. *FEBS J.* 2011 Jul;278(14):2500-10.
9. Chen B, Zhong Y, Peng W, Sun Y, Hu YJ, Yang Y, Kong WJ. Increased mitochondrial DNA damage and decreased base excision repair in the auditory cortex of D-galactose-induced aging rats. *Mol Biol Rep.* 2011 Aug;38(6):3635-42.
10. Liu Y, Zhou D, Zhang F, Tu Y, Xia Y, Wang H, Zhou B, Zhang Y, Wu J, Gao X, He Z, Zhai Q. Liver Ppt1 deficiency protects male mice from age-associated but not high-fat diet-induced hepatic steatosis. *J Lipid Res.* 2012 Mar;53(3):358-67.
11. Lu B, Xu Y, Xu L, Cong X, Yin L, Li H, Peng J. Mechanism investigation of dioscin against CCl₄-induced acute liver damage in mice. *Environ Toxicol Pharmacol.* 2012 Mar 30;31(2):127-135.
12. Feng Y, Zhang C, Luo Q, Wei X, Jiang B, Zhu H, Zhang L, Jiang L, Liu M, Xiao X. A novel WD-repeat protein, WDR26, inhibits apoptosis of cardiomyocytes induced by oxidative stress. *Free Radic Res.* 2012 Jun;46(6):777-84.
13. Zhao X, Cong X, Zheng L, Xu L, Yin L, Peng J. Dioscin, a natural steroid saponin, shows remarkable protective effect against acetaminophen-induced liver damage in vitro and in vivo. *Toxicol Lett.* 2012 Oct 2;214(1):69-80.
14. Zhong Y, Hu Y, Peng W, Sun Y, Yang Y, Zhao X, Huang X, Zhang H, Kong W. Age-related decline of the cytochrome c oxidase subunit expression in the auditory cortex of the mimetic aging rat model associated with the common deletion. *Hear Res.* 2012 Dec;294(1-2):40-8.
15. Wu L, Sun Y, Hu YJ, Yang Y, Yao LL, Zhou XX, Wang H, Zhang R, Huang X, Kong WJ. Increased p66Shc in the inner ear of D-galactose-induced aging mice with accumulation of mitochondrial DNA 3873-bp deletion: p66Shc and mtDNA damage in the inner ear during aging. *PLoS One.* 2012;7(11):e50483.
16. Tian X, Zhang L, Wang J, Dai J, Shen S, Yang L, Huang P. The protective effect of hyperbaric oxygen and Ginkgo biloba extract on Aβ25-35-induced oxidative stress and neuronal apoptosis in rats. *Behav Brain Res.* 2013 Apr 1;242:1-8.
17. Yin X, Yu B, Tang Z, He B, Ren J, Xiao X, Tang W. *Bifidobacterium infantis*-mediated HSV-TK/GCV suicide gene therapy induces both extrinsic and intrinsic apoptosis in a rat model of bladder cancer. *Cancer Gene Ther.* 2013 Feb;20(2):77-81.
18. Zhang S, Lu B, Han X, Xu L, Qi Y, Yin L, Xu Y, Zhao Y, Liu K, Peng J. Protection of the flavonoid fraction from *Rosa laevigata* Michx fruit against carbon tetrachloride-induced acute liver injury in mice. *Food Chem Toxicol.* 2013 May;55:60-9.
19. Zhang S, Zheng L, Dong D, Xu L, Yin L, Qi Y, Han X, Lin Y, Liu K, Peng J. Effects of flavonoids from *Rosa laevigata* Michx fruit against high-fat diet-induced non-alcoholic fatty liver disease in rats. *Food Chemistry.* 2013 Dec;141(3):2108-16.
20. Mo H, Chen Y, Huang L, Zhang H, Li J, Zhou W. Neuroprotective Effect of Tea Polyphenols on Oxyhemoglobin Induced Subarachnoid Hemorrhage in Mice. *Oxidative Medicine and Cellular Longevity.* Volume. 2013;2013:743938.
21. Huang L, Wan J, Chen Y, Wang Z, Hui L, Li Y, Xu D, Zhou W. Inhibitory effects of p38 inhibitor against mitochondrial dysfunction in the early brain injury after subarachnoid hemorrhage in mice. *Brain Res.* 2013 Jun 23;1517:133-40.
22. Guo Y, Wang W, Dong Y, Zhang Z, Zhou Y, Chen G. Carbon disulfide induces rat testicular injury via mitochondrial apoptotic pathway. *Chemosphere.* 2014 Aug;108:367-75.
23. Lan R, Zhang Y, Xiang J, Zhang W, Wang GH, Li WW, Xu LL, Cai DF, Xiao-Xu-Ming decoction preserves mitochondrial integrity and reduces apoptosis after focal cerebral ischemia and reperfusion via the mitochondrial p53 pathway. *J Ethnopharmacol.* 2014;151(1):307-16.
24. Li M, Sun M, Cao L, Gu JH, Ge J, Chen J, Han R, Qin YY, Zhou ZP, Ding Y, Qin ZH. A TIGAR-regulated metabolic pathway is critical for protection of brain ischemia. *J Neurosci.* 2014 May 28;34(22):7458-71.
25. Sun LY, Wang N, Ban T, Sun YH, Han Y, Sun LL, Yan Y, Kang XH, Chen S, Sun LH, Zhang R, Zhao YJ, Zhang H, Ai J, Yang BF. MicroRNA-23a mediates mitochondrial compromise in estrogen deficiency-induced concentric remodeling via targeting PGC-1α. *J Mol Cell Cardiol.* 2014 Oct;75:1-11.
26. Du LL, Chai DM, Zhao LN, Li XH, Zhang FC, Zhang HB, Liu LB, Wu K, Liu R, Wang JZ, Zhou XW. AMPK Activation Ameliorates Alzheimer's Disease-Like Pathology and Spatial Memory Impairment in a Streptozotocin-Induced Alzheimer's Disease Model in Rats. *J Alzheimers Dis.* 2015;43(3):775-84.
27. Zhang C, Zhang Z, Zhao Q, Wang X, Ji H, Zhang Y. (S)-ZJM-289 preconditioning induces a late phase protection against nervous injury induced by transient cerebral ischemia and oxygen-glucose deprivation. *Neurotox Res.* 2014 Jul;26(1):16-31.
28. Yu L, Wang S, Chen X, Yang H, Li X, Xu Y, Zhu X. Orientin alleviates cognitive deficits and oxidative stress in Aβ 1-42-induced mouse model of Alzheimer's disease. *Life Sci.* 2015 Jan 15;121:104-9.
29. Zhang Z, Huang Z, Dai H, Wei L, Sun S, Gao F. Therapeutic Efficacy of E-64-d, a Selective Calpain Inhibitor, in Experimental Acute Spinal Cord Injury. *Biomed Res Int.* 2015;2015:134242.
30. Li N, Tian Y, Wang C, Zhang P, You S. Protective effect of Lai Fu Cheng Qi decoction on severe acute pancreatitis-induced myocardial injury in a rat model. *Exp Ther Med.* 2015 Apr;9(4):1133-1140.
31. Xu Z, Zhang L, Li X, Jiang Z, Sun L, Zhao G, Zhou G, Zhang H, Shang J, Wang T. Mitochondrial fusion/fission process involved in the improvement of catalpol on high glucose-induced hepatic mitochondrial dysfunction. *Acta Biochim Biophys Sin(Shanghai).* 2015 Sep;47(9):730-40.
32. Xu X, He L, Zhang A, Li Q, Hu W, Chen H, Du J, Shen J. Toxoplasma gondii isolate with genotype Chinese 1 triggers trophoblast apoptosis through oxidative stress and mitochondrial dysfunction in mice. *Exp Parasitol.* 2015 Jul;154:51-61.
33. Su D, Zhang R, Zhang C, Huang F, Xiao J, Deng Y, Wei Z, Zhang Y, Chi J, Zhang M. Phenolic-rich lychee (*Litchi chinensis* Sonn.) pulp extracts offer hepatoprotection against restraint stress-induced liver injury in mice by modulating mitochondrial dysfunction. *Food Funct.* 2015 Nov 16. [Epub ahead of print].
34. Wang M, Sun GB, Zhang JY, Luo Y, Yu YL, Xu XD, Meng XB, Zhang MD, Lin WB, Sun XB. Elatoside C protects the heart from ischaemia/reperfusion injury through the modulation of oxidative stress and intracellular Ca²⁺ homeostasis. *Int J Cardiol.* 2015 Apr 15;185:167-76.
35. Liao XG, Li YL, Gao RF, Geng YQ, Chen XM, Liu XQ, Ding YB, Mu XY, Wang YX, He JL. Folate Deficiency Decreases Apoptosis of Endometrium Decidual Cells in Pregnant Mice via the Mitochondrial Pathway. *Nutrients.* 2015 Mar 13;7(3):1916-32.
36. Cao L, Chen J, Li M, Qin YY, Sun M, Sheng R, Han F, Wang G, Qin ZH. Endogenous level of TIGAR in brain is associated with vulnerability of neurons to ischemic injury. *Neurosci Bull.* 2015 Oct;31(5):527-40.
37. Du Z, Yang Q, Liu L, Li S, Zhao J, Hu J, Liu C, Qian D, Gao C. NADPH oxidase 2-dependent oxidative stress, mitochondrial damage and apoptosis in the ventral cochlear nucleus of d-galactose-induced aging rats. *J Ethnopharmacol.* 2015 Aug 22;172:124-32.
38. Wei W, Wang H, Wu Y, Ding K, Li T, Cong Z, Xu J, Zhou M, Huang L, Ding H, Wu H. Alpha lipoic acid inhibits neural apoptosis via a mitochondrial pathway in rats following traumatic brain injury. *Neurochem Int.* 2015 Aug;87:85-91.
39. Xing WM, Yuan TJ, Xu JD, Gu LL, Liang P, Lu H. Proteomic identification of mitochondrial targets involved in andrographolide sodium bisulfite-induced nephrotoxicity in a rat model. *Environ Toxicol Pharmacol.* 2015 Sep;40(2):592-9.
40. An FM, Chen S, Xu Z, Yin L, Wang Y, Liu AR, Yao WB, Gao XD. Glu cagon-like peptide-1 regulates mitochondrial biogenesis and tau phosphorylation against advanced glycation end product-induced neuronal insult: Studies in vivo and in vitro. *Neuroscience.* 2015 Aug 6;300:75-84.
41. Jin JL, Liou AK, Shi Y, Yin KL, Chen L, Li LL, Zhu XL, Qian L, Yang R, Chen J, Xu Y. CART treatment improves memory and synaptic structure in APP/PS1 mice. *Sci Rep.* 2015 May 11;5:10224.
42. Chen O, Ye Z, Cao Z, Manaenko A, Ning K, Zhai X, Zhang R, Zhang T, Chen X, Liu W, Sun X. Methane attenuates myocardial ischemia injury

- in rats through anti-oxidative, anti-apoptotic and anti-inflammatory actions. *Free Radic Biol Med.* 2016 Jan;90:1-11.
43. Xiong H, Du S, Ni J, Zhou J, Yao J. Mitochondria and nuclei dual-targeted heterogeneous hydroxyapatite nanoparticles for enhancing therapeutic efficacy of doxorubicin. *Biomaterials.* 2016 Jul;94:70-83.
 44. Li X, Wang H, Gao Y, Li L, Tang C, Wen G, Yang Y, Zhuang Z, Zhou M, Mao L, Fan Y. Quercetin induces mitochondrial biogenesis in experimental traumatic brain injury via the PGC-1α signaling pathway. *Am J Transl Res.* 2016 Aug 15;8(8):3558-66.
 45. Wang C, Duan X, Sun X, Liu Z, Sun P, Yang X, Sun H, Liu K, Meng Q. Protective effects of glycyrrhetic acid from edible botanical *glycyrrhiza* g labra against non-alcoholic steatohepatitis in mice. *Food Funct.* 2016 Sep 14;7(9):3716-23.
 46. Su D, Zhang R, Zhang C, Huang F, Xiao J, Deng Y, Wei Z, Zhang Y, Chi J, Zhang M. Phenolic-rich lychee (*Litchi chinensis* Sonn.) pulp extracts offer hepatoprotection against restraintstress-induced liver injury in mice by modulating mitochondrial dysfunction. *Food Funct.* 2016 Jan; 7(1): 508-15.
 47. Liu J, Wang Y, Zhuang Q, Chen M, Wang Y, Hou L, Han F. Protective effects of cyclosporine A and hypothermia on neuronal mitochondria in a rat asphyxialcardiac arrest model. *Am J Emerg Med.* 2016 Jun; 34(6): 1080-5.
 48. Li X, Wang H, Gao Y, Li L, Tang C, Wen G, Zhou Y, Zhou M, Mao L, Fan Y. Protective Effects of Quercetin on Mitochondrial Biogenesis in Experimental Traumatic Brain Injury via the Nrf2 Signaling Pathway. *PloS One.* 2016 Oct 25;11(10):e0164237.
 49. Jin W, Xu W, Chen J, Zhang X, Shi L, Ren C. Remote limb preconditioning protects against ischemia-induced neuronal death through amelioratingneuronal oxidative DNA damage and parthanatos. *J Neurol Sci.* 2016 Jul 15;366:8-17.
 50. Ma X, Xie Y, Chen Y, Han B, Li J, Qi S. Post-ischemia mdviv-1 treatment protects against ischemia/reperfusion-induced brain injury in a ratmodel. *Neurosci Lett.* 2016 Oct 6;632:23-32.
 51. Chen J, Wang Y, Hui C, Xi Y, Liu X, Qi F, Liu H, Wang Z, Niu S. Mechanisms of Heshouwuyin in regulating apoptosis of testicular cells in aging rats throughmitochondrial pathway. *BMC Complement Altern Med.* 2016 Sep 1;16:337.
 52. Zhai X, Ding Y, Wang Q, Zhang H, Li F. Rutin Acid Ameliorates Neural Apoptosis Induced by Traumatic Brain Injury via Mitochondrial Pathwaysin Mice. *Neuroimmunomodulation.* 2016;23(3):179-187.
 53. Chen O, Ye Z, Cao Z, Manaenko A, Ning K, Zhai X, Zhang R, Zhang T, Chen X, Liu W, Sun X. Methane attenuates myocardial ischemia injury in rats through anti-oxidative, anti-apoptotic and anti-inflammatory actions. *Free Radic Biol Med.* 2016 Jan;90:1-11.
 54. Huang HJ, Liu CW, Zhou X, Zhang CX, Bao YY. A mitochondrial membrane protein is a target for rice ragged stunt virus in its insect vector. *Virus Res.* 2017 Feb 2;229:48-56.
 55. Xu W, Jin W, Zhang X, Chen J, Ren C. Remote Limb Preconditioning Generates a Neuroprotective Effect by Modulating the Extrinsic Apoptotic Pathway and TRAIL-Receptors Expression. *Cell Mol Neurobiol.* 2017 Jan;37(1):169-182.
 56. Zhu Y, Di S, Hu W, Feng Y, Zhou Q, Gong B, Tang X, Liu J, Zhang W, Xi M, Jiang L, Guo C, Cao J, Fan C, Ma Z, Yang Y, Wen A. A new flavonoid glycoside (APG) isolated from *Clematis tangutica* attenuates myocardialischemia/reperfusion injury via activating PKCε signaling. *Biochim Biophys Acta.* 2017 Mar;1863(3):701-711.
 57. Wang X, Wang M, Xu J, Jia Z, Liu Z, Wang L, Song L. Soluble adenylyl cyclase mediates mitochondrial pathway of apoptosis and ATP metabolism in oyster *Crassostrea gigas* exposed to elevated CO₂. *Fish Shellfish Immun.* 2017 Jul;66:140-147.
 58. Xu J, Bian X, Liu Y, Hong L, Teng T, Sun Y, Xu Z. Adenosine A2 receptor activation ameliorates mitochondrial oxidative stress upon reperfusion throughthe posttranslational modification of NDUFV2 subunit of complex I in the heart. *Free Radical Bio Med.* 2017 May;106:208-218.
 59. Sun M, Wang R, Han Q. Inhibition of leukotriene B4 receptor 1 attenuates lipopolysaccharide-induced cardiaedysfunction: role of AMPK-regulated mitochondrial function. *Sci Rep-uk.* 2017 Mar 14;7:44352.
 60. Wang Z, Qiu Z, Gao C, Sun Y, Dong W, Zhang Y, Chen R, Qi Y, Li S, Guo Y, Piao Y, Li S, Piao F. 2,5-hexanedione downregulates nerve growth factor and induces neuron apoptosis inthe spinalcord of rats via inhibition ofthe PI3K/Akt signaling pathway. *Plos One.* 2017 Jun 27;12(6):e0179388.
 61. Ma ZN, Liu Z, Wang Z, Ren S, Tang S, Wang YP, Xiao SY, Chen C, Li W. Supplementation of American ginseng berry extract mitigated cisplatin-evoked nephrotoxicity by suppressing ROS-mediated activation of MAPK and NF- κ B signaling pathways. *Food Chem Toxicol.* 2017 Dec;110:62-73.
 62. Wang J, Wang Y, Shen L, Qian Y, Yang J, Wang F. Sulfated lentinan induced mitochondrial dysfunction leads to programmed cell death of tobaccoBY-2 cells. *Pestic Biochem Phys.* 2017 Apr;137:27-35.
 63. Yang RX, Lei J, Wang BD, Feng DY, Huang L, Li YQ, Li T, Zhu G, Li C, Lu FF, Nie TJ, Gao GD, Gao L. Pretreatment with Sodium Phenylbutyrate Alleviates Cerebral Ischemia/Reperfusion Injury by Upregulating DJ-1 Protein. *Front Neurol.* 2017 Jun 9;8:256.
 64. Zhao H, Liu YJ, Liu ZR, Tang DD, Chen XW, Chen YH, Zhou RN, Chen SQ, Niu HX. Role of mitochondrial dysfunction in renal fibrosis promoted by hypochlorite-modified albumin in a remnant kidney model and protective effects of antioxidant peptide SS-31. *Eur J Pharmacol.* 2017 Jun 5;804:57-67.
 65. Luan Y, Zhang F, Cheng Y, Liu J, Huang R, Yan M, Wang Y, He Z, Lai H, Wang H, Ying H, Guo F, Zhai Q. Hemin Improves Insulin Sensitivity and Lipid Metabolism in Cultured Hepatocytes and Mice Fed a High-Fat Diet. *Nutrients.* 2017 Jul 26;9(8). pii: E805.
 66. Peng X, Shang G, Wang W, Chen X, Lou Q, Zhai G, Li D, Du Z, Ye Y, Jin X, He J, Zhang Y, Yin Z. Fatty Acid Oxidation in Zebrafish Adipose Tissue Is Promoted by 1 α ,25(OH)2D3. *Cell Rep.* 2017 May 16;19(7):1444-1455.
 67. Huang HJ, Liu CW, Zhou X, Zhang CX, Bao YY. A mitochondrial membrane protein isa target for rice ragged stunt virus in its insect vector. *Virus Res.* 2017 Feb 2;229:48-56.
 68. Liu JD, Chen HJ, Wang DL, Wang H, Deng Q. Pim-1 Kinase Regulating Dynamics Related Protein 1 Mediates Sevoflurane Postconditioning-induced Cardioprotection. *Chinese Med J-peking.* 2017 Feb 5;130(3):309-317.
 69. Xu W, Jin W, Zhang X, Chen J, Ren C. Remote Limb Preconditioning Generates a Neuroprotective Effect by Modulating the ExtrinsicApoptotic Pathway and TRAIL-Receptors Expression. *Cell Mol Neurobiol.* 2017 Jan;37(1):169-182.
 70. Chen O, Cao Z, Li H, Ye Z, Zhang R, Zhang N, Huang J, Zhang T, Wang L, Han L, Liu W, Sun X. High-concentration hydrogen protects mouse heart against ischemia/reperfusion injury through activation of thePI3K/Akt1 pathway. *Sci Rep-uk.* 2017 Nov 1;7(1):14871.
 71. Peng K, Hu J, Xiao J, Dan G, Yang L, Ye F, Zou Z, Cao J, Sai YM. Mitochondrial ATP-sensitive potassium channel regulates mitochondrial dynamics to participate in neurodegeneration of Parkinson's disease. *Bba-Mol Basis Dis.* 2018 Apr;1864(4 Pt A):1086-1103.
 72. Li A, Liu Q, Li Q, Liu B, Yang Y, Zhang N. Berberine Reduces Pyruvate-driven Hepatic Glucose Production by Limiting Mitochondrial Import of Pyruvate through Mitochondrial Pyruvate Carrier 1. *Ebiomedicine.* 2018 Aug;34:243-255.
 73. Wang J, Dong B, Yu ZX, Yao CL. The impact of acute thermal stress on green mussel *Perna viridis*: Oxidative damage and responses. *Comp Biochem Physiol A Mol Integr Physiol.* 2018 Aug;222:7-15.
 74. Tang LX, Wang B, Wu ZK. Aerobic Exercise Training Alleviates Renal Injury by Interfering with Mitochondrial Function in Type-1 Diabetic Mice. *Med Sci Monitor.* 2018 Dec 14;24:9081-9089.
 75. Zhang J, Yu J, Chen Y, Liu L, Xu M, Sun L, Luo H, Wang Y, Meng G. Exogenous Hydrogen Sulfide Supplement Attenuates Isoproterenol-Induced Myocardial Hypertrophy in a Sirtuin 3-Dependent Manner. *Oxid Med Cell Longev.* 2018 Dec 17;2018:9396089.
 76. Zhong Y, Jin C, Wang X, Li X, Han J, Xue W, Wu P, Peng X, Xia X. Protective effects of apigenin against 3-MCPD-induced renal injury in rat. *Chem-Biol Interact.* 2018 Dec 25;296:9-17.
 77. Wang X, Lan YL, Xing JS, Lan XQ, Wang LT, Zhang B. Alantolactone plays neuroprotective roles in traumatic brain injury in rats via anti-inflammatory, anti-oxidative and anti-apoptosis pathways. *Am J Transl Res.* 2018 Feb 15;10(2):368-380. eCollection 2018.
 78. Fu C, Wang Q, Zhai X, Gao J. Sinomenine reduces neuronal cell apoptosis in mice after traumatic brain injury via its effect on mitochondrial pathway. *Drug Des Dev Ther.* 2018 Jan 5;12:77-84.
 79. Kang J, Jia Z, Ping Y, Liu Z, Yan X, Xing G, Yan W. Testosterone alleviates mitochondrial ROS accumulation and mitochondria-mediated apoptosis in the gastric mucosa of orchectomized rats. *Arch Biochem Biophys.* 2018 Jul 1;649:53-59.

80. Zhang X,Zhang Y,Tang S,Yu L,Zhao Y,Ren Q,Huang X,Xu W,Huang M,Peng J.Pien-Tze-Huang protects cerebral ischemic injury by inhibiting neuronal apoptosis in acute ischemic stroke rats. *J Ethnopharmacol*. 2018 Jun 12;219:117-125.
81. Zhou J,Wang H,Shen R,Fang J,Yang Y,Dai W,Zhu Y,Zhou MMitochondrial-targeted antioxidant MitoQ provides neuroprotection and reduces neuronal apoptosis in experimental traumatic brain injury possibly via the Nrf2-ARE pathway. *Am J Transl Res*. 2018 Jun 15;10(6):1887-1899. eCollection 2018.
82. Wei Q,Luo Q,Liu H,Chen L,Cui H,Fang J,Zuo Z,Deng J,Li Y,Wang X,Zhao L.The mitochondrial pathway is involved in sodium fluoride (NaF)-induced renal apoptosis in mice. *Toxicol Res-uk*. 2018 Jun 23;7(5):792-808.
83. Zhang Z,Liu J,Fan C,Mao L,Xie R,Wang S,Yang M,Yuan H,Yang X,Sun J,Wang J,Kong J,Huang S,Sun B.The GluN1/GluN2B NMDA receptor and metabotropic glutamate receptor 1 negative allosteric modulator has enhanced neuroprotection in a rat subarachnoid hemorrhage model. *Exp Neurol*. 2018 Mar;301(Pt A):13-25.
84. Wu H,Zhang Q,Gao J,Sun C,Wang J,Xia W,Cao Y,Hao Y,Wu L.Modulation of sphingosine 1-phosphate (S1P) attenuates spatial learning and memory impairments in the valproic acid rat model of autism. *Psychopharmacology*. 2018 Mar;235(3):873-886.
85. Shi Z,Li C,Yin Y,Yang Z,Xue H,Mu N,Wang Y,Liu M,Haerobic Interval Training Regulated SIRT3 Attenuates High-Fat-Diet-Associated Cognitive Dysfunction. *Biomed Res Int*. 2018 Mar 22;2018:2708491.
86. Zhang G,Sheng M,Wang J,Teng T,Sun Y,Yang Q,Xu ZZinc improves mitochondrial respiratory function and prevents mitochondrial ROS generation at reperfusion by phosphorylating STAT3 at Ser727. *J Mol Cell Cardiol*. 2018 May;118:169-182.
87. Liu ZR,Chen SQ,Zou YW,Wu XY,Li HY,Wang XQ,Shi Y,Niu HX. Hypochlorite modified albumins promote cell death in the tubule interstitium in rats via mitochondrial damage in obstructive nephropathy and the protective effects of antioxidant peptides. *Free Radical Res*. 2018 May;52(5):616-628.
88. Fang S,Zhuo Z,Yu X,Wang H,Feng J.Oral administration of liquid iron preparation containing excess iron induces intestine and liver injury, impairs intestinal barrier function and alters the gut microbiota in rats. *J Trace Elem Med Bio*. 2018 May;47:12-20.
89. Jia P,Liu C,Wu N,Jia D,Sun Y.Agomelatine protects against myocardial ischemia reperfusion injury by inhibiting mitochondrial permeability transition pore opening. *Am J Transl Res*. 2018 May 15;10(5):1310-1323. eCollection 2018.
90. Liu SM,Li XZ,Zhang SN,Yang ZM,Wang KX,Lu F,Wang CZ,Yuan CS. Acanthopanax senticosus Protects Structure and Function of Mesencephalic Mitochondria in A Mouse Model of Parkinson's Disease. *Chin J Integr Med*. 2018 Nov;24(11):835-843.
91. Xu J,Pan H,Xie X,Zhang J,Wang Y,Yang G. Inhibiting Succinate Dehydrogenase by Dimethyl Malonate Alleviates Brain Damage in a Rat Model of Cardiac Arrest. *Neuroscience*. 2018 Nov 21;393:24-32.
92. Shao X,Chen Q,Dou X,Chen L,Wu J,Zhang W,Shao H,Ling P,Liu F,Wang F. Lower range of molecular weight of xanthan gum inhibits cartilage matrix destruction via intrinsic bax-mitochondria cytochrome c-caspase pathway. *Carbohydr Polym*. 2018 Oct 15;198:354-363..
93. Zhou W,Chen X,Zhao G,Xu D,Jiang Z,Zhang L,Wang T. Psoralen Induced Liver Injury by Attenuating Liver Regenerative Capability. *Front Pharmacol*. 2018 Oct 22;9:1179.
94. Xu C,Sun Y,Cai X,You T,Zhao H,Li Y,Zhao H. Medial Habenula-Interpeduncular Nucleus Circuit Contributes to Anhedonia-Like Behavior in a Rat Model of Depression. *Front Behav Neurosci*. 2018 Oct 9;12:238.
95. Che Y,Wang ZP,Yuan Y,Zhang N,Jin YG,Wan CX,Tang QZ. Role of autophagy in a model of obesity: A long-term high fat diet induces cardiac dysfunction. *Mol Med Rep*. 2018 Sep;18(3):3251-3261.
96. Sun Y,Zhang T,Zhang Y,Li J,Jin L,Sun Y,Shi N,Liu K,Sun X. Ischemic Postconditioning Alleviates Cerebral Ischemia-Reperfusion Injury Through Activating Autophagy During Early Reperfusion in Rats. *Neurochem Res*. 2018 Sep;43(9):1826-1840.
97. Xiong H,Du S,Zhang P,Jiang Z,Zhou J,Yao J. Primary tumor and pre-metastatic niches co-targeting "peptides-lego" hybrid hydroxyapatite nanoparticles for metastatic breast cancer treatment. *Biomater Sci-uk*. 2018 Sep 25;6(10):2591-2604.
98. Fang J,Zhu Y,Wang H,Cao B,Fei M,Niu W,Zhou Y,Wang X,Li X,Zhou M. Baicalin Protects Mice Brain From Apoptosis in Traumatic Brain Injury Model Through Activation of Autophagy. *Front Neurosci-Switz*. 2019 Jan 9;12:1006.
99. Ye JS,Chen L,Lu YY,Lei SQ,Peng M,Xia ZYSIRT3 activator honokiol ameliorates surgery/anesthesia-induced cognitive decline in mice through anti-oxidative stress and anti-inflammatory in hippocampus. *CNS Neurosci Ther*. 2019 Mar;25(3):355-366.
100. Su X,Wu Z,Mai F,Fan Z,Du S,Qian H,Zhu J.'Governor vessel-unblocking and mind-regulating' acupuncture therapy ameliorates cognitive dysfunction in a rat model of middle cerebral artery occlusion. *Int J Mol Med*. 2019 Jan;43(1):221-232.
101. Wen D,Cui C,Duan W,Wang W,Wang Y,Liu Y,Li Z,Li C.The role of insulin-like growth factor 1 in ALS cell and mouse models: A mitochondrial protector. *Brain Res Bull*. 2019 Jan;144:1-13.
102. Miao J,Huang Z,Liu S,Li X,Jia P,Guo Y,Wu N,Jia D.Hydroxytyrosol protects against myocardial ischemia reperfusion injury by inhibiting mitochondrial permeability transition pore opening. *Exp Ther Med*. 2019 Jan;17(1):671-678.
103. Qin J,Wang P,Li Y,Yao L,Liu Y,Yu T,Lin J,Fang X,Huang Z.Activation of Sigma-1 Receptor by Cutame sine Attenuates Neuronal Apoptosis by Inhibiting Endoplasmic Reticulum Stress and Mitochondrial Dysfunction in a Rat Model of Asphyxia Cardiac Arrest. *Shock*. 2019 Jan;51(1):105-113.
104. Hui F,Qin X,Zhang Q,Li R,Liu M,Ren T,Zhao M,Zhao Q.Alpinia oxyphylla oil induces apoptosis of hepatocellular carcinoma cells via PI3K/Akt pathway in vitro and in vivo. *Biomed Pharmacother*. 2019 Jan;109:2365-2374.
105. Song J,Lin C,Yang X,Xie Y,Hu P,Li H,Zhu W,Hu H.Mitochondrial targeting nanodrugs self-assembled from 9-O-octadecyl substituted berberine derivative for cancer treatment by inducing mitochondrial apoptosis pathways. *J Control Release*. 2019 Jan 28;294:27-42.
106. Mao J,Li Y,Li S,Li J,Tian Y,Feng S,Liu X,Bian Q,Li J,Hu Y,Zhang L,JI H.Bufei Jianpi Granules Reduce Quadriceps Muscular Cell Apoptosis by Improving Mitochondrial Function in Rats with Chronic Obstructive Pulmonary Disease.EVID-BASED COMPL ALT. 2019 Jan 6;2019:1216305.
107. Li K,Li W,Yin H,Cheong YK,Ren G,Yang Z.Pretreatment-Etidronate Alleviates CoCl₂ Induced-SH-SY5Y Cell Apoptosis via Decreased HIF-1 α and TRPC5 Channel Proteins. *Neurochem Res*. 2019 Feb;44(2):428-440.
108. Wei L,Zhao W,Hu Y,Wang X,Liu X,Zhang P,Han F.Exploration of the optimal dose of HOE-642 for the protection of neuronal mitochondrial function after cardiac arrest in rats. *Biomed Pharmacother*. 2019 Feb;110:818-824.
109. Zuo E,Zhang C,Mao J,Gao C,Hu S,Shi X,Piao F.2,5-Hexanedione mediates neuronal apoptosis through suppression of NGF via PI3K/Akt signaling in the rat sciatic nerve. *BIOSCIENCE REP*. 2019 Feb 12;39(2). pii: BSR20181122.
110. Ye JS,Chen L,Lu YY,Lei SQ,Peng M,Xia ZY.Honokiol-Mediated Mitophagy Ameliorates Postoperative Cognitive Impairment Induced by Surgery/Sevoflurane via Inhibiting the Activation of NLRP3 Inflammasome in the Hippocampus. *Oxid Med Cell Longev*. 2019 Feb 24;2019:8639618.
111. Diao J,Wei J,Yan R,Fan G,Lin L,Chen M..Effects of resveratrol on regulation on UCP2 and cardiac function in diabetic rats. *J Physiol Biochem*. 2019 Feb;75(1):39-51.
112. Ye JS,Chen L,Lu YY,Lei SQ,Peng M,Xia ZY.SIRT3 activator honokiol ameliorates surgery/anesthesia-induced cognitive decline in mice through anti-oxidative stress and anti-inflammatory in hippocampus. *CNS Neurosci Ther*. 2019 Mar;25(3):355-366.
113. Wang M,Bi Y,Zeng S,Liu Y,Shao M,Liu K,Deng Y,Wen G,Sun X,Zeng P,Jing L,Lv Z.Modified Xiaoyaosan ameliorates depressive-like behaviors by triggering autophagosome formation to alleviate neuronal apoptosis. *Biomed Pharmacother*. 2019 Mar;111:1057-1065.
114. Du ZD,Yu S,Qi Y,Qu TF,He L,Wei W,Liu K,Gong SS.NADPH oxidase inhibitor apocynin decreases mitochondrial dysfunction and apoptosis in the ventral cochlear nucleus of D-galactose-induced aging model in rats. *Neurochem Int*. 2019 Mar;124:31-40.
115. Xu X,Wang D,Zheng C,Gao B,Fan J,Cheng P,Liu B,Yang L,Luo Z.Progerin accumulation in nucleus pulposus cells impairs mitochondrial function and induces intervertebral disc degeneration and therapeutic effects of sulforaphane. *Theranostics*. 2019 Apr 12;9(8):2252-2267.
116. Geng J,Wei M,Yuan X,Liu Z,Wang X,Zhang D,Luo L,Wu J,Guo W,Qin ZH.TIGAR regulates mitochondrial functions through SIRT1-PGC1 α

- pathway and translocation of TIGAR into mitochondria in skeletal muscle.FASEB J. 2019 May;33(5):6082-6098
- 117.Xie M,Jiang L,Dun Y,Zhang W,Liu S.Trimetazidine combined with exercise improves exercise capacity and anti-fatal stress ability through enhancing mitochondrial quality control.Life Sci. 2019 May 1;224:157-168
- 118.Shi N,He J,Guo Q,Liu T,Han J.Liraglutide protects against diabetes mellitus complicated with focal cerebral ischemic injury by activating mitochondrial ATP-sensitive potassium channels.Neuroreport. 2019 May 1;30(7):479-484
- 119.Xian W,Li T,Li L,Hu L,Cao J.Maresin 1 attenuates the inflammatory response and mitochondrial damage in mice with cerebral ischemia/reperfusion in a SIRT1-dependent manner.Brain Res. 2019 May 15;1711:83-90
- 120.Hu L,Ding M,Tang D,Gao E,Li C,Wang K,Qi B,Qiu J,Zhao H,Chang P,Fu F,Li Y.Targeting mitochondrial dynamics by regulating Mfn2 for therapeutic intervention in diabetic cardiomyopathy.Theranostics. 2019 May 31;9(13):3687-3706
- 121.Hasnat M,Yuan Z,Naveed M,Khan A,Raza F,Xu D,Ullah A,Sun L,Zhang L,Jiang Z.Drp1-associated mitochondrial dysfunction and mitochondrial autophagy: a novel mechanism in triptolid-induced hepatotoxicity.Cell Biol Toxicol. 2019 Jun;35(3):267-280
- 122.Xu M,Zhang Y,Wang M,Zhang H,Chen Y,Adcock IM,Chung KF,Mo J,Zhang Y,Li F,TRPV1 and TRPA1 in Lung Inflammation and Airway Hyperresponsiveness Induced by Fine Particulate Matter.Oxid Med Cell Longev. 2019 Jun 2;2019:7450151
- 123.Yang HX,Wang P,Wang NN,Li SD,Yang MH.Tongxinluo Ameliorates Myocardial Ischemia-Reperfusion Injury Mainly via Activating Parkin-Mediated Mitophagy and Downregulating Ubiquitin-Proteasome System.Clin J Integr Med. 2019 Jun 21
- 124.Wang S,Zhao Z,Fan Y,Zhang M,Feng X,Lin J,Hu J,Cheng Z,Sun C,Liu T,Xiong Z,Yang Z,Wang H,Sun D.Mst1 inhibits Sirt3 expression and contributes to diabetic cardiomyopathy through inhibiting Parkin-dependent mitophagy.BBA-MOL BASIS DIS. 2019 Jul 1;1865(7):1905-1914
- 125.Wang XY,Yu J,Zhang FY,Liu KJ,Xiang B.Phenylephrine Alleviates 131I Radiation Damage in Submandibular Gland Through Maintaining Mitochondrial Homeostasis.INT J RADIAT ONCOL. 2019 Jul 1;104(3):644-655
- 126.Li X,Jia P,Huang Z,Liu S,Miao J,Guo Y,Wu N,Jia D.Lycopene protects against myocardial ischemia-reperfusion injury by inhibiting mitochondrial permeability transition pore opening.DRUG DES DEV THER. 2019 Jul 11;13:2331-2342
- 127.Wang Y,Ni J,Gao C,Xie L,Zhai L,Cui G,Yin X.Mitochondrial transplantation attenuates lipopolysaccharide-induced depression-like behaviors.PROG NEURO-PSYCHOPH. 2019 Jul 13;93:240-249
- 128.Zhang H,Hu Y,Wang H,Tian L,Li W,Han L,Xu H,Ma J,Mo K,Xiao B,Chen L.Cytoplasmic upregulation of Cyto c and AIF serve as biomarkers of mechanical asphyxia death.Am J Transl Res. 2019 Jul 15;11(7):4568-4583. eCollection 2019
- 129.Du ZD,Wei W,Yu S,Song QL,Liu K,Gong SS.NADPH Oxidase 2-Mediated Insult in the Auditory Cortex of Zucker Diabetic Fatty Rats.Neural Plast. 2019 Jul 30;2019:3591605
- 130.Feibo Xu,Peiyan Wang,Qiucheng Yao,Bing Shao,Hongyan Yu,Kaiyuan Yu,Yanfei Li.Lycopene Alleviates AFB 1-induced Immunosuppression by Inhibiting Oxidative Stress and Apoptosis in the Spleen of Mice Food Funct. 2019 Jul 17;10(7):3868-3879.;doi: 10.1039/c8fo02300j
- 131.Ren Y,Chen X,Li P,Zhang H,Su C,Zeng Z,Wu Y,Xie X,Wang Q,Han J,Guo S,Liu B,Wang W.Si-Miao-Yong-An decoction ameliorates cardiac function through restoring the equilibrium of SOD and NOX2 in heart failure mice.Pharmacol Res. 2019 Aug;146:104318
- 132.Xu H,Li Y,Liu R,Wu L,Zhang C,Ding N,Ma A,Zhang J,Xie X.Protective effects of ghrelin on brain mitochondria after cardiac arrest and resuscitation.Neuropeptides. 2019 Aug;76:101936
- 133.Lin Q,Li S,Jiang N,Shao X,Zhang M,Jin H,Zhang Z,Shen J,Zhou Y,Zhou W,Gu L,Lu R,Ni Z.PINK1-parkin pathway of mitophagy protects against contrast-induced acute kidney injury via decreasing mitochondrial ROS and NLRP3 inflammasome activation.Redox Biol. 2019 Sep;26:101254
- 134.Xiao X,Chen Q,Zhu X,Wang Y.ABAD/17 β -HSD10 reduction contributes to the protective mechanism of huperzine a on the cerebral mitochondrial function in APP/PS1 mice.Neurobiol Aging. 2019 Sep;81:77-87
- 135.Yu J,Chen Y,Xu M,Sun L,Luo H,Bao X,Meng G,Zhang W.Ca²⁺/Calmodulin-Dependent Protein Kinase II Regulation by Inhibitor 1 of Protein Phosphatase 1 Protects Against Myocardial Ischemia-Reperfusion Injury.J Cardiovasc Pharmacol Ther. 2019 Sep;24(5):460-473.
- 136.He J,Huang Y,Du G,Wang Z,Xiang Y,Wang Q.Lasting spatial learning and memory deficits following chronic cerebral hypoperfusion are associated with hippocampal mitochondrial aging in rats.Neuroscience. 2019 Sep 1;415:215-229
- 137.Li R,Chen J.Salidroside Protects Dopaminergic Neurons by Enhancing PINK1/Parkin-Mediated Mitophagy.Oxid Med Cell Longev. 2019 Sep 10;2019:9341018
- 138.Yu Y,Xu C,Zhen L,Yang S,Zhou J,Yao J.Bio-inspired drug-dominated supramolecular nanocomplex based on low molecular weight heparin for progressive tumor therapy.CARBOHYD POLYM. 2019 Sep 15;220:30-42
- 139.Sun J,Zhou Q,Hu X.Integrating multi-omics and regular analyses identifies the molecular responses of zebrafish brains to graphene oxide: Perspectives in environmental criteria.ECOTOX ENVIRON SAFE. 2019 Sep 30;180:269-279
- 140.Li X,Luo S,Zhang J,Yuan Y,Jiang W,Zhu H,Ding X,Zhan L,Wu H,Xie Y,Song R,Pan Z,Lu Y.lncRNA H19 Alleviated Myocardial I/Ri via Suppressing miR-877-3p/Bcl-2-Mediated Mitochondrial Apoptosis.MOL THER-NUCL ACIDS. 2019 Sep 6;17:297-309
- 141.Wang M,Yang L,Yang J,Zhou Y,Wang C.Magnesium lithospermate B attenuates renal injury in 5/6 renal ablation/infarction rats by mitochondrial pathway of apoptosis.Biomed Pharmacother. 2019 Oct;118:109316
- 142.Liu D,Wu L,Wu Y,Wei X,Wang W,Zhang S,Yi M,Li J,Liu H,Ma X.Heat shock factor 1-mediated transcription activation of Omi/HtrA2 induces myocardial mitochondrial apoptosis in the aging heart.AGING-US. 2019 Oct 18;11(20):8982-8997
- 143.Yang X,Lu GP,Cai XD,Lu ZJ,Kissoon N.Alterations of complex IV in the tissues of a septic mouse model.Mitochondrion. 2019 Nov;49:89-96
- 144.Zhang Z,Wang J,Zhu Y,Zhang H,Wang H.Astragaloside IV alleviates myocardial damage induced by type 2 diabetes via improving energy metabolism.Mol Med Rep. 2019 Nov;20(5):4612-4622
- 145.Shan S,Shen Z,Zhang C,Kou R,Xie K,Song F.Mitophagy protects against acetaminophen-induced acute liver injury in mice through inhibiting NLRP3 inflammasome activation.Biochem Pharmacol. 2019 Nov;169:113643
- 146.Wu Y,Wang Y,Yin D,Wu W,Sun X,Zhang Y,Guo X,Chen J,Yuan J.Effect of supplementation of nicotinamide and sodium butyrate on the growth performance, liver mitochondrial function and gut microbiota of broilers at high stocking density.Food Funct. 2019 Nov 1;10(11):7081-7090
- 147.Zhao W,Xu Z,Cao J,Fu Q,Wu Y,Zhang X,Long Y,Zhang X,Yang Y,Li Y,Mi W.Elamipretide (SS-31) improves mitochondrial dysfunction, synaptic and memory impairment induced by lipopolysaccharide in mice.J NEUROINFLAMM. 2019 Nov 20;16(1):230
- 148.XR,Rong JB,Lu HS,Daugherty A,Shi P,Ke CL,Zhang ZC,Xu YC,Wang JA.Angiotensinogen in hepatocytes contributes to Western diet-induced liver steatosis.J Lipid Res. 2019 Dec;60(12):1983-1995
- 149.Liu K,Jing MJ,Liu C,Yan DY,Ma Z,Wang C,Deng Y,Liu W,Xu B.Effect of trehalose on manganese-induced mitochondrial dysfunction and neuronal cell damage in mice.BASIC CLIN PHARMACOL. 2019 Dec;125(6):536-547
- 150.Wang J,Dong ZH,Gui MT,Yao L,Li JH,Zhou XJ,Fu DY,Huo Xue Qian Yang Qu Tan Recipe attenuates left ventricular hypertrophy in obese hypertensive rats by improving mitochondrial function through SIRT1/PGC-1 α deacetylation pathway.BIOSCIENCE REP. 2019 Dec 20;39(12). pii: BSR20192909
- 151.He X,Gao J,Hou H,Qi Z,Chen H,Zhang XX.Inhibition of Mitochondrial Fatty Acid Oxidation Contributes to Development of Nonalcoholic Fatty Liver Disease Induced by Environmental Cadmium Exposure.Environ Sci Technol. 2019 Dec 3;53(23):13992-14000
- 152.Li S,Lin Q,Shao X,Zhu X,Wu J,Wu B,Zhang M,Zhou W,Zhou Y,Jin H,Zhang Z,Qi C,Shen J,Mou S,Gu L,Ni Z.Drp1-regulated PARK2-dependent mitophagy protects against renal fibrosis in unilateral ureteral obstruction.FREE RADICAL BIO MED. 2019 Dec 9;pii: S0891-5849(19)31684-3
- 153.Wanyue Huang,Zheng Cao,Qiucheng Yao,Qiang Ji,Jian Zhang,Yanfei Li.Mitochondrial damage are involved in Aflatoxin B 1-induced testicular damage and spermatogenesis disorder in mice Sci Total Environ. 2020 Jan 20;701:135077.;doi: 10.1016/j.scitotenv.2019.135077

- 154.Wen Xu,Yue Fu,Longyuan Jiang,Zhengfei Yang,Yue Wang,Wanchun Tang,Xiangshao Fang.Cardiopulmonary resuscitation ameliorates myocardial mitochondrial dysfunction in a cardiac arrest rat model Am J Emerg Med. 2020 Jan;38(1):65-72.;doi: 10.1016/j.ajem.2019.04.024
- 155.Yan Cao,Xiaotong Han,Hongwei Pan,Yu Jiang,Xiang Peng,Weiwei Xiao,Jingjing Rong,Fang Chen,Jin He,Lianhong Zou,Yi Tang,Yanfang Pei,Jiao Zheng,Jia Wang,Jie Zhong,Xiuqing Hong,Zhengyu Liu,Zhaofen Zheng.Emerging protective roles of shengmai injection in septic cardiomyopathy in mice by inducing myocardial mitochondrial autophagy via caspase-3/Beclin-1 axis Inflamm Res. 2020 Jan;69(1):41-50.;doi: 10.1007/s00011-019-01292-2
- 156.Long Yang,Peng Xie,Jianjiang Wu,Jin Yu,Xin Li,Haiping Ma,Tian Yu,Haiying Wang,Jianrong Ye,Jiang Wang,Hong Zheng.Deferoxamine Treatment Combined With Sevoflurane Postconditioning Attenuates Myocardial Ischemia-Reperfusion Injury by Restoring HIF-1/BNIP3-Mediated Mitochondrial Autophagy in GK Rats Front Pharmacol. 2020 Feb 19;11:6.;doi: 10.3389/fphar.2020.00006
- 157.Zhenyu Zhu,Huihui Li,Wanli Chen,Yameng Cui,Anan Huang,Xin Qi.Perindopril Improves Cardiac Function by Enhancing the Expression of SIRT3 and PGC-1 α in a Rat Model of Isoproterenol-Induced Cardiomyopathy Front Pharmacol. 2020 Feb 21;11:94.;doi: 10.3389/fphar.2020.00094
- 158.Changchun Cao,Jun Zhou,Xiaoli Wu,Yuanyuan Qian,Yali Hong,Junyu Mu,Lai Jin,Chao Zhu,Shengnan Li.Activation of CRHR1 contributes to cerebral endothelial barrier impairment via cPLA2 phosphorylation in experimental ischemic stroke Cell Signal. 2020 Feb;66:109467.;doi: 10.1016/j.cellsig.2019.109467
- 159.Tingjuan Ni,Na Lin,Xingxiao Huang,Wenqiang Lu,Zhenzhu Sun,Jie Zhang,Hui Lin,Jufang Chi,Handyuan Guo.Icarin Ameliorates Diabetic Cardiomyopathy Through Apelin/Sirt3 Signalling to Improve Mitochondrial Dysfunction Front Pharmacol. 2020 Mar 19;11:256.;doi: 10.3389/fphar.2020.00256
- 160.Wenfang Zhang,Yinchuan Jin,Dong Wang,Jingjing Cui.Neuroprotective effects of leptin on cerebral ischemia through JAK2/STAT3/PGC-1-mediated mitochondrial function modulation Brain Res Bull. 2020 Mar;156:118-130.;doi: 10.1016/j.brainresbull.2020.01.002
- 161.Weiwei Li,Jiancheng Yang,Qiufeng Lyu,Gaofeng Wu,Shumei Lin,Qunhui Yang,Jiammin Hu.Taurine prevents cardiomyocyte apoptosis by inhibiting the calpain-1/cytochrome c pathway during RVH in broilers Amino Acids. 2020 Mar;52(3):453-463.;doi: 10.1007/s00726-020-0284-5
- 162.Xue Bai,Tian-Yang Tan,Yun-Xin Li,Yue Li,Ya-Fei Chen,Ru Ma,Shu-Yan Wang,Qiang Li,Zhen-Quan Liu.The protective effect of cordyceps sinensis extract on cerebral ischemic injury via modulating the mitochondrial respiratory chain and inhibiting the mitochondrial apoptotic pathway Biomed Pharmacother. 2020 Apr;124:109834.;doi: 10.1016/j.biopha.2020.109834
- 163.Wei Xue,Xin Wang,Hong Tang,Fanfan Sun,Huaqing Zhu,Dake Huang,Liuyi Dong.Vitexin attenuates myocardial ischemia/reperfusion injury in rats by regulating mitochondrial dysfunction induced by mitochondrial dynamics imbalance Biomed Pharmacother. 2020 Apr;124:109849.;doi: 10.1016/j.biopha.2020.109849
- 164.Hao Zhang,Yanan Chen,Yueping Chen,Peilu Jia,Shuli Ji,Jianxiong Xu,Yue Li,Tian Wang.Comparison of the effects of resveratrol and its derivative pterostilbene on hepatic oxidative stress and mitochondrial dysfunction in piglets challenged with diquat Food Funct. 2020 May 1;11(5):4202-4215.;doi: 10.1039/d0fo00732c
- 165.Na Zhao,Qing-Wei Yan,Jie Xia,Xian-Liang Zhang,Bai-Xia Li,Ling-Yu Yin,Bo Xu.Treadmill Exercise Attenuates A β -Induced Mitochondrial Dysfunction and Enhances Mitophagy Activity in APP/PS1 Transgenic Mice Neurochem Res. 2020 May;45(5):1202-1214.;doi: 10.1007/s11064-020-03003-4
- 166.Bingchao Qi,Linjie He,Ya Zhao,Ling Zhang,Yuanfang He,Jun Li,Congye Li,Bo Zhang,Qichao Huang,Jinliang Xing,Fei Li,Yan Li,Lele Ji.Akap1 deficiency exacerbates diabetic cardiomyopathy in mice by NDUFS1-mediated mitochondrial dysfunction and apoptosis Diabetologia. 2020 May;63(5):1072-1087.;doi: 10.1007/s00125-020-05103-w
- 167.Wen Wang,Fengrui Zhang,Xiaoyu Li,Juan Luo,Yang Sun,Jing Wu,Maojuan Li,Yunling Wen,Hao Liang,Kunhua Wang,Junkun Niu,Yinglei Miao.Heat shock transcription factor 2 inhibits intestinal epithelial cell apoptosis through the mitochondrial pathway in ulcerative colitis BIOCHEM BIOPH RES CO. 2020 Jun 18;527(1):173-179.;doi: 10.1016/j.bbrc.2020.04.103
- 168.Wen-Yuan Ling,Ying Cui,Jun-Ling Gao,Xiao-Hua Jiang,Kai-Jie Wang,Yan-Xia Tian,Hua-Xin Sheng,Jian-Zhong Cui.Long-term chemogenetic activation of M1 glutamatergic neurons attenuates the behavioral and cognitive deficits caused by intracerebral hemorrhage BIOCHEM BIOPH RES CO. 2020 Jun 18;527(1):22-28.;doi: 10.1016/j.bbrc.2020.04.083
- 169.Xiaoling Zhang,Wenmin Huang,Qianhang Shao,Yuan Yang,Zhengxin Xu,Jing Chen,Xiaoyan Zhang,Xiaoqun Ge.Drp1, a potential therapeutic target for Parkinson's disease, is involved in olfactory bulb pathological alteration in the Rotenone-induced rat model Toxicol Lett. 2020 Jun 1;325:1-13.;doi: 10.1016/j.toxlet.2020.02.009
- 170.Xiao-Hui Zhang,Jia-Xin Wu,Jun-Zhou Sha,Bo Yang,Jia-Rui Sun,En-Dong Bao.Heat shock protein 90 relieves heat stress damage of myocardial cells by regulating Akt and PKM2 signaling in vivo Int J Mol Med. 2020 Jun;45(6):1888-1908.;doi: 10.3892/ijmm.2020.4560
- 171.Shuhuan Liu,Jing Yang,Xudong Peng,Jingjing Li,Cunjing Zhu.The Natural Product Fucoidan Inhibits Proliferation and Induces Apoptosis of Human Ovarian Cancer Cells: Focus on the PI3K/Akt Signaling Pathway Cancer Manag Res. 2020 Jul 23;12:6195-6207.;doi: 10.2147/CMAR.S254784
- 172.Yan Xiao,Wanying Chen,Zehao Zhong,Liang Ding,Hua Bai,Hao Chen,Hongru Zhang,Yihuang Gu,Shengfeng Lu.Electroacupuncture preconditioning attenuates myocardial ischemia-reperfusion injury by inhibiting mitophagy mediated by the mTORC1-ULK1-FUND1 pathway Biomed Pharmacother. 2020 Jul;127:110148.;doi: 10.1016/j.biopha.2020.110148
- 173.Xinwei Li,Zhen Shi,Yiwei Zhu,Taiyu Shen,Heyuan Wang,Guangzhou Shui,Juan J Loor,Zhiyuan Fang,Meng Chen,Xinghui Wang,Zhicheng Peng,Yuxiang Song,Zhe Wang,Xiliang Du,Guowen Liu.Cyanidin-3-O-glucoside improves non-alcoholic fatty liver disease by promoting PINK1-mediated mitophagy in mice BRIT J PHARMACOL. 2020 Aug;177(15):3591-3607.;doi: 10.1111/bph.15083
- 174.Pin Zhao,Li Zhang,Longfei Gao,Qian Ding,Qian Yang,Jianke Kuai.Ulinastatin attenuates lipopolysaccharide-induced cardiac dysfunction by inhibiting inflammation and regulating autophagy Exp Ther Med. 2020 Aug;20(2):1064-1072.;doi: 10.3892/etm.2020.8755
- 175.Jinxin Zhou,Ya Zhang,Shiyue Li,Qian Zhou,Yuanfu Lu,Jingshan Shi,Jie Liu,Qin Wu,Shaoyu Zhou.Dendrobium nobile Lindl. alkaloids-mediated protection against CCl4-induced liver mitochondrial oxidative damage is dependent on the activation of Nrf2 signaling pathway Biomed Pharmacother. 2020 Sep;129:110351.;doi: 10.1016/j.biopha.2020.110351
- 176.Sha Yang,Xiaoning Wang,Cancan Duan,Jianyong Zhang.A novel approach combining metabolomics and molecular pharmacology to study the effect of Gei Herba on mouse hematopoietic function Biomed Pharmacother. 2020 Sep;129:110437.;doi: 10.1016/j.biopha.2020.110437
- 177.Ming Qi,Jing Wang,Bie Tan,Simeng Liao,Cimin Long,Yulong Yin.Postnatal growth retardation is associated with intestinal mucosa mitochondrial dysfunction and aberrant energy status in piglets J Cell Mol Med. 2020 Sep;24(17):10100-10111.;doi: 10.1111/jcmcm.15621
- 178.Tao Zhang,Chun-Feng Liu,Tie-Ning Zhang,Ri Wen,Wen-Liang Song.Overexpression of Peroxisome Proliferator-Activated Receptor γ Coactivator 1 α Protects Cardiomyocytes from Lipopolysaccharide-Induced Mitochondrial Damage and Apoptosis Inflammation. 2020 Oct;43(5):1806-1820.;doi: 10.1007/s10753-020-01255-4
- 179.Shuo Huang,Yanyan Tang,Tianjun Liu,Ning Zhang,Xueyan Yang,Dingwei Yang,Ge Hong.A Novel Antioxidant Protects Against Contrast Medium-Induced Acute Kidney Injury in Rats Front Pharmacol. 2020 Nov 27;11:599577.;doi: 10.3389/fphar.2020.599577
- 180.Sheng Zhang,Yanping Xu,Jinyuan Zhu,Jinlan Ma,Qingsheng Niu,Xiaohong Wang.Carbon monoxide attenuates LPS-induced myocardial dysfunction in rats by regulating the mitochondrial dynamic equilibrium Eur J Pharmacol. 2020 Dec 15;889:173726.;doi: 10.1016/j.ejphar.2020.173726
- 181.Xuefei Yu,Yanli Sun,Qing Cai,Xinyi Zhao,Ziyun Liu,Xindong Xue,Jianhua Fu.Hyperoxia exposure arrests alveolarization in neonatal rats via PTEN-induced putative kinase 1-Parkin and Nip3-like protein X-mediated mitophagy disorders Int J Mol Med. 2020 Dec;46(6):2126-2136.;doi: 10.3892/ijmm.2020.4766
- 182.Lujing Jiang,Xiaojian Yin,Ya-Hui Chen,Yan Chen,Wei Jiang,Hao Zheng,Feng-Qing Huang,Baolin Liu,Wei Zhou,Lian-Wen Qi,Jia Li.Proteomic analysis reveals ginsenoside Rb1 attenuates myocardial ischemia/reperfusion injury through inhibiting ROS production from mitochondrial complex I Theranostics. 2021 Jan 1;11(4):1703-1720.;doi: 10.7150/thno.43895
- 183.Peiven Zhang,Shen Chen,Huanwen Tang,Wanjun Fang,Ke Chen,Xu Chen.CoQ10 protects against acetaminophen-induced liver injury by

- enhancing mitophagy TOXICOL APPL PHARM. 2021 Jan 1;410:115355.;doi: 10.1016/j.taap.2020.115355
- 184.Xinyu Feng,Shanjie Wang,Xingjun Yang,Jie Lin,Wanrong Man,Yuan Dong,Yan Zhang,Zhijing Zhao,Haichang Wang,Dongdong Sun.Mst1 Knockout Alleviates Mitochondrial Fission and Mitigates Left Ventricular Remodeling in the Development of Diabetic Cardiomyopathy Front Cell Dev Biol. 2021 Jan 21;8:628842.;doi: 10.3389/fcell.2020.628842
- 185.Xue Han,Jielong Guo,Yunxiao Gao,Jicheng Zhan,Yilin You,Weidong Huang.Gentisic acid prevents diet-induced obesity in mice by accelerating the thermogenesis of brown adipose tissue Food Funct. 2021 Feb 15;12(3):1262-1270.;doi: 10.1039/d0fo02474k
- 186.Wanyue Huang,Menglin Liu,Bonan Xiao,Jian Zhang,Miao Song,Yanfei Li,Zheng Cao.Aflatoxin B 1 disrupts blood-testis barrier integrity by reducing junction protein and promoting apoptosis in mice testes Food Chem Toxicol. 2021 Feb;148:111972.;doi: 10.1016/j.fct.2021.111972
- 187.Lina Wu,Qingzhu Wang,Feng Guo,Xiaojun Ma,Jiao Wang,Yanyan Zhao,Yushan Yan,Guijun Qin.Involvement of miR-27a-3p in diabetic nephropathy via affecting renal fibrosis, mitochondrial dysfunction, and endoplasmic reticulum stress J Cell Physiol. 2021 Feb;236(2):1454-1468.;doi: 10.1002/jcp.29951
- 188.Wenxing Cui,Xun Wu,Yingwu Shi,Wei Guo,Jianing Luo,Haixiao Liu,Longlong Zheng,Yong Du,Ping Wang,Qiang Wang,Dayun Feng,Shunnan Ge,Yan Qu.20-HETE synthesis inhibition attenuates traumatic brain injury-induced mitochondrial dysfunction and neuronal apoptosis via the SIRT1/PGC-1 α pathway: A translational study CELL PROLIFERAT. 2021 Feb;54(2):e12964.;doi: 10.1111/cpr.12964
- 189.Yan Zhang,Xiaomei Yang,Shuang Wang,Shuang Song,Xiudong Yang.Gentiopicroside prevents alcoholic liver damage by improving mitochondrial dysfunction in the rat model Phytother Res. 2021 Apr;35(4):2230-2251.;doi: 10.1002/ptr.6981

Version 2021.09.01