碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

BeyoFast™ SYBR Green qPCR Mix (2X)

产品编号	产品名称	包装
D7260-1ml	BeyoFast™ SYBR Green qPCR Mix (2X)	1ml
D7260-5ml	BeyoFast™ SYBR Green qPCR Mix (2X)	5ml
D7260-25ml	BeyoFast™ SYBR Green qPCR Mix (2X)	25ml

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1. Ma Y,Yang H,Huang J.Icariin ameliorates dexamethasone-induced bone deterioration in an experimental mouse model via activation of microRNA-186 inhibition of cathepsin K.Mol Med Rep. 2018 Jan;17(1):1633-1641
2. Zhong L, Sun S, Yao S, Han X, Gu M, Shi J.Histone deacetylase 5 promotes the proliferation and invasion of lung cancer cells.2018 Oct;
3. Liu Y,Long YH,Wang SQ,Zhang YY,Li YF,Mi JS,Yu CH,Li DY,Zhang JH,Zhang XJ.JMJD6 regulates histone H2A.X phosphorylation and promotes autophagy in triple-negative breast cancer cells via a novel tyrosine kinase activity.Oncogene. 2019 Feb;38(7):980-997
4. Wei Y,Shao J,Shen H,Wang Y,Cao G,Chen W,Zhang J,Yin L.The inhibitory role of recombinant P-selectin glycoprotein ligand immunoglobulin G on portal vein thrombosis based on a novel rat model.Exp Ther Med. 2019 May;17(5):3589-3597
5. Gu Y,Zhang J,Zhang X,Liang G,Xu T,Niu W.Three-dimensional Printed Mg-Doped β -TCP Bone Tissue Engineering Scaffolds: Effects of Magnesium Ion Concentration on Osteogenesis and Angiogenesis In Vitro.Tissue Eng Regen Med. 2019 Jun 17;16(4):415-429
6. Chen Z,Du Y,Liu X,Chen H,Weng X,Guo J,Wang M,Wang X,Wang L.EZH2 inhibition suppresses bladder cancer cell growth and metastasis via the JAK2/STAT3 signaling pathway.Oncol Lett. 2019 Jul;18(1):907-915
7. Lou K,Huang P,Ma H,Wang X,Xu H,Wang W.Orlistat increases arsenite tolerance in THP-1 derived macrophages through the up-regulation of ABCA1.Drug Chem Toxicol. 2019 Oct 31:1-9
8. Lingchen Kong,Peng Wu,Jianzhong Li.miR-331 inhibits CLDN2 expression and may alleviate the vascular endothelial injury induced by sepsis Exp Ther Med. 2020 Aug;20(2):1343-1352.;doi: 10.3892/etm.2020.8854

Version 2021.09.01