碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

BeyoFast™ SYBR Green qPCR Mix (2X, Low ROX)

产品编号	产品名称	包装
D7262-1ml	BeyoFast™ SYBR Green qPCR Mix (2X, Low ROX)	1ml
D7262-5ml	BeyoFast™ SYBR Green qPCR Mix (2X, Low ROX)	5ml
D7262-25ml	BeyoFast™ SYBR Green qPCR Mix (2X, Low ROX)	25ml

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1. Yingchao Cui, Jingfeng Gao, Da Zhang, Yifan Zhao, Yuwei Wang. Rapid start-up of partial nitrification process using benzethonium chloride-a novel nitrite oxidation inhibitor. *BIORESOURCE TECHNOL.* 2020 Nov;315:123860.
2. Shaoyuan Yang, Hairun Pei, Xiaoying Zhang, Qiang Wei, Jia Zhu, Jimin Zheng, Zongchao Jia. Characterization of DicB by partially masking its potent inhibitory activity of cell division. *Open Biol.* 2016 Jul;6(7):160082.
3. Zhang C, Wang Y. Metformin attenuates cells stemness and epithelial-mesenchymal transition in colorectal cancer cells by inhibiting the Wnt3a/ β -catenin pathway. *Mol Med Rep.* 2019 Feb 19(2):1203-1209.
4. Shi Q, Zhao L, Xu C, Zhang L, Zhao H. High Molecular Weight Hyaluronan Suppresses Macrophage M1 Polarization and Enhances IL-10 Production in PM2.5-Induced Lung Inflammation. *Molecules.* 2019 May 7 24(9). pii: E1766.
5. Jing Jing Fu, Xiao Hui Lv, Lin Xiang Wang, Xiu He, Yuan Li, Ling Yu, Chang Ming Li. Cutting and Bonding Parafilm® to Fast Prototyping Flexible Hanging Drop Chips for 3D Spheroid Cultures. *Cell Mol Bioeng.* 2020 Oct 27;14(2):187-199.
6. Jianbo He, Hongjuan Fu, Cancan Li, Zhihui Deng, Hui Chang. Eriodictyol inhibits breast carcinogenesis by targeting circ_0007503 and repressing PI3K/Akt pathway. *Phytomedicine.* 2022 Jul 20:102:154159.
7. Jiong Jiang, Yan Cheng, Shejiao Dai, Baicang Zou, Xiaoyan Guo. Suppression of rhomboid domain-containing 1 produces anticancer effects in pancreatic adenocarcinoma through affection of the AKT/GSK-3 β / β -catenin pathway. *Environ Toxicol.* 2022 Aug;37(8):1944-1956.
8. Xueyan Wang, Chengwei Duan, Yu Li, Hongjian Lu, Kongwei Guo, Xin Ge, Tianpeng Chen, Yanxing Shang, Hongbin Liu, Dongmei Zhang. Sodium butyrate reduces overnutrition-induced microglial activation and hypothalamic inflammation. *Int Immunopharmacol.* 2022 Oct:111:109083.

Version 2024.03.12