碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301 或 800-8283301  
 订货 e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

## BeyoGel™ SDS-PAGE预制胶(Tris-Gly, 8%, 12孔)

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1. Dan Wu, Lushan Xu, Wen-Min Cai, Shi-Yu Zhan, Guoqiang Wan, Yun Xu, Yun Stone Shi . A splicing-dependent ER retention signal regulates surface expression of the mechanosensitive TMEM63B cation channel J Biol Chem. 2023 Jan;299(1):102781.
2. Zhiqiang Ma, Yan Zuo, Wei Wang . Ginsenoside Rg3 inhibits renal cell carcinoma cell migration, invasion, colony formation, and tube formation and enhances apoptosis through promoting the DNA demethylation and histone acetylation J Pharm Pharmacol. 2023 Jan 31;75(1):76-86.
3. Lan Zhan, Zhuang Mu, Hao Jiang, Shicun Zhang, Yu Pang, Hongwei Jin, Jing Chen, Cuiying Jia, Hongyan Guo . MiR-21-5p protects against ischemic stroke by targeting IL-6R Ann Transl Med. 2023 Jan 31;11(2):101.
4. Xiaoyun Bin, Yu Luo, Zefeng Sun, Chaoqun Lin, Peng Huang, Zhenbo Tu, Ling Li, Cong Qu, Jiamin Long, Sufang Zhou . The Role of H2-Calponin Antigen in Cancer Metastasis: Presence of Autoantibodies in Liver Cancer Patients Int J Mol Sci. 2023 Jun 7;24(12):9864.
5. Chen Chen, Jianan Ma, Chenchen Pi, Wei Huang, Tao Zhang, Cong Fu, Wentao Liu, Yong-Guang Yang . PPAR $\delta$  inhibition blocks the induction and function of tumor-induced IL-10+ regulatory B cells and enhances cancer immunotherapy Cell Discov. 2023 Jun 8;9(1):54.
6. Xi Wen, Yixiang Wang, Yan Gu . Transferrin promotes chondrogenic differentiation in condylar growth through inducing autophagy via ULK1-ATG16L1 axis Clin Sci (Lond). 2023 Sep 27;137(18):1431-1449.
7. Lihui Lin, Yuting Liang, Tianyu Cao, Yuji Huang, Weize Li, Jia Li, Juan Wang, Xia Peng, Yiqin Ge, Yanning Li, Li Li . Transcriptome profiling and ceRNA network of small extracellular vesicles from resting and degranulated mast cells Epigenomics. 2023 Sep;15(17):845-862.
8. Wenjie Han, Haojun Liang, Jianqiang Bao . Efficient Large DNA Fragment Knock-in by Long dsDNA with 3'-Overhangs Mediated CRISPR Knock-in (LOCK) in Mammalian Cells Bio Protoc. 2023 Oct 20;13(20):e4853.
9. Wenmin Ji, Wenyan Wang, Yanfang Wei . TP53TG1/STAT axis is involved in the development of colon cancer through affecting PD-L1 expression and immune escape mechanism of tumor cells Am J Cancer Res. 2023 Nov 15;13(11):5218-5235.
10. Yin Liang, Xiao-Dan Xu, Xi Xu, Yang-Bo Cai, Zi-Xian Zhu, Lin Zhu, Kun Ren . Linc00657 promoted pyroptosis in THP-1-derived macrophages and exacerbated atherosclerosis via the miR-106b-5p/TXNIP/NLRP3 axis Int J Biol Macromol. 2023 Dec 31;253(Pt 4):126953.
11. Zirui He, Chuanhao Sun, Yifan Ma, Xi Chen, Ying Wang, Kai Chen, Fangru Xie, Yan Zhang, Yuan Yuan, Changsheng Liu . Rejuvenating Aged Bone Repair through Multihierarchy Reactive Oxygen Species-Regulated Hydrogel Adv Mater. 2024 Mar;36(9):e2306552.
12. Wenzhi Chen, Shishi Jiang, Shu Li, Cheng Li, Renshi Xu . OSMR is a potential driver of inflammation in amyotrophic lateral sclerosis Neural Regen Res. 2024 Nov 1;19(11):2513-2521.

Version 2025.04.03