碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

BeyoMag™ Protein G Plus磁珠

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++++, Strong Binding
 +++~++++, Medium Binding
 +, Weak Binding
 +/-, Weak or No Binding
 -, No Binding

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Characteristics	Description
Product content	10mg/ml magnetic beads in specific protective buffer

Beads size	~200nm
Magnetization	Superparamagnetic
Coupled protein	Recombinant Protein G Plus
M.W. of protein	~25kDa (Protein G Plus)
Antibody concentration	≥0.6mg Protein G Plus per ml beads
Binding capacity	≥ 0.9mg human IgG per ml beads
Specificity	Antibodies from many different species, including mouse, human, rabbit, cow, goat and sheep
Elution method	Elution with acid, competing peptide or SDS-PAGE loading buffer
Application	IP, Co-IP, Ch-IP, antibody purification

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Problem	Possible Causes	Solution
Very few or no target protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the target protein by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the transfer was adequate by using prestained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
Background is too high.	Proteins bind nonspecifically to the antibody, insufficient washing on magnetic beads, or the microcentrifuge tubes.	1. Pre-clear lysate with Normal IgG to remove nonspecific binding proteins. 2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before separation.
	Washes are insufficient.	1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Increase the salt and/or detergent concentrations in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	Purify the target protein at lower temperature, such as 4°C.
	Protein degradation due to proteases activity during purification process.	Add protease inhibitors to cell lysate.
	Non-specific binding.	1. Prepare cell lysate again. 2. Add additional wash steps.

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产品编号	产品名称	包装
P2102-1ml	BeyoMag™ Protein A磁珠	1ml
P2102-5ml	BeyoMag™ Protein A磁珠	5ml
P2105-1ml	BeyoMag™ Protein G磁珠	1ml
P2105-5ml	BeyoMag™ Protein G磁珠	5ml
P2106-1ml	BeyoMag™ Protein G Plus磁珠	1ml
P2106-5ml	BeyoMag™ Protein G Plus磁珠	5ml
P2108-1ml	BeyoMag™ Protein A+G磁珠	1ml
P2108-5ml	BeyoMag™ Protein A+G磁珠	5ml

Version 2021.10.05