## Mouse IL-6 ELISA Kit

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- Wang J, Zhang Y . Neuroprotective effect of berberine agonist against impairment of learning and memory skills in severe traumatic brain injury via Sirt1/p38 MAPK expression. Mol Med Rep. 2018 May;17(5):6881-6886.
- Jin Y, Tu Q, Liu M . MicroRNA-125b regulates Alzheimer's disease through SphK1 regulation. Mol Med Rep. 2018 Aug;18(2):2373-2380.
- Shao L, Liu C, Wang S, Liu J, Wang L, Lv L, Zou Y . The impact of exogenous CO releasing molecule CORM-2 on inflammation and signaling of orthotopic lung cancer. Oncol Lett. 2018 Sep;16(3):3223-3230.
- Sun J, Huang N, Ma W, Zhou H, Lai K . Protective effects of metformin on lipopolysaccharide-induced airway epithelial cell injury via NF- $\kappa$ B signaling inhibition. Mol Med Rep. 2019 Mar 19(3):1817-1823.
- Chang X, Li S, Xue XD, Chang F . Propranolol regulates ERK1/2 signaling pathway and promotes chronic wound healing in diabetic rats. EUR REV MED PHARMACO. 2019 May 23(10):4498-4506.
- Rao G, Zhang W, Song S . MicroRNA-217 inhibition relieves cerebral ischemia/reperfusion injury by targeting SIRT1. Mol Med Rep. 2019 Aug 20(2):1221-1229.
- Sui C, Zhang L, Hu Y . MicroRNA-let-7a inhibition inhibits LPS-induced inflammatory injury of chondrocytes by targeting IL6R. Mol Med Rep. 2019 Sep 20(3):2633-2640.
- Ma X, Liu X, Feng J, Zhang D, Huang L, Li D, Yin L, Li L, Wang XZ . Fraxin Alleviates LPS-Induced ARDS by Downregulating Inflammatory Responses and Oxidative Damages and Reducing Pulmonary Vascular Permeability. Inflammation. 2019 Oct 42(5):1901-1912.
- Chen L, Liu Z . Downregulation of FSTL1 attenuates the inflammation injury during Streptococcus pneumoniae infection by inhibiting the NLRP3 and TLR4/NF- $\kappa$ B signaling pathway. Mol Med Rep. 2019 Dec 20(6):5345-5352.
- Shao-Qing Chen, Yan-Qing Song, Cheng Wang, Shan Tao, Fang-Ying Yu, Hai-Ya Lou, Fu-Qiang Hu, Hong Yuan . Chitosan-modified lipid nanodrug delivery system for the targeted and responsive treatment of ulcerative colitis Carbohydr Polym. 2020 Feb 15;230:115613.
- Yuting Pan, Weihao Li, Ying Feng, Jing Xu, Huifang Cao . Edaravone attenuates experimental asthma in mice through induction of HO-1 and the Keap1/Nrf2 pathway Exp Ther Med. 2020 Feb;19(2):1407-1416.
- Yufeng Lou, Zhenrong Huang . microRNA-15a-5p participates in sepsis by regulating the inflammatory response of macrophages and targeting TNIP2 Exp Ther Med. 2020 Apr;19(4):3060-3068.
- Zengtao Hou, Jian Zhang, Ke Yu, Feng Song . Irisin ameliorates the postoperative depressive-like behavior by reducing the surface expression of epidermal growth factor receptor in mice Neurochem Int. 2020 May;135:104705.
- Hui Chen, Yao Wang, Fang-Zhou Jiao, Fan Yang, Xun Li, Lu-Wen Wang . Sinomenine Attenuates Acetaminophen-Induced Acute Liver Injury by Decreasing Oxidative Stress and Inflammatory Response via Regulating TGF- $\beta$ /Smad Pathway in vitro and in vivo Drug Des Devel Ther. 2020 Jun 17;14:2393-2403.
- Shaohong Wu, Xiaoping Huang, Wei Sun, Li Chen, Yan Huang, Yan

- Wang, Erping Luo, Ailan Qin, Weifeng Zhao, Jianhe Gan . Role of the microRNA-214/Bax axis in the progression of acute liver failure *Mol Med Rep.* 2020 Jul;22(1):117-126.
16. Lei Shi, Shun Guo, Song Zhang, Xiaobo Gao, An Liu, Qinhui Wang, Tian Zhang, Yan Zhang, Aidong Wen . Glycyrrhetic acid attenuates disturbed vitamin a metabolism in non-alcoholic fatty liver disease through AKR1B10 *Eur J Pharmacol.* 2020 Sep 15;883:173167.
  17. Zijing Yang, Li Wang, Yingjun Hu, Feixiang Wang . Butorphanol protects PC12 cells against OGD/R-induced inflammation and apoptosis *Mol Med Rep.* 2020 Sep;22(3):1969-1975.
  18. Xianliang Zeng, Ping Chen, Li Zhao, Sizeng Chen . Acylated and unacylated ghrelin relieve cancer cachexia in mice through multiple mechanisms *Chin J Physiol.* Sep-Oct 2020;63(5):195-203.
  19. Fei Liao, Xiaohong Lu, Weiguo Dong . Exosomes derived from T regulatory cells relieve inflammatory bowel disease by transferring miR-195a-3p *IUBMB Life.* 2020 Oct 26.
  20. Peng-Cheng Xing, Peng An, Guo-Yong Hu, Dong-Lian Wang, Min-Jie Zhou . LncRNA MIAT Promotes Inflammation and Oxidative Stress in Sepsis-Induced Cardiac Injury by Targeting miR-330-5p/TRAF6/NF- $\kappa$ B Axis *Biochem Genet.* 2020 Oct;58(5):783-800.
  21. Cong Fu, Qiancheng Xu, Shengxing Tang, Yuhan Cao, Can Liu, Yihua Wang, Yan Qian, Fei Shi, Jingmin Gui, Qun Fan, Yang Ling . The mobilization of splenic reservoir myeloid-derived suppressor cells in sepsis-induced myocardial injury *Am J Transl Res.* 2020 Nov 15;12(11):7114-7126.
  22. Xin Xu, Quanli Zhu, Guoliang Li, Junjian Ma, Zhijian Pan, Wei Wu . Protective Role of Somatostatin in Sepsis-Induced Intestinal Barrier Dysfunction through Inhibiting the Activation of NF-  $\kappa$  B Pathway *Gastroenterol Res Pract.* 2020 Dec 10;2020:2549486.
  23. Yunfu Yu, Rui Yan, Xiaozhen Chen, Tao Sun, Jifeng Yan . Paeonol suppresses the effect of ox-LDL on mice vascular endothelial cells by regulating miR-338-3p/TET2 axis in atherosclerosis *Mol Cell Biochem.* 2020 Dec;475(1-2):127-135.
  24. Zhongyang Lv, Xingquan Xu, Ziyang Sun, Yannick Xiaofan Yang, Hu Guo, Jiawei Li, Kuoyang Sun, Rui Wu, Jia Xu, Qing Jiang, Shiro Ikegawa, Dongquan Shi. TRPV1 alleviates osteoarthritis by inhibiting M1 macrophage polarization via Ca<sup>2+</sup>/CaMKII/Nrf2 signaling pathway. *Cell Death Dis.* 2021 May 18;12(6):504.
  25. Yuanyuan Lei, Li Tang, Shuang Liu, Shiping Hu, Lingyi Wu, Yaojiang Liu, Min Yang, Shengjie Huang, Xuefeng Tang, Tao Tang, Xiaoyan Zhao, Israel Vlodavsky, Shuo Zeng, Bo Tang, Shiming Yang. Parabacteroides produces acetate to alleviate heparanase-exacerbated acute pancreatitis through reducing neutrophil infiltration. *Microbiome.* 2021 May 20;9(1):115.
  26. Bin Ye, Liqun Lai. Yu Shi An Chang Fang Ameliorates TNBS-Induced Colitis in Mice by Reducing Inflammatory Response and Protecting the Intestinal Mucosal Barrier. *Evid Based Complement Alternat Med.* 2021 May 4;2021:8870901.
  27. Di Zhao, Bo Yang, Chen Ye, Shaoyi Zhang, Xiaoqiong Lv, Qiyi Chen. Enteral nutrition ameliorates the symptoms of Crohn's disease in mice via activating special pro-resolving mediators through innate lymphoid cells. *Innate Immun.* 2021 Oct;27(7-8):533-542.
  28. Huajie Mao, Bin Han, Hao Li, Yiqing Tao, Weigang Wu. FABP4 knockdown suppresses inflammation, apoptosis and extracellular matrix degradation in IL-1 $\beta$ -induced chondrocytes by activating PPAR  $\gamma$  to regulate the NF- $\kappa$ B signaling pathway. *Mol Med Rep.* 2021 Dec;24(6):855.
  29. Biwei Pei, Na Zhang, Tingting Pang, Gengyun Sun. Linagliptin ameliorates pulmonary fibrosis in systemic sclerosis mouse model via inhibition of endothelial-to-mesenchymal transition. *Mol Cell Biochem.* 2022 Apr;477(4):995-1007.
  30. Haifeng Guo, Xia Guo, Shiting Jiang. Long non-coding RNA lincRNA-erythroid prosurvival (EPS) alleviates cerebral ischemia/reperfusion injury by maintaining high-temperature requirement protein A1 (Htra1) stability through recruiting hete. *Bioengineered.* 2022 May;13(5):12248-12260.
  31. Rongxin Chen, Hailin Xu, Zebin Guo, Peng Zhang, Jianxia Chen, Zheng Chen. CID16020046, a GPR55 antagonist, attenuates sepsis-induced acute kidney injury. *Mol Med Rep.* 2022 May;25(5):155.
  32. Dongbei Shen, Caie Wu, Gongjian Fan, Tingting Li, Jinfeng Dou, Jinpeng Zhu, Chunmei Li, Xiuying Kou. Jujube peel polyphenols synergistically inhibit lipopolysaccharide-induced inflammation through multiple signaling pathways in RAW 264.7 cells. *Food Chem Toxicol.* 2022 Jun;164:113062.
  33. Shucheng Zhang, Wenhua Jin, Weibo Zhang, Fazheng Ren, Pengjie Wang, Ning Liu. Pea Albumin Attenuates Dextran Sulfate Sodium-Induced Colitis by Regulating NF- $\kappa$ B Signaling and the Intestinal Microbiota in Mice. *Nutrients.* 2022 Sep 1;14(17):3611.
  34. Jing Tian, Dapeng Zhou, Liangbi Xiang, Bing Xie, Baichuan Wang, Yang Li, Xinwei Liu. Calycosin represses AIM2 inflammasome-mediated inflammation and pyroptosis to attenuate monosodium urate-induced gouty arthritis through NF- $\kappa$ B and p62-Keap1 pathways. *Drug Dev Res.* 2022 Nov;83(7):1654-1672.
  35. Jun Sun, Ping Zhao, Xufeng Ding, Fang Li, Jie Jiang, Hua Huang, Lijiang Ji. Cayratia japonica Prevents Ulcerative Colitis by Promoting M2 Macrophage Polarization through Blocking the TLR4/MAPK/NF-  $\kappa$  B Pathway. *Mediators Inflamm.* 2022 Dec 30;2022:1108569.
  36. Ciyu Yao, Wenzhen Ye, Mengxue Chen. Inhibition of Mast Cell Degranulation in Atopic Dermatitis by Celastrol through Suppressing MRGPRX2. *Dis Markers.* 2023 Jan 18;2023:9049256.
  37. Qiu-Jiao Zhu, Jian Wang, Ying Li, Zhen-Jiang Bai, Xu-Bei Guo, Tao Pan. PRKCA Promotes Mitophagy through the miR-15a-5p/PDK4 Axis to Relieve Sepsis-Induced Acute Lung Injury. *Infect Immun.* 2023 Jan 24;91(1):e0046522.
  38. Pan Huang, Fan Hu, Zhi-Bo Yang, Yi Pan, Rong Zhou, Yi-Ning Yan, Hai-Zhen Wang, Chang Wang. Matrine regulates Th1/Th2 inflammatory responses by inhibiting the Hsp90/NF- $\kappa$ B signaling axis to alleviate atopic dermatitis. *Kaohsiung J Med Sci.* 2023 May;39(5):501-510.
  39. Zhiwen Zhang, Quan Guo, Zhenzhou Zhao, Ming Nie, Qingbo Shi, En Li, Kaiyuan Liu, Haosen Yu, Lixin Rao, Muwei Li. DNMT3B activates FGFR3-mediated endoplasmic reticulum stress by regulating PTPN2 promoter methylation to promote the development of atherosclerosis. *FASEB J.* 2023 Aug;37(8):e23085.
  40. Lin Lin, Kuimiao Deng, Zongrong Gong, Huifeng Fan, Dongwei Zhang, Gen Lu. Sinensetin Attenuated Macrophagic NLRP3 Inflammasomes Formation via SIRT1-NRF2 Signaling. *ACS Omega.* 2023 Sep 7;8(37):33514-33525.
  41. Wei Chen, Juan Shi, Wei Zhang. Administration of Akebia Saponin D Improved Blood Lipid Levels and Pregnancy Outcomes in Mice with Gestational Diabetes Mellitus. *Balkan Med J.* 2023 Oct 20;40(6):422-429.
  42. Zheqian Wu, Yong Wang, Shijie Lu, Lili Yin, Lihua Dai. SIRT3 alleviates sepsis-induced acute lung injury by inhibiting pyroptosis via regulating the deacetylation of FoxO3a. *Pulm Pharmacol Ther.* 2023 Oct;82:102244.
  43. Shuxiang Chen, Xiaotong Zhu, Wenhuan Ou, Le Kang, Jian Situ, Zhipeng Liao, Li Huang, Weizhong Qi, Songjia Ni. ETS2 overexpression ameliorates cartilage injury in osteoarthritis by the ETS2/miR-155/STAT1/DNMT1 feedback loop pathway. *Biochim Biophys Acta Gene Regul Mech.* 2023 Dec;1866(4):194965.
  44. Jingying Liu, Meijuan Miao, Fujiang Wei. Angelicin Alleviates Maternal Isoflurane Exposure-Induced Offspring Cognitive Defects

Through the Carbonic Anhydrase 4/Aquaporin-4 Pathway. Mol Biotechnol. 2024 Jan;66(1):34-43.

45. Qian Yan, Zixing Liu, Yujing Chen, Xinxin Zhang, Wenjiang Zheng, Xiaohong Liu, Huiting Huang, Qiong Liu, Yong Jiang, Shaofeng Zhan, Xiufang Huang. ITGAM-macrophage modulation as a potential

strategy for treating neutrophilic Asthma: insights from bioinformatics analysis and in vivo experiments. Apoptosis. 2024 Apr;29(3-4):393-411.

Version 2025.02.07