



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301 或 800-8283301
订货 e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

Beyozol (总RNA抽提试剂)

产品编号	产品名称	包装
R0011	Beyozol (总RNA抽提试剂)	100ml

产品简介:

- 碧云天生产的Beyozol是一种用于细胞或组织总RNA抽提的试剂。本产品采用和Invitrogen公司的TRIzol完全相同的原理和方法, 抽提的方法和步骤完全相同。
- Beyozol的颜色和TRIzol相同, 加入氯仿后上层呈无色, 下层呈紫红色, 便于吸取上层水相。
- Beyozol对动植物细胞或组织及细菌的总RNA抽提均适用。
- Beyozol可以抽提长达15 kb的RNA, 也可以抽提microRNA等小RNA。抽提小RNA时宜-70°C沉淀过夜。
- Beyozol抽提所得RNA无DNA和蛋白污染。一般所得RNA溶于DEPC水后的A260/280值为1.8-2.0。
- 裂解细胞或和组织共匀浆时, Beyozol可以保持样品中RNA的完整性, 即可以有效抑制RNA的降解。
- 每一百万细胞用Beyozol抽提可得5-15µg RNA; 每毫克组织用Beyozol抽提可得1-10µg RNA。产量因细胞和组织不同而异。
- 抽提两个样品约需一小时。
- Beyozol抽提所得RNA可直接用于Northern, 点杂交, 纯化mRNA, 体外翻译, RNase protection assay, cDNA克隆, 以及RT-PCR; 也可以用于基因表达芯片分析、高通量测序(deep sequencing)等对RNA质量要求较高的情况。
- 碧云天生产的Beyozol(R0011)和Trizol(R0016)的成分有细微差别, 实际抽提效果无任何显著差异。
- 每100ml Beyozol可以抽提100个六孔板中的样品或100个50-80mg的组织样品。

包装清单:

产品编号	产品名称	包装
R0011	Beyozol	100ml
—	说明书	1份

保存条件:

4°C保存, 一年有效。

注意事项:

- 需自备氯仿, 异丙醇, DEPC, 75%乙醇(DEPC水配制), 和DEPC水。DEPC(ST036), DEPC水(R0021)可向碧云天订购。
- 所有离心管, 枪头及相关溶液都必须无RNA酶污染。耐高温器物可150°C烘烤4小时以去除RNA酶, 其它器物去除RNA酶可考虑用0.01%的DEPC水浸泡过夜, 然后灭菌, 烘干。溶液需用DEPC水配制。加0.01%(体积比) diethylpyrocarbonate(DEPC)至重蒸水或Milli-Q级水中, 处理过夜, 灭菌即成DEPC水。
- 使用冻存的细胞或组织抽提总RNA的效果通常比新鲜的细胞或组织差一些。因为在细胞或组织冻融过程中一些细胞或组织内的RNase会被释放出来并剪切样品。如果不能及时抽提RNA, 推荐先加入适量Beyozol, 并裂解样品后冻存。
- 必须戴一次性手套操作, 且尽量不要对着RNA样品呼气或说话, 以防RNA酶污染。建议戴一次性口罩操作。
- Beyozol含有毒物质苯酚, 避免接触皮肤或吸入。为防止溅入眼睛, 请戴防护眼镜或使用透明保护屏。如皮肤接触Beyozol, 请立即用大量去垢剂和水冲洗, 如仍有不适, 请听取医生意见。
- Beyozol抽提总RNA的同时, 理论上也可抽提蛋白和DNA, 但未经测试。有兴趣者可按Invitrogen公司的TRIzol的操作步骤进行操作。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. 细胞裂解或组织匀浆。

a. 贴壁细胞

吸尽培养液, 每10平方厘米细胞加入1ml Beyozol。一般六孔板每孔加1ml Beyozol, 12孔板每孔加0.5ml Beyozol。晃动3-5下, 再用枪吹打2-3下, 确保全部裂解, 然后吸至离心管中。

b. 悬浮细胞

离心收集细胞, 吸尽液体, 每五百万至一千万动植物或酵母细胞, 或一千万细菌, 加入1ml Beyozol。用枪吹打或适当vortex, 确保全部裂解。某些酵母和细菌如裂解不充分, 可用匀浆器匀浆, 确保全部裂解。

c. 组织

- 先将组织剪切成小块，放入普通玻璃匀浆器内。每50mg-80mg组织加入1ml Beyozol，匀浆。对于RNA完整性要求比较高的情况，推荐先液氮冷冻组织块，然后在低温下用研钵研碎组织，随后再加入Beyozol进行总RNA抽提。
- 对于某些蛋白，多糖或脂含量很高的细胞或组织，Beyozol裂解后可能会有不溶物或油脂状漂浮物。需12,000g 4°C离心10分钟，然后吸取澄清的Beyozol裂解产物至一新的离心管中。
 - 室温放置5分钟，使样品充分裂解。
 - 每毫升 Beyozol加入0.2ml氯仿，vortex混匀或剧烈晃动15秒，室温放置2-3分钟。
 - 12,000g 4°C离心15分钟，然后吸取含总RNA的上层无色水相至一新的离心管中，每毫升Beyozol约可吸取0.5-0.55ml。
 - 按每毫升最初的Beyozol加入0.5ml异丙醇，颠倒数次混匀，室温沉淀10分钟。如果希望提取microRNA等小RNA，推荐-70°C沉淀过夜。
 - 12,000g 4°C离心10分钟，在管底可见RNA沉淀，弃上清。
 - 每毫升最初的Beyozol加入1ml 75%乙醇(DEPC水配制)，vortex或颠倒混匀。
 - 7,500g 4°C离心5分钟，弃上清。再用离心机甩一下 (>5,000rpm，离心1秒)，小心吸尽液体。
 - 待RNA略干后，加入20μl DEPC水溶解，-70°C冻存。注意，切勿让RNA过分干燥，否则将极难溶解，且测出的A260/280值会低于1.6。

使用本产品的文献：

- Xiong L, Liu Y, Zhu F, Lin J, Wen D, Wang Z, Bai J, Ge G, Xu C, Gu Y, Xu Y, Zhou J, Geng D . Acetyl-11-keto-β-boswellic acid attenuates titanium particle-induced osteogenic inhibition via activation of the GSK-3β/β-catenin signaling pathway. *Theranostics*. 2019 Sep 23 9(24):7140-7155.
- Haiqin Huang, Kang Sha, Hanitrimalala Veroniaina, Ziheng Wu, Zhenghong Wu, Xiaole Qi . Ca²⁺ participating self-assembly of an apoferritin nanostructure for nucleic acid drug delivery *Nanoscale*. 2020 Apr 3;12(13):7347-7357.
- Peng Zhao, Jianhua Cheng, Bin Li, Ding Nie, Chuzhong Li, Songbai Gui, Hongyun Wang, Yazhuo Zhang . Up-regulation of the expressions of MiR-149-5p and MiR-99a-3p in exosome inhibits the progress of pituitary adenomas *Cell Biol Toxicol*. 2021 Jan 5. doi: 10.1007/s10565-020-09570-0.
- Haiqin Huang, Shirui Yuan, Zhuo Ma, Peng Ji, Xiaonan Ma, Zhenghong Wu, Xiaole Qi . Genetic recombination of poly(l-lysine) functionalized apoferritin nanocages that resemble viral capsid nanometer-sized platforms for gene therapy *BIOMATER SCI-UK*. 2020 Mar 17;8(6):1759-1770.
- Qing Wang, Gaoran Ge, Xiaolong Liang, Jiayang Bai, Wei Wang, Wenhao Zhang, Kai Zheng, Sen Yang, Minggang Wei, Huilin Yang, Yaozeng Xu, Boling Liu, Dechun Geng . Punicalagin ameliorates wear-particle-induced inflammatory bone destruction by bi-directional regulation of osteoblastic formation and osteoclastic resorption *BIOMATER SCI-UK*. 2020 Sep 21;8(18):5157-5171.
- Ruilian Li, Fangyan Jia, Kaixuan Ren, Mai Luo, Xiaoyun Min, Shengxiang Xiao, Yumin Xia . Fibroblast growth factor inducible 14 signaling facilitates anti-dsDNA IgG penetration into mesangial cells *J Cell Physiol*. 2021 Jan;236(1):249-259.
- Yang D, Lin F, Huang Y, Ye J, Xiao M . Separation, purification, structural analysis and immune-enhancing activity of sulfated polysaccharide isolated from sea cucumber viscera. *Int J Biol Macromol*. 2019 Nov 9 pii: S0141-8130(19)35632-6.
- Li H, Zhang YY, Tan HW, Jia YF, Li D . Therapeutic effect of tripterine on adjuvant arthritis in rats. *J Ethnopharmacol*. 2008 Aug 13;118(3):479-84.
- ZHANG ZY, LIU XJ, LIU T, YAN L, WANG YJ, WANG CB . Polypeptide from *Chlamys farreri* inhibits UVB-induced apoptosis of HaCaT cells via iNOS/NO and HSP90. *Chinese Journal of Oceanology and Limnology*. 2009;27(3):594-9.
- Longxi P, Buwu F, Yuan W, Sinan G . Expression of p53 in the effects of artesunate on induction of apoptosis and inhibition of proliferation in rat primary hepatic stellate cells. *PLoS One*. 2011;6(10):e26500.
- Zhao YQ, Feng HW, Jia T, Chen XM, Zhang H, Xu AT, Zhang HL, Fan XL . Antiproliferative effects of celecoxib in Hep-2 cells through telomerase inhibition and induction of apoptosis. *ASIAN PAC J CANCER P* . 2014;15(12):4919-23.
- Wang X, Jiang Q, Wang W, Su L, Han Y, Wang C . Molecular mechanism of polypeptides from *Chlamys farreri* (PCF)'s anti-apoptotic effect in UVA-exposed HaCaT cells involves HSF1/HSP70, JNK, XO, iNOS and NO/ROS. *J PHOTOCH PHOTOBIO B*. 2014 Jan 5;130:47-56.
- Jin Q, Jiang S, Wu YL, Bai T, Yang Y, Jin X, Lian LH, Nan JX . Hepatoprotective effect of cryptotanshinone from *Salvia miltiorrhiza* in D-galactosamine/lipopolysaccharide-induced fulminant hepatic failure. *Phytomedicine*. 2014 Jan 15;21(2):141-7.
- Ye J, Han Y, Chen X, Xie J, Liu X, Qiao S, Wang C . l-Carnitine attenuates H₂O₂-induced neuron apoptosis via inhibition of endoplasmic reticulum stress. *Neurochem Int*. 2014 Sep 9;78C:86-95.
- Yu B, Dai CQ, Chen J, Deng L, Wu XL, Wu S, Zhao CL, Jiang ZY, Chen XY . Dysbiosis of gut microbiota induced the disorder of helper T cells in influenza virus-infected mice. *HUM VACC IMMUNOTHER*. 2015;11(5):1140-6.
- Zhao S, Wang Z, Chen J, Chen J . Preparation of heparan sulfate-like polysaccharide and application in stem cell chondrogenic differentiation. *CARBOHYD RES*. 2015 Jan 12;401:32-8.
- Zhao S, Deng C, Wang Z, Teng L, Chen J . Heparan sulfate 6-O-sulfotransferase 3 is involved in bone marrow mesenchymal stromal cell osteogenic differentiation. *BIOCHEMISTRY-MOSCOW+*. 2015 Mar;80(3):379-89.
- Wu Y, Wang M, Feng H, Peng Y, Sun J, Qu X, Li C . Lactate induces osteoblast differentiation by stabilization of HIF1α. *Mol Cell Endocrinol*. 2017 Sep 5;452:84-92.
- Chen X, Han R, Hao P, Wang L, Liu M, Jin M, Kong D, Li X . Nepetin inhibits IL-1β induced inflammation via NF-κB and MAPKs signaling pathways in ARPE-19 cells. *Biomed Pharmacother*. 2018 May;101:87-93.
- Zhou XQ, Yao ZW, Peng Y, Mao SS, Xu D, Qin XF, Zhang RJ . PQQ ameliorates D-galactose induced cognitive impairments by reducing glutamate neurotoxicity via the GSK-3β/Akt signaling pathway in mouse. *SCI REP-UK*. 2018 Jun 11;8(1):8894.
- Wang X, Hu H, Liu H . RNA binding protein Lin28B confers gastric cancer cells stemness via directly binding to NRP-1. *Biomed Pharmacother*. 2018 Aug;104:383-389.

22. Xiao H, Yang YJ, Lin YZ, Peng S, Lin S, Song ZY . Transcription factor Tbx18 induces the differentiation of c-kit+ canine mesenchymal stem cells (cMSCs) into SAN-like pacemaker cells in a co-culture model in vitro. *Am J Transl Res*. 2018 Aug 15;10(8):2511-2528. eCollection 2018.
23. Pang P, Yu B, Shi Y, Deng L, Xu H, Wu S, Chen X . Alteration of Intestinal Flora Stimulates Pulmonary microRNAs to Interfere with Host Antiviral Immunity in Influenza. *Molecules*. 2018 Nov 30;23(12). pii: E3151.
24. Yi J, Mei L, Li X, Zheng W, Li Y, Zhao Z . Effects of continuous and intermittent parathyroid hormone administration on midpalatal suture expansion in rats. *Arch Oral Biol*. 2019 Mar 99:161-168.
25. Zhou J, Ma W, Wang X, Liu H, Miao Y, Wang J, Du P, Chen Y, Zhang Y, Liu Z . Matrine Suppresses Reactive Oxygen Species (ROS)-Mediated MKKs/p38-Induced Inflammation in Oxidized Low-Density Lipoprotein (ox-LDL)-Stimulated Macrophages. *MED SCI MONITOR*. 2019 Jun 3 25:4130-4136.
26. Lyu Z, Gao X, Wang W, Dang J, Yang L, Yan M, Ali SA, Liu Y, Liu B, Yu M, Du L, Liu K . mTORC1-Sch9 regulates hydrogen sulfide production through the transsulfuration pathway. *AGING-US*. 2019 Oct 3 11(19):8418-8432.
27. Li J, Xiang X, Xu H, Shi Y . Cilostazol Promotes Angiogenesis and Increases Cell Proliferation After Myocardial Ischemia-Reperfusion Injury Through a cAMP-Dependent Mechanism. *Cardiovasc Eng Technol*. 2019 Dec 10(4):638-647.
28. Wei Huang, Chao Cheng, Wen-Shan Shan, Zhen-Fei Ding, Fu-En Liu, Wei Lu, Wei He, Jie-Gou Xu, Zong-Sheng Yin . Knockdown of SGK1 alleviates the IL-1 β -induced chondrocyte anabolic and catabolic imbalance by activating FoxO1-mediated autophagy in human chondrocytes *FEBS J*. 2020 Jan;287(1):94-107.
29. Shancheng Zhao, Zhen Wang . Changes in heparan sulfate sulfotransferases and cell-surface heparan sulfate during SKM-1 cells granulocytic differentiation and A549 cells epithelial-mesenchymal transition *GLYCOCONJUGATE J*. 2020 Apr;37(2):151-164.
30. Ying Peng, Dong Xu, Shishi Mao, Xingqin Zhou . Neurotoxicity and apoptosis induced by pyrroloquinoline quinone and its ester derivative on primary cortical neurons *Neurotoxicology*. 2020 May;78:47-56.
31. Wen-Liang Huang, Shu-Fen Wu, Song-Tao Xu, Yong-Chao Ma, Rui Wang, Shaoju Jin, Shan Zhou . Alicin enhances the radiosensitivity of colorectal cancer cells via inhibition of NF- κ B signaling pathway *J Food Sci*. 2020 Jun;85(6):1924-1931.
32. Zhang Feng, Liu Jie, Lv Guimin, Wang Xi . Mixed Lineage Leukemia 1 Promoted Neuron Apoptosis in Ischemic Penumbra via Regulating ASK-1/TNF- α Complex *Front Neuroanat*. 2020 Jul 24;14:36.
33. Jingyun Chen, Fangyan Jia, Kaixuan Ren, Mai Luo, Xiaoyun Min, Ping Wang, Shengxiang Xiao, Yumin Xia . Inhibition of suppressor of cytokine signaling 1 mediates the profibrotic effect of TWEAK/Fn14 signaling on kidney cells *Cell Signal*. 2020 Jul;71:109615.
34. Xiaoning Zhang, Xia Zhao, Yanqing Wang, Liqun Xing . Long Non-Coding RNA LINC00491 Contributes to the Malignancy of Non-Small-Cell Lung Cancer via Competitively Binding to microRNA-324-5p and Thereby Increasing Specificity Protein 1 Expression *Cancer Manag Res*. 2020 Aug 6;12:6779-6793.
35. Qin hong Xu, Lijing Zhou, Ganghua Yang, Fandi Meng, Yong Wan, Lin Wang, Lei Zhang . CircIL4R facilitates the tumorigenesis and inhibits ferroptosis in hepatocellular carcinoma by regulating the miR-541-3p/GPX4 axis *Cell Biol Int*. 2020 Nov;44(11):2344-2356.
36. Jianzhao Liao, Fan Yang, Wenlan Yu, Na Qiao, Hui Zhang, Qingyue Han, Lianmei Hu, Ying Li, Jianying Guo, Jiaqiang Pan, Zhaoxin Tang . Copper induces energy metabolic dysfunction and AMPK-mTOR pathway-mediated autophagy in kidney of broiler chickens *ECOTOX ENVIRON SAFE*. 2020 Dec 15;206:111366.
37. Jianguo Lin, Ying Peng, Qingzhu Liu, Ke Li, Gaochao Lv, Yann Seimbille, Gang Huang, Feng Gao, Ling Qiu . Pharmacological evaluation of imidazole-derived bisphosphonates on receptor activator of nuclear factor- κ B ligand-induced osteoclast differentiation and function *Chem Biol Drug Des*. 2021 Jan;97(1):121-133.
38. Xiaoxi Fan, Zhongwei Zhao, Jingjing Song, Dengke Zhang, Fazong Wu, Jianfei Tu, Min Xu, Jiansong Ji . LncRNA-SNHG6 promotes the progression of hepatocellular carcinoma by targeting miR-6509-5p and HIF1A *Cancer Cell Int*. 2021 Mar 4;21(1):150.
39. Xin He, Juan Wang, Ru Zhou, Shanshan Yu, Jue Jiang, Qi Zhou . Kinesin family member 23 exerts a protumor function in breast cancer via stimulation of the Wnt/ β -catenin pathway *Toxicol Appl Pharmacol*. 2022 Jan 15;435:115834.
40. Xiaoyu Zhou, Haitao Liu, Yingying Pang, Muqun Wang, Shengming Liu . UTMD-mediated delivery of miR-21-5p inhibitor suppresses the development of lung cancer *Tissue Cell*. 2022 Feb;74:101719.
41. Yongqiang Li, Lihong Yi, Xuelian Li, Ying Li . Benzodiazepines safeguards nerve cells from the toxicity of lidocaine via miR-133a-3p/EGFR pathway *Transpl Immunol*. 2022 Apr;71:101510.
42. Jie Pan, Xiaoping Cai, Xiao Zheng, Xiaoyu Zhu, Jihong Feng, Xiaoqi Wang . Luteolin inhibits viability, migration, angiogenesis and invasion of non-small cell lung cancer vascular endothelial cells via miR-133a-3p/purine rich element binding protein B-mediated MAPK and PI3K/ *Tissue Cell*. 2022 Apr;75:101740.
43. Hongmei Zhao, Lian Shi, Xiaohong Wang, Xiuli Yu, Danfeng Wang . Sp1 transcription factor represses transcription of phosphatase and tensin homolog to aggravate lung injury in mice with type 2 diabetes mellitus-pulmonary tuberculosis *Bioengineered*. 2022 Apr;13(4):9928-9944.
44. Hanxun Yao, Xiaobin Li, Xuyan Pan, Jie Xu, Shufa Zhao, Zhongzhou Su, Sheng Qiu . Fraxetin exerts anticancer effect in glioma by suppressing MiR-21-3p *Drug Dev Res*. 2022 Apr;83(2):501-511.
45. Qi Li, Dong-Na Zhou, Yi-Qing Tu, Xin-Wei Wu, Da-Qing Pei, Yun Xiong . Ketamine administration ameliorates anesthesia and surgery-induced cognitive dysfunction via activation of TRPV4 channel opening *Exp Ther Med*. 2022 Jun 1;24(1):478.
46. Danping Luan, Chun Jiang . The mechanism of lncRNA TALNEC2 regulating miR-19a-3p/JNK to alleviate cerebral ischemia injury in rats with acute cerebral infarction *Cell Mol Biol (Noisy-le-grand)*. 2022 Jun 30;68(6):17-24.
47. Nanzheng Chen, Qifei Wu, Guangjian Zhang, Junke Fu, Qianqian Geng, Yong Zhang . Deactivation of AKT/GSK-3 β -mediated Wnt/ β -catenin pathway by silencing of KIF26B weakens the malignant behaviors of non-small cell lung cancer *Tissue Cell*. 2022 Jun;76:101750.
48. Haiyan Zheng, Qingyuan Lin, Yamin Rao . A-Kinase Interacting Protein 1 Knockdown Restores Chemosensitivity via Inactivating PI3K/AKT and

- β -Catenin Pathways in Anaplastic Thyroid Carcinoma *Front Oncol.* 2022 Jul 28;12:854702.
49. Lei Zhang, Donglai Ma, Fujun Li, Gongcai Qiu, Dongsheng Sun, Zhaolin Zeng . Lnc-PKD2-2-3/miR-328/GPAM ceRNA Network Induces Cholangiocarcinoma Proliferation, Invasion and 5-FU Chemoresistance *Front Oncol.* 2022 Jul 29;12:871281.
50. Qingxiang Mao, Lixia Tian, Jianxiong Wei, Xiaoqiong Zhou, Hong Cheng, Xuan Zhu, Xiang Li, Zihao Gao, Xi Zhang, Lingli Liang . Transcriptome analysis of microRNAs, circRNAs, and mRNAs in the dorsal root ganglia of paclitaxel-induced mice with neuropathic pain *Front Mol Neurosci.* 2022 Aug 31;15:990260.
51. Yi-Ren Qin, Chi-Qian Ma, Jian-Hua Jiang, Da-Peng Wang, Quan-Quan Zhang, Mei-Rong Liu, Hong-Ru Zhao, Qi Fang, Yang Liu . Artesunate restores mitochondrial fusion-fission dynamics and alleviates neuronal injury in Alzheimer's disease models *J Neurochem.* 2022 Aug;162(3):290-304.
52. Weiwei Li, Ting Cheng, Tingting Jiang, Mengyue Zhou, Bowen Gong, Guode Zhao, Jing Li, Rong Tan, Xiaojing Yang, Kandarp Joshi, Ying Peng, Maosheng Cheng, Ting Liu, Dan Ohtan Wang, Jiang Zheng . Hepatic RNA adduction derived from metabolic activation of retrorsine in vitro and in vivo *Chem Biol Interact.* 2022 Sep 25;365:110047.
53. Li-Min Zhang, Yue Xin, Zhi-You Wu, Rong-Xin Song, Hui-Tao Miao, Wei-Chao Zheng, Yan Li, Dong-Xue Zhang, Xiao-Chun Zhao . STING mediates neuroinflammatory response by activating NLRP3-related pyroptosis in severe traumatic brain injury *J Neurochem.* 2022 Sep;162(5):444-462.
54. Yuling Guo, Yiping Gao, Shasha Liu . lncRNA XIST is associated with preeclampsia and mediates trophoblast cell invasion via miR-340-5p/KCNJ16 signaling pathway *Transpl Immunol.* 2022 Oct;74:101666.

Version 2024.03.12