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R0125	RNase, DNase and DNA Away	250ml
R0127	RNase, DNase, RNA and DNA Away	250ml
ST036	DEPC	10g

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1. Jiachuan Xiong, Li Ran, Yingguo Zhu, Yaqin Wang, Shaobo Wang, Yue Wang, Qigang Lan, Wenhao Han, Yong Liu, Yinghui Huang, Ting He, Yan Li, Li Liu, Jinghong Zhao, Ke Yang . DUSP2-mediated inhibition of tubular epithelial cell pyroptosis confers nephroprotection in acute kidney injury Theranostics. 2022 Jul 4;12(11):5069-5085.
2. Jinzhi Meng, Huawei Du, Haiyuan Lv, Jinfeng Lu, Jia Li, Jun Yao . Identification of the osteoarthritis signature gene PDK1 by machine learning and its regulatory mechanisms on chondrocyte autophagy and apoptosis Front Immunol. 2023 Jan 6;13:1072526.
3. Yu-Wei Yang, Xin Meng, Yuan-Yuan Meng, Hai-Kang Tang, Ming-Hui Cheng, Zi-Qi Zhang, Wen-Qing Xu, Wei Long . ceRNA regulatory network of FIH inhibitor as a radioprotector for gastrointestinal toxicity by activating the HIF-1 pathway Mol Ther Nucleic Acids. 2021 May 19;25:173-185.
4. Jiefu Zhu, Yafei Zhang, Lang Shi, Yao Xia, Hongchu Zha, Huimin Li, Zhixia Song . RP105 protects against ischemic and septic acute kidney injury via suppressing TLR4/NF- κ B signaling pathways Int Immunopharmacol. 2022 Aug;109:108904.
5. Wendan He, Xianlong Xie, Chenxi Li, Huang Ding, Jishi Ye . Adenosine A2A Receptor Antagonist Improves Cognitive Impairment by Inhibiting Neuroinflammation and Excitatory Neurotoxicity in Chronic Periodontitis Mice Molecules. 2022 Sep 23;27(19):6267.
6. Linxin Pan, Ying Hu, Cheng Qian, Yan Yao, Shuxian Wang, Wanrong Shi, Tao Xu . RNF2 mediates pulmonary fibroblasts activation and proliferation by regulating mTOR and p16-CDK4-Rb1 signaling pathway Inflamm Res. 2022 Nov;71(10-11):1283-1303.
7. Huangming Zhuang, Bin Li, Ting Xie, Changgeng Xu, Xunshan Ren, Fuze Jiang, Tianrun Lei, Panghu Zhou . Indole-3-aldehyde alleviates chondrocytes inflammation through the AhR-NF- κ B signalling pathway Int Immunopharmacol. 2022 Dec;113(Pt A):109314.
8. Liqun Wang, Wentao Zhang, Ruiyan Cen, Chenda Yue, Tianzhen Xiao, Yumeng Deng, Lingfei Li, Kedai Sun, Xia Lei . ALA-PDT regulates macrophage M1 polarization via ERK/MAPK-NLRP3 pathway to promote the early inflammatory response Lasers Surg Med. 2022 Dec;54(10):1309-1320.

Version 2024.03.12