



碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301 或 800-8283301
 订货 e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

ATP检测试剂盒

产品编号	产品名称	包装
S0026B	ATP检测试剂盒	200次

产品简介:

- ATP检测试剂盒(ATP Assay Kit)可以用于检测普通溶液、细胞或组织内的ATP(adenosine 5'-triphosphate)水平。
- ATP, 作为最重要的能量分子在细胞的各种生理、病理过程中起着重要作用。ATP水平的改变, 会影响细胞的功能。通常细胞在凋亡、坏死或处于一些毒性状态下, ATP水平会下降, 而高葡萄糖刺激等对于一些细胞可以上调细胞内ATP水平。通常ATP水平的下降表明线粒体的功能受损或下降, 在细胞凋亡时ATP水平的下降通常和线粒体的膜电位下降同时发生。
- 样品制备非常简单。本试剂盒提供了可以用于细胞和组织裂解的ATP检测裂解液, 简单裂解后即可用于ATP检测。无需进行高氯酸或三氯乙酸(TCA)抽提, 或样品裂解后的煮沸等繁琐操作。
- 本试剂盒根据萤火虫萤光素酶(firefly luciferase, 也称荧光素酶)催化萤光素产生萤光时需要ATP提供能量研制而成。当萤火虫萤光素酶和萤光素都过量时, 在一定的浓度范围内荧光的产生和ATP的浓度成正比, 这样就可以高灵敏地检测溶液中的ATP浓度。本试剂盒在样品体积为100微升时可以检测浓度低至5nM的ATP。而常规的细胞或组织裂解液中ATP的浓度仅为0.1-1 μ M, 一些常见细胞的细胞内ATP水平约为10nmol/mg蛋白。并且本试剂盒的检测浓度范围非常大, 检测上限可以高达10 μ M, 并在10nM-10 μ M范围内可以形成良好的标准曲线。
- 使用本试剂盒中的ATP检测裂解液裂解获得的细胞或组织样品, 不仅可以用于ATP检测, 还可用于检测蛋白浓度、进行SDS-PAGE或一些常规的较易溶解蛋白的Western检测。
- 本试剂盒的使用方便快捷, 通常10-20个样品可以在30-60分钟内测定完毕。
- 碧云天的三款ATP检测试剂盒的主要特点和差异如下:

产品编号	S0026	S0026B	S0027
产品名称	ATP检测试剂盒	ATP检测试剂盒	增强型ATP检测试剂盒
检测灵敏度	++++	+++	+++++
特殊样品兼容性	++++	+++++	++++
信号值	++++	+++	+++++
信号稳定性	+++++	++++	+++++
检测下限(ATP)	1nM	5nM	0.1nM
检测上限(ATP)	10 μ M	10 μ M	10 μ M
线性范围(ATP)	1nM-10 μ M	10nM-10 μ M	0.1nM-10 μ M
推荐指数	++++	+++	+++++

- 国际顶级期刊Cell发表的论文中注明使用了本产品(Cell. 2014 Jul 31;158(3):607-19)。
- 一个包装的本试剂盒至少可以检测200个样品。

包装清单:

产品编号	产品名称	包装
S0026B-1	ATP检测试剂	0.2 ml(50 μ l/管, 共4管)
S0026B-2	ATP检测试剂稀释液	20ml
S0026B-3	ATP标准溶液 (0.5mM)	0.1ml
S0026B-4	ATP检测裂解液	100ml
—	说明书	1份

保存条件:

-20 $^{\circ}$ C保存, 半年有效。-70 $^{\circ}$ C保存, 一年有效。ATP检测试剂需避光保存。

注意事项:

- ATP检测试剂中含有萤光素酶, 反复冻融会导致其逐渐失活。尽管经测试ATP检测试剂反复冻融5次对于其检测效果无明显影响, 为取得较好的使用效果, 建议用户使用时的冻融次数不宜超过3次。ATP检测试剂稀释成ATP检测工作液后, 最好一次用完, 不宜冻存后再使用。
- ATP, 特别是裂解后样品中的ATP在室温不太稳定, 需在4 $^{\circ}$ C或冰上操作。ATP在冰上可以稳定长达6小时。

- 本试剂盒需使用luminometer，即化学发光仪(检测荧光素酶报告基因时所用的仪器)。如果没有luminometer，也可以使用液闪仪。液闪仪的测定效果取决于液闪仪的检测灵敏度和检测精度。
- 使用可检测化学发光的多功能酶标仪时，推荐使用孔和孔之间不透光的96孔白板或黑板。如使用普通的透明96孔板，须特别注意在检测孔之间设置间隔孔，以减少邻近孔之间的相互干扰。对于透明96孔板，一个发光孔可以使上下或左右邻近孔的RLU值升高该孔的10-20%左右，使上下或左右间隔一个孔的邻近孔的RLU值升高该孔的1%-4%左右；对于相同的样品，底部不透光的96孔白板的化学发光读数可以达到透明96孔板的5-10倍左右，达到底部透光96孔白板读数的3倍左右(实测数据会因96孔板、检测仪器和样品等的不同而存在差异)。
- 本试剂盒提供的ATP检测裂解液可以有效裂解并释放常见的培养细胞和组织中的ATP。对于一些特殊的组织或样品，如果发现检测出来的ATP水平显著低于预期水平，可以在裂解样品后并且在离心前，取部分样品煮沸2分钟以充分释放ATP。煮沸后样品中的蛋白会变性，从而会在后续的离心步骤中被沉淀，因此煮沸的样品不能用于蛋白浓度测定、SDS-PAGE和Western检测。可以使用剩余的部分样品进行蛋白浓度测定、SDS-PAGE和Western检测。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明：

1. 样品测定的准备：(注意：样品裂解需在4°C或冰上操作)

对于贴壁细胞：

吸除培养液，按照6孔板每孔加入200微升裂解液的比例(即相当于细胞培养液量2毫升的1/10)加入裂解液，裂解细胞。裂解细胞时为了裂解充分，可以使用移液器进行反复吹打或晃动培养板使裂解液充分接触并裂解细胞。通常细胞在接触裂解液后会立即裂解。裂解后4°C 12000g离心5分钟，取上清，用于后续的测定。

对于悬浮细胞：

用离心管离心沉淀细胞，弃上清，轻轻弹散细胞，按照6孔板每孔的细胞量加入200微升裂解液的比例加入裂解液，裂解细胞。裂解细胞时为了裂解充分可以弹击离心管管底或适当Vortex使裂解液充分接触并裂解细胞。通常细胞在接触裂解液后会立即裂解。裂解后4°C 12000g离心5分钟，取上清，用于后续的测定。

对于组织样品：

按照每20毫克组织加入约100-200微升裂解液的比例加入裂解液，然后用玻璃匀浆器或其它匀浆设备进行匀浆。充分匀浆可以确保组织被完全裂解。裂解后4°C 12000g离心5分钟，取上清，用于后续的测定。

2. 标准曲线测定的准备：

冰浴上溶解待用试剂，把ATP标准溶液用ATP检测裂解液稀释成适当的浓度梯度。具体的浓度需根据样品中ATP的浓度而定。初次检测可以检测0.1、1和10微摩尔/升这几个浓度，在后续的实验中，可以根据样品中ATP的浓度对标准品的浓度进行适当调节。

3. ATP检测工作液的配制：

按照每个样品或标准品需100微升ATP检测工作液的比例配制适当量的ATP检测工作液。把待用试剂在冰浴上溶解。取适量的ATP检测试剂，按照1:100的比例用ATP检测试剂稀释液稀释ATP检测试剂。例如50微升ATP检测试剂可以加入5毫升ATP检测试剂稀释液配制成5毫升ATP检测工作液。稀释后的ATP检测试剂即为用于后续实验的ATP检测工作液。ATP检测工作液可在冰浴上暂时保存。

4. ATP浓度的测定：

- a. 加100微升ATP检测工作液到检测孔或检测管内。室温放置3-5分钟，以使本底性的ATP全部被消耗掉，从而降低本底。可以一次性把10-20个检测孔或检测管分别加上100微升ATP检测工作液，从而节省时间。
- b. 在检测孔或检测管内加上20微升样品或标准品，迅速用枪(微量移液器)混匀，至少间隔2秒后，用luminometer或液闪仪测定RLU值或CPM。(注：样品的体积可以自行在10-100微升范围内调节。如果样品中的ATP浓度比较低则可以加入100微升样品，如果样品中ATP浓度比较高则可以加入较小体积的样品，同时标准品也需要使用相同的体积。如果样品中ATP的浓度特别高，可以用ATP检测裂解液稀释样品后再测定。本试剂盒在加入100微升标准品时，大致在10nmol/L-10 μ mol/L的浓度范围内有很好的线性关系。)
- c. 根据标准曲线计算出样品中ATP的浓度。
- d. 为了消除样品制备时由于蛋白量的差异而造成的误差，可以用碧云天生产的BCA蛋白浓度测定试剂盒(P0009/P0010/P0010S/P0011/P0012/P0012S)测定样品中的蛋白浓度。然后把ATP的浓度换算成nmol/mg蛋白的形式。

常见问题：

1. Luminometer和荧光分光光度计有何不同？

荧光分光光度计检测的样品本身不能发光，样品需要由特定波长的激发光激发，然后才能产生荧光并被荧光分光光度计检测。Luminometer，即化学发光检测仪，检测的样品本身可以发光，不需要激发光进行激发。也就是说Luminometer是检测化学发光(荧光)的仪器。有些型号的荧光分光光度计也具有luminometer的功能，即也可以检测化学发光。您所使用的荧光分光光度计能否用于化学发光的测定请仔细阅读该仪器的说明书。

2. 可以进行荧光素酶报告基因检测的仪器是否就可以用于本试剂盒的ATP检测？

是。荧光素酶报告基因的检测原理和本ATP检测试剂盒的原理相同，可以用相同的仪器测定，并且可以选择相同的测定参数，例如检测前等待时间为2秒，检测时间为10秒等。

3. 多功能酶标仪是否可以用于本试剂盒的ATP检测？

不一定。如果该多功能酶标仪具有检测化学发光的功能，即具有luminometer的功能，就可以用于本试剂盒的检测，否则就不能了。

相关产品：

产品编号	产品名称	包装
S0026	ATP检测试剂盒	200次
S0026B	ATP检测试剂盒	200次
S0027	增强型ATP检测试剂盒	200次

使用本产品的文献：

- Zhu CH, Lu FP, He YN, Han ZL, Du LX. Regulation of avilamycin biosynthesis in *Streptomyces viridochromogenes*: effects of glucose, ammonium ion, and inorganic phosphate. *APPL MICROBIOL BIOT* 2007 Jan;73(5):1031-8.
- Wang YM, Pu P, Le WD. ATP depletion is the major cause of MPP+ induced dopamine neuronal death and worm lethality in alpha-synuclein transgenic *C. elegans*. *Neurosci Bull* 2007 Nov;23(6):329-35.
- Sheng B, Gong K, Niu Y, Liu L, Yan Y, Lu G, Zhang L, Hu M, Zhao N, Zhang X, Tang P, Gong Y. Inhibition of gamma-secretase activity reduces Aβ production, reduces oxidative stress, increases mitochondrial activity and leads to reduced vulnerability to apoptosis: Implications for the treatment of Alzheimer's disease. *FREE RADICAL BIO MED* 2009 May 15;46(10):1362-75.
- Chu L, Wang J, Wang B, Xing XH, Yan S, Sun X, Jurcik B. Changes in biomass activity and characteristics of activated sludge exposed to low ozone dose. *Chemosphere* 2009 Sep;77(2):269-72.
- Liu Q, Zhang Y, Lin Z, Shen H, Chen L, Hu L, Jiang H, Shen X. Danshen extract 15,16-dihydroartemisinin I functions as a potential modulator against metabolic syndrome through multi-target pathways. *J STEROID BIOCHEM* 2010;120(4-5):155-63.
- Yang Q, Gong ZJ, Zhou Y, Yuan JQ, Cheng J, Tian L, Li S, Lin XD, Xu R, Zhu ZR, Mao C. Role of *Drosophila* alkaline ceramidase (Dacer) in *Drosophila* development and longevity. *Cell Mol Life Sci* 2010 May;67(9):1477-90.
- Zhang W, Shi Y, Chen Y, Yu S, Hao J, Luo J, Sha X, Fang X. Enhanced antitumor efficacy by paclitaxel-loaded pluronic P123/F127 mixed micelles against non-small cell lung cancer based on passive tumor targeting and modulation of drug resistance. *Eur J Pharm Biopharm* 2010 Aug;75(3):341-53.
- Zhu SS, Ren Y, Zhang M, Cao JQ, Yang Q, Li XY, Bai H, Jiang L, Jiang Q, He ZG, Chen Q. Wld(S) protects against peripheral neuropathy and retinopathy in an experimental model of diabetes in mice. *Diabetologia* 2011 Sep;54(9):2440-50.
- Zhang W, Shi Y, Chen Y, Ye J, Sha X, Fang X. Multifunctional Pluronic P123/F127 mixed polymeric micelles loaded with paclitaxel for the treatment of multidrug resistant tumors. *Biomaterials* 2011 Apr;32(11):2894-906.
- Cao X, Yang M, Wei RC, Zeng Y, Gu JF, Huang WD, Yang DQ, Li HL, Ding M, Wei N, Zhang KJ, Xu B, Liu XR, Qian QJ, Liu XY. Cancer targeting Gene-Viro-Therapy of liver carcinoma by dual-regulated oncolytic adenovirus armed with TRAIL gene. *Gene Ther* 2011 Aug;18(8):765-77.
- Chang WQ, Wu XZ, Cheng AX, Zhang L, Ji M, Lou HX. Retigeric acid B exerts antifungal effect through enhanced reactive oxygen species and decreased cAMP. *BBA-BIOMEMBRANES* 2011 May;1810(5):569-76.
- Cui Y, Zhao Y, Tian Y, Zhang W, Lü X, Jiang X. The molecular mechanism of action of bactericidal gold nanoparticles on *Escherichia coli*. *Biomaterials* 2012 Mar;33(7):2327-33.
- Feng Q, Yu AF, Chu LB, Chen HZ, Xing XH. Mechanistic Study of On-site sludge Reduction in a Baffled Bioreactor Consisting of three series of Alternating Aerobic and Anaerobic Compartments. *Biochemical Engineering Journal* 2012 Aug;67(15):45-51.
- Chen Y, Yin Q, Ji X, Zhang S, Chen H, Zheng Y, Sun Y, Qu H, Wang Z, Li Y, Wang X, Zhang K, Zhang L, Shi J. Manganese oxide-based multifunctionalized mesoporous silica nanoparticles for pH-responsive MRI, ultrasonography and circumvention of MDR in cancer cells. *Biomaterials* 2012 Oct;33(29):7126-37.
- Peng Y, Dong D, Jiang C, Yu B, Wang X, Ji Y. Relationship between respiration deficiency and azole resistance in clinical *Candida glabrata*. *FEMS Yeast Res* 2012 Sep;12(6):719-27.
- Yu SJ, Liu HC, Ling-Ling E, Wang DS, Zhu GX. Proliferation and differentiation of osteoblasts from the mandible of osteoporotic rats. *EXP BIOL MED* 2012 Apr;237(4):395-406.
- Liu D, Liu S, Jing X, Li X, Li W, Huang Y. Necrosis of cervical carcinoma by dichloroacetate released from electrospun polylactide mats. *Biomaterials* 2012 Jun;33(17):4362-9.
- Li XY, Jing CQ, Zang XY, Yang S, Wang JJ. Toxic cytological alteration and mitochondrial dysfunction in PC12 cells induced by 1-octyl-3-methylimidazolium chloride. *Toxicol In Vitro* 2012 Oct;26(7):1087-92.
- Song GX, Shen YH, Liu YQ, Sun W, Miao LP, Zhou LJ, Liu HL, Yang R, Kong XQ, Cao KJ, Qian LM, Sheng YH. Overexpression of FABP3 promotes apoptosis through inducing mitochondrial impairment in embryonic cancer cells. *J Cell Biochem* 2012 Dec;113(12):3701-8.
- Zhu S, Wang Y, Chen M, Jin J, Qiu Y, Huang M, Huang Z. Protective Effect of Schisandrin B Against Cyclosporine A-Induced Nephrotoxicity In Vitro and In Vivo. *AM J CHINESE MED* 2012;40(3):551-66.
- Yu Ping Li, Hu Su, Xiao Fang Pi, Yan Chun Gong, Xiang Yuan Xiong, Guang Jie Wu, Zi Ling Li. The Stimulatory Activities of Baicalein and Baicalin Compounds Derived from *Scutellaria baicalensis* on Insulin Secretion In Vitro. *ADV MATER* 2012;554-556:1673.
- Huang HY, Cai JF, Liu QC, Hu GH, Xia B, Mao JY, Wu DS, Liu JJ, Zhuang ZX. Role of poly(ADP-ribose) glycohydrolase in the regulation of cell fate in response to benzo(a)pyrene. *Exp Cell Res* 2012 Mar 10;318(5):682-90.
- Yuan Y, Ming Z, Gong-Hua H, Lan G, Lu D, Peng L, Feng J, Cai-Gao Z. Cr(VI) induces the decrease of ATP level and the increase of apoptosis rate mediated by ROS or VDAC1 in L-02 hepatocytes. *ENVIRON TOXICOL CHEM* 2012 Jul 7;34(2):579-587.
- Liu YQ, Cheng X, Guo LX, Mao C, Chen YJ, Liu HX, Xiao QC, Jiang S, Yao ZJ, Zhou GB. Identification of an annonaceous acetogenin mimetic, AA005, as an AMPK activator and autophagy inducer in colon cancer cells. *PLoS One* 2012;7(10):e47049.
- Yue R, Hu H, Yiu KH, Luo T, Zhou Z, Xu L, Zhang S, Li K, Yu Z. Lycopene protects against hypoxia/reoxygenation-induced apoptosis by preventing mitochondrial dysfunction in primary neonatal mouse cardiomyocytes. *PLoS One* 2012;7(11):e50778.
- Peng Y, Yang J, Zhang E, Sun H, Wang Q, Wang T, Su Y, Shi C. Human positive coactivator 4 is a potential novel therapeutic target in non-small cell lung cancer. *Cancer Gene Ther* 2012 Oct;19(10):690-6.
- Wei Z, Yuan S, Hao J, Fang X. Mechanism of inhibition of P-glycoprotein mediated efflux by Pluronic P123/F127 block copolymers: Relationship between copolymer concentration and inhibitory activity. *Eur J Pharm Biopharm* 2013 Feb;83(2):266-74.
- Ma C, Tang Z, Wang K, Yang X, Tan W. A novel sensitive and selective ligation-based ATP assay using a molecular beacon. *Analyst* 2013 Apr 22;138(10):3013-7.
- Tang Y, Wang F, Jin C, Liang H, Zhong X, Yang Y. Mitochondrial injury induced by nanosized titanium dioxide in A549 cells and rats. *ENVIRON TOXICOL CHEM* 2013 Jul;36(1):66-72.

30. Wang X, Bai H, Zhang X, Liu J, Cao P, Liao N, Zhang W, Wang Z, Hai C. Inhibitory effect of oleanolic acid on hepatocellular carcinoma via ERK-p53-mediated cell cycle arrest and mitochondrial-dependent apoptosis. *Carcinogenesis* 2013 Jun;34(6):1323-30.
31. Qiu L, Chen Y, Gao M, Zheng C, Zhao Q. Phagocytic uptake and ROS-mediated cytotoxicity in human hepatic cell line of polyphosphazene nanoparticles. *J Biomed Mater Res A* 2013 Jan;101(1):285-97.
32. Wang K, Gao F, Ji Y, Liu Y, Dan Z, Yang P, Zhu Y, Li S. ORFH79 impairs mitochondrial function via interaction with a subunit of electron transport chain complex III in Honglian cytoplasmic male sterile rice. *New Phytol* 2013 Apr;198(2):408-18.
33. Wang J, Chen Y, Zhang W, Zheng G, Meng S, Che H, Ke T, Yang J, Chen J, Luo W. Akt Activation Protects Liver Cells from Apoptosis in Rats during Acute Cold Exposure. *INT J BIOL SCI* 2013; 9(5):509-17.
34. Li N, Zhang CX, Wang XX, Zhang L, Ma X, Zhou J, Ju RJ, Li XY, Zhao WY, Lu WL. Development of targeting lonidamine liposomes that circumvent drug-resistant cancer by acting on mitochondrial signaling pathways. *Biomaterials* 2013 Apr;34(13):3366-80.
35. Dong LY, Jin J, Lu G, Kang XL. Astaxanthin Attenuates the Apoptosis of Retinal Ganglion Cells in db/db Mice by Inhibition of Oxidative Stress. *Mar Drugs* 2013 Mar 21;11(3):960-74.
36. Qi M, Fan S, Yao G, Li Z, Zhou H, Tashiro S, Onodera S, Xia M, Ikejima T. Pseudolaric acid B-induced autophagy contributes to senescence via enhancement of ROS generation and mitochondrial dysfunction in murine fibrosarcoma L929 cells. *J Pharmacol Sci* 2013;121(3):200-11.
37. Yin Q, Shen J, Zhang Z, Yu H, Chen L, Gu W, Li Y. Multifunctional Nanoparticles Improve Therapeutic Effect for Breast Cancer by Simultaneously Antagonizing Multiple Mechanisms of Multidrug Resistance. *Biomacromolecules* 2013 Jul 8;14(7):2242-52.
38. Zhang X, Zuo X, Yang B, Li Z, Xue Y, Zhou Y, Huang J, Zhao X, Zhou J, Yan Y, Zhang H, Guo P, Sun H, Guo L, Zhang Y, Fu XD. MicroRNA directly enhances mitochondrial translation during muscle differentiation. *Cell* 2014 Jul 31;158(3):607-19.
39. Cui J, He W, Yi B, Zhao H, Lu K, Ruan H, Ma D. mTOR pathway is involved in ADP-evoked astrocyte activation and ATP release in the spinal dorsal horn in a rat neuropathic pain model. *Neuroscience* 2014 Sep 5;275:395-403.
40. Luo DJ, Feng Q, Wang ZH, Sun DS, Wang Q, Wang JZ, Liu GP. Knockdown of phosphotyrosyl phosphatase activator induces apoptosis via mitochondrial pathway and the attenuation by simultaneous tau hyperphosphorylation. *J Neurochem* 2014 Sep;130(6):816-25.
41. Meng G, Xia M, Xu C, Yuan D, Schnurr M, Wei J. Multifunctional antitumor molecule 5'-triphosphate siRNA combining glutaminase silencing and RIG-I activation. *Int J Cancer* 2014 Apr 15;134(8):1958-71.
42. Sun S, Han Y, Liu J, Fang Y, Tian Y, Zhou J, Ma D, Wu P. Trichostatin A targets the mitochondrial respiratory chain, increasing mitochondrial reactive oxygen species production to trigger apoptosis in human breast cancer cells. *PLoS One* 2014 Mar 13;9(3):e91610.
43. Xie Y, Li Q, Yang Q, Yang M, Zhang Z, Zhu L, Yan H, Feng R, Zhang S, Huang C, Liu Z, Wen T. Overexpression of DCF1 inhibits glioma through destruction of mitochondria and activation of apoptosis pathway. *SCI REP-UK* 2014 Jan 15;4:3702.
44. Su B, Ji YS, Sun XL, Liu XH, Chen ZY. Brain-derived neurotrophic factor (BDNF)-induced mitochondrial motility arrest and presynaptic docking contribute to BDNF-enhanced synaptic transmission. *J Biol Chem* 2014 Jan 17;289(3):1213-26.
45. Jia WC, Liu G, Zhang CD, Zhang SP. Formononetin attenuates hydrogen peroxide (H₂O₂)-induced apoptosis and NF- κ B activation in RGC-5 cells. *EUR REV MED PHARMACO* 2014;18(15):2191-7.
46. Zhang H, Sun R, Liu XY, Shi XM, Wang WF, Yu LG, Guo XL. A tetramethylpyrazine piperazine derivative CXK137 prevents cell injury in SH-SY5Y cells and improves memory dysfunction of rats with vascular Dementia. *Neurochem Res* 2014 Feb;39(2):276-86.
47. Xu P, Yu H, Zhang Z, Meng Q, Sun H, Chen X, Yin Q, Li Y. Hydrogen-bonded and reduction-responsive micelles loading atorvastatin for therapy of breast cancer metastasis. *Biomaterials* 2014 Aug;35(26):7574-87.
48. Chen G, Dai J, Tan S, Meng S, Liu Z, Li M, Cui Q, Yu M. MTERF1 regulates the oxidative phosphorylation activity and cell proliferation in HeLa cells. *ACTA BIOCH BIOPH SIN* 2014 Jun;46(6):512-21.
49. Li Q, Zhang T, Wang J, Zhang Z, Zhai Y, Yang GY, Sun X. Rapamycin attenuates mitochondrial dysfunction via activation of mitophagy in experimental ischemic stroke. *BIOCHEM BIOPH RES CO* 2014 Feb 7;444(2):182-8.
50. He Z, Liu J, Du L. The unexpected effect of PEGylated gold nanoparticles on the primary function of erythrocytes. *Nanoscale* 2014 Aug 7;6(15):9017-24.
51. Xia M, Meng G, Jiang A, Chen A, Dahlhaus M, Gonzalez P, Beltinger C, Wei J. Mitophagy switches cell death from apoptosis to necrosis in NSCLC cells treated with oncolytic measles virus. *ONCOTARGET* 2014 Jun 15;5(11):3907-18.
52. Shen Y, Yang J, Li J, Shi X, Ouyang L, Tian Y, Lu J. Carnosine inhibits the proliferation of human gastric cancer SGC-7901 cells through both of the mitochondrial respiration and glycolysis pathways. *PLoS One* 2014 Aug 12;9(8):e104632.
53. Jiang XR, Yu XY, Fan JH, Guo L, Zhu C, Jiang W, Lu SH. RFT2 is overexpressed in esophageal squamous cell carcinoma and promotes tumorigenesis by sustaining cell proliferation and protecting against cell death. *Cancer Lett* 2014 Oct 10;353(1):78-86.
54. Li J, Ke W, Zhou Q, Wu Y, Luo H, Zhou H, Yang B, Guo Y, Zheng Q, Zhang Y. Tumour necrosis factor- α promotes liver ischaemia-reperfusion injury through the PGC-1 α /Mfn2 pathway. *J Cell Mol Med* 2014 Sep;18(9):1863-73.
55. Li X, Xu Z, Jiang Z, Sun L, Ji J, Miao J, Zhang X, Li X, Huang S, Wang T, Zhang L. Hypoglycemic effect of catalpol on high-fat diet/streptozotocin-induced diabetic mice by increasing skeletal muscle mitochondrial biogenesis. *ACTA BIOCH BIOPH SIN* 2014 Sep;46(9):738-48.
56. Liu B, Tang M, Han Z, Li J, Zhang J, Lu P, Song N, Wang Z, Yin C, Zhang W. Co-Incubation of Human Spermatozoa with Anti-VDAC Antibody Reduced Sperm Motility. *CELL PHYSIOL BIOCHEM* 2014;33(1):142-50.
57. Yan H, Zhang D, Hao S, Li K, Hang CH. Role of Mitochondrial Calcium Uniporter in Early Brain Injury After Experimental Subarachnoid Hemorrhage. *Mol Neurobiol* 2015 Dec;52(3):1637-47.
58. Du LL, Chai DM, Zhao LN, Li XH, Zhang FC, Zhang HB, Liu LB, Wu K, Liu R, Wang JZ, Zhou AMPK Activation Ameliorates Alzheimer's Disease-Like Pathology and Spatial Memory Impairment in a Streptozotocin-Induced Alzheimer's Disease Model in Rats. *J Alzheimers Dis* 2015;43(3):775-84.
59. Chen Y, Zhang Z, Yang K, Du J, Xu Y, Liu S. Myeloid zinc-finger 1 (MZF-1) suppresses prostate tumor growth through enforcing ferroportin-conducted iron egress. *Oncogene* 2015 Jul;34(29):3839-47.
60. Lu L, Han H, Tian Y, Li W, Zhang J, Feng M, Li Y. Aurora kinase A mediates c-Myc's oncogenic effects in hepatocellular carcinoma. *MOL CARCINOGEN* 2015 Nov;54(11):1467-79.
61. Ji X, Yi B, Xu Y, Zhao Y, Zhong H, Ding C. A novel fluorescent biosensor for Adenosine Triphosphate detection based on the polydopamine nanospheres integrating with enzymatic recycling amplification. *Talanta* 2017 Jul 1;169:8-12.
62. Du X, Fu X, Yao K, Lan Z, Xu H, Cui Q, Yang E. Bcl-2 delays cell cycle through mitochondrial ATP and ROS. *Cell Cycle* 2017 Apr 3;16(7):707-713.
63. Shen J, Xu S, Zhou H, Liu H, Jiang W, Hao J, Hu Z. IL-1 β induces apoptosis and autophagy via mitochondria pathway in human degenerative nucleus pulposus cells. *SCI REP-UK* 2017 Jan 25;7:41067.
64. Qiu M, Ke L, Zhang S, Zeng X, Fang Z, Liu J. JS-K, a GST-activated nitric oxide donor prodrug, enhances chemo-sensitivity in renal carcinoma cells and prevents cardiac myocytes toxicity induced by Doxorubicin. *CANCER CHEMOTH PHARM* 2017 Aug;80(2):275-286.

65. Li J, Song P, Zhu L, Aziz N, Zhou Q, Zhang Y, Xu W, Feng L, Chen D, Wang X, Jin H. Synthetic lethality of glutaminolysis inhibition, autophagy inactivation and asparagine depletion in colon cancer. *ONCOTARGET* 2017 Jun 27;8(26):42664-42672.
66. Liao X, Song L, Zhang L, Wang H, Tong Q, Xu J, Yang G, Yang S, Zheng H. LAMP3 regulates hepatic lipid metabolism through activating PI3K/Akt pathway. *Mol Cell Endocrinol* 2017 Oct 19. pii: S0303-7207(17)30547-6.
67. He H, Lai Y, Hao Y, Liu Y, Zhang Z, Liu X, Guo C, Zhang M, Zhou H, Wang N, Luo XG, Huo L, Ma W, Zhang TC. Selective p300 inhibitor C646 inhibited HPV E6-E7 genes, altered glucose metabolism and induced apoptosis in cervical cancer cells. *Eur J Pharmacol* 2017 Oct 5;812:206-215.
68. Zhang J, Cai Q, Jiang M, Liu Y, Gu H, Guo J, Sun H, Fang J, Jin L. Mesencephalic astrocyte-derived neurotrophic factor alleviated 6-OHDA-induced cell damage via ROS-AMPK/mTOR mediated autophagic inhibition. *Exp Gerontol* 2017 Mar;89:45-56.
69. Hu LT, Zhu BL, Lai YJ, Long Y, Zha JS, Hu XT, Zhang JH, Chen GJ. HMGCS2 promotes autophagic degradation of the amyloid- β precursor protein through ketonebody-mediated mechanisms. *BIOCHEM BIOPH RES CO* 2017 Apr 29;486(2):492-498.
70. Ju L, Tong W, Qiu M, Shen W, Sun J, Zheng S, Chen Y, Liu W, Tian J. Antioxidant MMCC ameliorates catch-up growth related metabolic dysfunction. *ONCOTARGET* 2017 Oct 23;8(59):99931-99939.
71. Pan T, Zhang M, Zhang F, Yan G, Ru Y, Wang Q, Zhang Y, Wei X, Xu X, Shen L, Zhang J, Wu K, Yao L, Li X. NDRG2 overexpression suppresses hepatoma cells survival during metabolic stress through disturbing the activation of fatty acid oxidation. *BIOCHEM BIOPH RES CO* 2017 Feb 5;483(2):860-866.
72. Bai F, Guo F, Jiang T, Wei H, Zhou H, Yin H, Zhong H, Xiong L, Wang Q. Arachidonyl-2-Chloroethylamide Alleviates Cerebral Ischemia Injury Through Glycogen Synthase Kinase-3 β -Mediated Mitochondrial Biogenesis and Functional Improvement. *Mol Neurobiol* 2017 Mar;54(2):1240-1253.
73. Zhu P, Xue J, Zhang ZJ, Jia YP, Tong YN, Han D, Li Q, Xiang Y, Mao XH, Tang B. Helicobacter pylori VacA induces autophagic cell death in gastric epithelial cells via the endoplasmic reticulum stress pathway. *Cell Death Dis* 2017 Dec 13;8(12):3207.
74. Chen E, Chen X, Yuan X, Wei S, Zhou L, Zhou J, Shen J. One-pot method to prepare a theranostic nanosystem with magnetic resonance imaging function and anticancer activity through multiple mechanisms. *DALTON T* 2017 Apr 19;46(16):5151-5158.
75. Sun R, Zhang J, Wei H, Meng X, Ding Q, Sun F, Cao M, Yin L, Pu Y. Acetyl-l-carnitine partially prevents benzene-induced hematotoxicity and oxidative stress in C3H/He mice. *ENVIRON TOXICOL CHEM* 2017 Apr;36(1):108-113.
76. Yin K, Zhu R, Wang S, Zhao RC. Low-Level Laser Effect on Proliferation, Migration, and Antiapoptosis of Mesenchymal Stem Cells. *Stem Cells Dev* 2017 May 15;26(10):762-775.
77. Huang Q, Sun M, Li M, Zhang D, Han F, Wu JC, Fukunaga K, Chen Z, Qin ZH. Combination of NAD⁺ and NADPH Offers Greater Neuroprotection in Ischemic Stroke Models by Relieving Metabolic Stress. *Mol Neurobiol* 2017 Nov 21.
78. He W, Cui L, Zhang C, Zhang X, He J, Xie Y, Chen Y. Sonic hedgehog promotes neurite outgrowth of cortical neurons under oxidative stress: Involving of mitochondria and energy metabolism. *Exp Cell Res* 2017 Jan 1;350(1):83-90.
79. Xie K, Zhu M, Xiang P, Chen X, Kasimimali A, Lu R, Wang Q, Mou S, Ni Z, Gu L, Pang H. Protein Kinase A/CREB Signaling Prevents Adriamycin-Induced Podocyte Apoptosis via Upregulation of Mitochondrial Respiratory Chain Complexes. *Mol Cell Biol* 2017 Dec 13;38(1). pii: e00181-17.
80. Yang X, Wang H, Huang C, He X, Xu W, Luo Y, Huang K. Zinc enhances the cellular energy supply to improve cell motility and restore impaired energetic metabolism in a toxic environment induced by OTA. *SCI REP-UK* 2017 Nov 7;7(1):14669.
81. Khan MR, Xiang S, Song Z, Wu M. The p53-inducible long noncoding RNA TRINGS protects cancer cells from necrosis under glucose starvation. *EMBO J* 2017 Dec 1;36(23):3483-3500.
82. Xiao N, Yang LL, Yang YL, Liu LW, Li J, Liu B, Liu K, Qi LW, Li P. Ginsenoside Rg5 Inhibits Succinate-Associated Lipolysis in Adipose Tissue and Prevents Muscle Insulin Resistance. *Front Pharmacol* 2017 Feb 14;8:43.
83. Yao GD, Yang J, Li XX, Song XY, Hayashi T, Tashiro SI, Onodera S, Song SJ, Ikejima T. Blocking the utilization of glucose induces the switch from senescence to apoptosis in pseudolaric acid B-treated human lung cancer cells in vitro. *Acta Pharmacol Sin* 2017 Oct;38(10):1401-1411.
84. Wang Y, Gao P, Wei C, Li H, Zhang L, Zhao Y, Wu B, Tian Y, Zhang W, Wu L, Wang R, Xu C. Calcium sensing receptor protects high glucose-induced energy metabolism disorder via blocking gp78-ubiquitin proteasome pathway. *Cell Death Dis* 2017 May 18;8(5):e2799.
85. Zhang JJ, Liu WQ, Peng JJ, Ma QL, Peng J, Luo XJ. miR-21-5p/203a-3p promote ox-LDL-induced endothelial cell senescence through down-regulation of mitochondrial fission protein Drp1. *Mech Ageing Dev* 2017 Jun;164:8-19.
86. Yang X, Liu T, Chen B, Wang F, Yang Q, Chen X. Grape seed proanthocyanidins prevent irradiation-induced differentiation of human lung fibroblasts by ameliorating mitochondrial dysfunction. *SCI REP-UK* 2017 Mar 3;7(1):62.
87. Wang C, Qi S, Liu C, Yang A, Fu W, Quan C, Duan P, Yu T, Yang K. Mitochondrial Dysfunction and Ca²⁺ Overload in Injured Sertoli Cells Exposed to Bisphenol A. *Environ Toxicol* 2017 Mar;32(3):823-831.
88. Luo C, Li Y, Guo L, Zhang F, Liu H, Zhang J, Zheng J, Zhang J, Guo S. Graphene Quantum Dots Downregulate Multiple Multidrug-Resistant Genes via Interacting with Their C-Rich Promoters. *Adv Healthc Mater* 2017 Nov;6(21).
89. Tian X, He W, Yang R, Liu Y. Di-3-n-butylphthalide protects the heart against ischemic injury and H9c2 cardiomyoblasts against oxidative stress: involvement of mitochondrial function and biogenesis. *J Biomed Sci* 2017 Jun 15;24(1):38.
90. Jiang Y, Wang X, Hu D. Furanodienone induces G0/G1 arrest and causes apoptosis via the ROS/MAPKs-mediated caspase-dependent pathway in human colorectal cancer cells: a study in vitro and in vivo. *Cell Death Dis* 2017 May 25;8(5):e2815.
91. Zhang L, He D, Li K, Liu H, Wang B, Zheng L, Li J. Emodin targets mitochondrial cyclophilin D to induce apoptosis in HepG2 cells. *Biomed Pharmacother* 2017 Jun;90:222-228.
92. Zhao Y, Yang L, He J, Yang H. STYK1 promotes Warburg effect through PI3K/AKT signaling and predicts a poor prognosis in nasopharyngeal carcinoma. *TUMOR BIOL* 2017 Jul;39(7):1010428317711644.
93. Liu Q, Wang Q, Xu C, Shao W, Zhang C, Liu H, Jiang Z, Gu A. Organochloride pesticides impaired mitochondrial function in hepatocytes and aggravated disorders of fatty acid metabolism. *SCI REP-UK* 2017 Apr 11;7:46339.
94. Li H, Han W, Wang H, Ding F, Xiao L, Shi R, Ai L, Huang Z. Tanshinone IIA Inhibits Glutamate-Induced Oxidative Toxicity through Prevention of Mitochondrial Dysfunction and Suppression of MAPK Activation in SH-SY5Y Human Neuroblastoma Cells. *Oxid Med Cell Longev* 2017;2017:4517486.
95. Li J, Yang YL, Li LZ, Zhang L, Liu Q, Liu K, Li P, Liu B, Qi LW. Succinate accumulation impairs cardiac pyruvate dehydrogenase activity through GRP91-dependent and independent signaling pathways: Therapeutic effects of ginsenoside Rb1. *BBA-BIOMEMBRANES* 2017 Nov;1863(11):2835-2847.
96. Wu LY, Ye ZN, Zhou CH, Wang CX, Xie GB, Zhang XS, Gao YY, Zhang ZH, Zhou ML, Zhuang Z, Liu JP, Hang CH, Shi JX. Roles of Pannexin-1 Channels in Inflammatory Response through the TLRs/NF- κ B Signaling Pathway Following Experimental Subarachnoid Hemorrhage in Rats. *Front Mol Neurosci* 2017 Jun 6;10:175.
97. Xu Y, Wang Y, Wang G, Ye X, Zhang J, Cao G, Zhao Y, Gao Z, Zhang Y, Yu B, Kou J. YiQiFuMai Powder Injection Protects against Ischemic Stroke via Inhibiting Neuronal Apoptosis and PKC δ /Drp1-Mediated Excessive Mitochondrial Fission. *Oxid Med Cell Longev* 2017;2017:1832093.
98. Su Y, Yu QH, Wang XY, Yu LP, Wang ZF, Cao YC, Li JD. JMJD2A promotes the

- Warburg effect and nasopharyngeal carcinoma progression by transactivating LDHA expression. *BMC Cancer* 2017 Jul 11;17(1):477.
99. Xue C, Wang C, Sun Y, Meng Q, Liu Z, Huo X, Sun P, Sun H, Ma X, Ma X, Peng J, Liu K. Targeting P-glycoprotein function, p53 and energy metabolism: Combination of metformin and 2-deoxyglucose reverses the multidrug resistance of MCF-7/Dox cells to doxorubicin. *ONCOTARGET* 2017 Jan 31;8(5):8622-8632.
 100. Zhou W, Tian D, He J, Zhang L, Tang X, Zhang L, Wang Y, Li L, Zhao J, Yuan X, Peng S. Exposure scenario: Another important factor determining the toxic effects of PM2.5 and possible mechanisms involved. *Environ Pollut* 2017 Jul;226:412-425.
 101. Qian G, Liu D, Hu J, Gan F, Hou L, Chen X, Huang K. Ochratoxin A-induced autophagy in vitro and in vivo promotes porcine circovirus type 2 replication. *Cell Death Dis* 2017 Jun 29;8(6):e2909.
 102. Zhou H, Hu S, Jin Q, Shi C, Zhang Y, Zhu P, Ma Q, Tian F, Chen Y. Mff-Dependent Mitochondrial Fission Contributes to the Pathogenesis of Cardiac Microvasculature Ischemia/Reperfusion Injury via Induction of mROS-Mediated Cardiolipin Oxidation and HK2/VDAC1 Disassociation-Involved mPTP Opening. *J Am Heart Assoc* 2017 Mar 13;6(3).
 103. Wang Y, Hu L, Zhang X, Zhao H, Xu H, Wei Y, Jiang H, Xie C, Zhou Y, Zhou F. Downregulation of Mitochondrial Single Stranded DNA Binding Protein (SSBP1) Induces Mitochondrial Dysfunction and Increases the Radiosensitivity in Non-Small Cell Lung Cancer Cells. *J Cancer* 2017 May 12;8(8):1400-1409.
 104. Xue Q, Pei H, Liu Q, Zhao M, Sun J, Gao E, Ma X, Tao L. MICU1 protects against myocardial ischemia/reperfusion injury and its control by the importer receptor Tom70. *Cell Death Dis* 2017 Jul 13;8(7):e2923.
 105. Sun Y, Zong L, Gao Z, Zhu S, Tong J, Cao Y. Mitochondrial DNA damage and oxidative damage in HL-60 cells exposed to 900 MHz radiofrequency fields. *MUTAT RES-REV MUTAT* 2017 Mar;797-799:7-14.
 106. Shan P, Fan G, Sun L, Liu J, Wang W, Hu C, Zhang X, Zhai Q, Song X, Cao L, Cui Y, Zhang S, Wang C. SIRT1 Functions as a Negative Regulator of Eukaryotic Poly(A)RNA Transport. *Curr Biol* 2017 Aug 7;27(15):2271-2284.e5.
 107. Wu D, Wu J, Liu H, Yu M, Tao L, Dong S, Tong X. Role of Pannexin1 channels in the resistance of I-10 testicular cancer cells to cisplatin mediated by ATP/IP3 pathway. *Biomed Pharmacother* 2017 Oct;94:514-522.
 108. Huang XT, Li C, Peng XP, Guo J, Yue SJ, Liu W, Zhao FY, Han JZ, Huang YH, Yang-Li, Cheng QM, Zhou ZG, Chen C, Feng DD, Luo ZQ. An excessive increase in glutamate contributes to glucose-toxicity in β -cells via activation of pancreatic NMDA receptors in rodent diabetes. *SCI REP-UK* 2017 Mar 17;7:44120.
 109. Qian X, Xu W, Xu J, Shi Q, Li J, Weng Y, Jiang Z, Feng L, Wang X, Zhou J, Jin H. Enolase 1 stimulates glycolysis to promote chemoresistance in gastric cancer. *ONCOTARGET* 2017 Jul 18;8(29):47691-47708.
 110. Huang HJ, Liu CW, Zhou X, Zhang CX, Bao YY. A mitochondrial membrane protein is a target for rice ragged stunt virus in its insect vector. *Virus Res* 2017 Feb 2;229:48-56.
 111. Zhang X, Du L, Zhang W, Yang Y, Zhou Q, Du G. Therapeutic effects of baicalein on rotenone-induced Parkinson's disease through protecting mitochondrial function and biogenesis. *SCI REP-UK* 2017 Aug 30;7(1):9968.
 112. Zou L, Wang D, Hu Y, Fu C, Li W, Dai L, Yang L, Zhang J. Drug resistance reversal in ovarian cancer cells of paclitaxel and borneol combination therapy mediated by PEG-PAMAM nanoparticles. *ONCOTARGET* 2017 Jul 31;8(36):60453-60468.
 113. Guo XD, Sun GL, Zhou TT, Wang YY, Xu X, Shi XF, Zhu ZY, Rukachaisirikul V, Hu LH, Shen X. LX2343 alleviates cognitive impairments in AD model rats by inhibiting oxidative stress-induced neuronal apoptosis and tauopathy. *Acta Pharmacol Sin* 2017 Aug;38(8):1104-1119.
 114. Huang XT, Yue SJ, Li C, Huang YH, Cheng QM, Li XH, Hao CX, Wang LZ, Xu JP, Ji M, Chen C, Feng DD, Luo ZQ. A Sustained Activation of Pancreatic NMDARs Is a Novel Factor of β -Cell Apoptosis and Dysfunction. *Endocrinology* 2017 Nov 1;158(11):3900-3913.
 115. Li HJ, Sun XM, Li ZK, Yin QW, Pang H, Pan JJ, Li X, Chen W. LncRNA UCA1 Promotes Mitochondrial Function of Bladder Cancer via the MiR-195/ARL2 Signaling Pathway. *CELL PHYSIOL BIOCHEM* 2017;43(6):2548-2561.
 116. Sun GL, Li Z, Wang WZ, Chen Z, Zhang L, Li Q, Wei S, Li BW, Xu JH, Chen L, He ZY, Ying K, Zhang X, Xu H, Zhang DC, Xu ZK. miR-324-3p promotes gastric cancer development by activating Smad4-mediated Wnt/ β -catenin signaling pathway. *J Gastroenterol* 2017 Nov 4.
 117. Luo X, Dan Wang, Luo X, Zhu X, Wang G, Ning Z, Li Y, Ma X, Yang R, Jin S, Huang Y, Meng Y, Li X. Caveolin 1-related autophagy initiated by aldosterone-induced oxidation promotes liver sinusoidal endothelial cells defenestration. *Redox Biol* 2017 Oct;13:508-521.
 118. Pei HF, Hou JN, Wei FP, Xue Q, Zhang F, Peng CF, Yang Y, Tian Y, Feng J, Du J, He L, Li XC, Gao EH, Li, Yang YJ. Melatonin attenuates postmyocardial infarction injury via increasing Tom70 expression. *J Pineal Res* 2017 Jan;62(1).
 119. Yan X, Zhang G, Bie F, Lv Y, Ma Y, Ma M, Wang Y, Hao X, Yuan N, Jiang X. Eugenol inhibits oxidative phosphorylation and fatty acid oxidation via downregulation of c-Myc/PGC-1 β /ERR α signaling pathway in MCF10A-ras cells. *SCI REP-UK* 2017 Oct 10;7(1):12920.
 120. Liu Y, Han M, Li X, Wang H, Ma M, Zhang S, Guo Y, Wang S, Wang Y, Duan N, Xu B, Yin J, Yao Y. Age-related changes in the mitochondria of human mural granulosa cells. *Hum Reprod* 2017 Dec 1;32(12):2465-2473.
 121. Cao G, Wang Q, Huang W, Tong J, Ye D, He Y, Liu Z, Tang X, Cheng H, Wen Q, Li D, Chau HT, Wen Y, Zhong H, Meng Z, Liu H, Wu Z, Zhao L, Flavell RA, Zhou H, Xu A, Yang H, Yin Z. Long-term consumption of caffeine-free high sucrose cola beverages aggravates the pathogenesis of EAE in mice. *Cell Discov* 2017 Jun 20;3:17020.
 122. Liu X, Yang S, Wang F, Dai X, Yang Y, Bai Z. Comparative analysis of the *Corynebacterium glutamicum* transcriptome in response to changes in dissolved oxygen levels. *J IND MICROBIOL BIOT* 2017 Feb;44(2):181-195.
 123. Wang B, Xiong S, Lin S, Xia W, Li Q, Zhao Z, Wei X, Lu Z, Wei X, Gao P, Liu D, Zhu Z. Enhanced Mitochondrial Transient Receptor Potential Channel, Canonical Type 3-Mediated Calcium Handling in the Vasculature From Hypertensive Rats. *J Am Heart Assoc* 2017 Jul 15;6(7). pii: e005812.
 124. Feng J, Gu Y, Quan Y, Gao W, Dang Y, Cao M, Lu X, Wang Y, Song C, Wang S. Construction of energy-conserving sucrose utilization pathways for improving poly- γ -glutamic acid production in *Bacillus amyloliquefaciens*. *Microb Cell Fact* 2017 Jun 6;16(1):98.
 125. Souza CF, Baldissera MD, Bottari NB, Moreira KLS, da Rocha MIUM, da Veiga ML, Santos RCV, Baldisserotto B. Purinergic signaling modulates the cerebral inflammatory response in experimentally infected fish with *Streptococcus agalactiae*: an attempt to improve the immune response. *Mol Cell Biochem* 2017 Oct 31.
 126. Tang G, Zhang C, Ju Z, Zheng S, Wen Z, Xu S, Chen Y, Ma Z. The mitochondrial membrane protein FgLetm1 regulates mitochondrial integrity, production of endogenous reactive oxygen species and mycotoxin biosynthesis in *Fusarium graminearum*. *Mol Plant Pathol* 2017 Oct 27.
 127. Zou X, Li P, Lou J, Zhang H. Surface coating-modulated toxic responses to silver nanoparticles in *Wolffia globosa*. *Aquat Toxicol* 2017 Aug;189:150-158.
 128. Zhang S, Lin X, Li G, Shen X, Niu D, Lu G, Fu X, Chen Y, Cui M, Bai Y. Knockout of *Eva1a* leads to rapid development of heart failure by impairing autophagy. *Cell Death Dis* 2017 Feb 2;8(2):e2586.
 129. Wang H, Zhang L, Guo X, Bai Y, Li YX, Sha J, Peng C, Wang YL, Liu M. MiR-195 modulates oxidative stress-induced apoptosis and mitochondrial energy production in human trophoblasts via flavin adenine dinucleotide-dependent oxidoreductase domain-containing protein 1 and pyruvate dehydrogenase phosphatase regulatory subunit. *J Hypertens* 2017 Aug 30.
 130. Li H, He F, Zhao X, Zhang Y, Chu X, Hua C, Qu Y, Duan Y, Ming L. YAP Inhibits the Apoptosis and Migration of Human Rectal Cancer Cells via Suppression of

- JNK-Drp1-Mitochondrial Fission-HtrA2/Omi Pathways. *CELL PHYSIOL BIOCHEM* 2017;44(5):2073-2089.
131. Shi Y, Qin W, Nie F, Wen H, Lu K, Cui J. Ulinastatin attenuates neuropathic pain via the ATP/P2Y2 receptor pathway in rat models. *Gene* 2017 Sep 5;627:263-270.
132. Chen W, Hu H, Zhang C, Huang F, Zhang D, Zhang H. Adaptation response of *Pseudomonas fragi* on refrigerated solid matrix to a moderate electric field. *BMC Microbiol* 2017 Feb 10;17(1):32.
133. Castaneda JM, Hua R, Miyata H, Oji A, Guo Y, Cheng Y, Zhou T, Guo X, Cui Y, Shen B, Wang Z, Hu Z, Zhou Z, Sha J, Prunskaitė-Hyyryläinen R, Yu Z, Ramirez-Solis R, Ikawa M, Matzuk MM, Liu M. TCTE1 is a conserved component of the dynein regulatory complex and is required for motility and metabolism in mouse spermatozoa. *P NATL ACAD SCI USA* 2017 Jul 3;114(27):E5370-E5378.
134. Xu JD, Jiang HS, Wei TD, Zhang KY, Wang XW, Zhao XF, Wang JX. Interaction of the Small GTPase Cdc42 with Arginine Kinase Restricts White Spot Syndrome Virus in Shrimp. *J Virol* 2017 Feb 14;91(5). pii: e01916-16.
135. Wang B, Wang H, Xiong J, Zhou Q, Wu H, Xia L, Li L, Yu Z. A Proteomic Analysis Provides Novel Insights into the Stress Responses of *Caenorhabditis elegans* towards Nematicidal Cry6A Toxin from *Bacillus thuringiensis*. *SCI REP-UK* 2017 Oct 26;7(1):14170.
136. Shen Z, Zheng Y, Wu J, Chen Y, Wu X, Zhou Y, Yuan Y, Lu S, Jiang L, Qin Z, Chen Z, Hu W, Zhang X. PARK2-dependent mitophagy induced by acidic postconditioning protects against focal cerebral ischemia and extends the reperfusion window. *Autophagy* 2017 Mar 4;13(3):473-485.
137. Wang C, Lv H, Yang W, Li T, Fang T, Lv G, Han Q, Dong L, Jiang T, Jiang B, Yang G, Wang H. SVCT-2 determines the sensitivity to ascorbate-induced cell death in cholangiocarcinoma cell lines and patient derived xenografts. *Cancer Lett* 2017 Jul 10;398:1-11.
138. Xiao J, Zhang R, Huang F, Liu L, Deng Y, Ma Y, Wei Z, Tang X, Zhang Y, Zhang M. Lychee (*Litchi chinensis* Sonn.) Pulp Phenolic Extract Confers a Protective Activity against Alcoholic Liver Disease in Mice by Alleviating Mitochondrial Dysfunction. *J AGR FOOD CHEM* 2017 Jun 21;65(24):5000-5009.
139. Qiu M, Shi F, Dai F, Song R, Wang S, You Y, Zhao B. A reactive oxygen species activation mechanism contributes to Sophoridine-induced apoptosis in rat liver BRL-3A cells. *J Ethnopharmacol* 2018 Mar 1;213:376-383.
140. Wu X, Luo P, Rao W, Dai S, Zhang L, Ma W, Pu J, Yu Y, Wang J, Fei Z. Homer1a Attenuates Hydrogen Peroxide-Induced Oxidative Damage in HT-22 Cells through AMPK-Dependent Autophagy. *FRONT NEUROSCI-SWITZ* 2018 Feb 9;12:51.
141. Zhang JJ, Zhang YZ, Peng JJ, Li NS, Xiong XM, Ma QL, Luo XJ, Liu B, Peng J. Atorvastatin exerts inhibitory effect on endothelial senescence in hyperlipidemic rats through a mechanism involving down-regulation of miR-21-5p/203a-3p. *Mech Ageing Dev* 2018 Jan;169:10-18.
142. Liu MX, Jin L, Sun SJ, Liu P, Feng X, Cheng ZL, Liu WR, Guan KL, Shi YH, Yuan HX, Xiong Y. Metabolic reprogramming by PCK1 promotes TCA cataplerosis, oxidative stress and apoptosis in liver cancer cells and suppresses hepatocellular carcinoma. *Oncogene* 2018 Jan 16.
143. Li J, Chen L, Liu Q, Tang M, Wang Y, Yu J. Buformin suppresses proliferation and invasion via AMPK/S6 pathway in cervical cancer and synergizes with paclitaxel. *Cancer Biol Ther* 2018 Feb 5:1-11.
144. Zhang Z, Xu S, Wang Y, Yu Y, Li F, Zhu H, Shen Y, Huang S, Guo S. Near-infrared triggered co-delivery of doxorubicin and quercetin by using gold nanocages with tetradecanol to maximize anti-tumor effects on MCF-7/ADR cells. *J COLLOID INTERF SCI* 2018 Jan 1;509:47-57.
145. Jin Q, Li R, Hu N, Xin T, Zhu P, Hu S, Ma S, Zhu H, Ren J, Zhou H. DUSP1 alleviates cardiac ischemia/reperfusion injury by suppressing the Mff-required mitochondrial fission and Bnip3-related mitophagy via the JNK pathways. *Redox Biol* 2018 Apr;14:576-587.
146. Chen Z, Li Y, Wang Y, Qian J, Ma H, Wang X, Jiang G, Liu M, An Y, Ma L, Kang L, Jia J, Yang C, Zhang G, Chen Y, Gao W, Fu M, Huang Z, Tang H, Zhu Y, Ge J, Gong H, Zou Y. Cardiomyocyte-Restricted Low Density Lipoprotein Receptor-Related Protein 6 (LRP6) Deletion Leads to Lethal Dilated Cardiomyopathy Partly Through Drp1 Signaling. *Theranostics* 2018 Jan 1;8(3):627-643.
147. Li Y, Liu Q, Sun J, Wang J, Liu X, Gao J. Mitochondrial protective mechanism of simvastatin protects against amyloid β peptide-induced injury in SH-SY5Y cells. *Int J Mol Med* 2018 Feb 5.
148. Zheng S, Chang W, Zhang M, Shi H, Lou H. Chiloscyphenol A derived from Chinese liverworts exerts fungicidal action by eliciting both mitochondrial dysfunction and plasma membrane destruction. *SCI REP-UK* 2018 Jan 10;8(1):326.
149. Zhao Q, Ye M, Yang W, Wang M, Li M, Gu C, Zhao L, Zhang Z, Han W, Fan W, Meng Y. Effect of Mst1 on Endometriosis Apoptosis and Migration: Role of Drp1-Related Mitochondrial Fission and Parkin-Required Mitophagy. *CELL PHYSIOL BIOCHEM* 2018;45(3):1172-1190.
150. Bao H, Zhang Q, Du Y, Zhang C, Xu H, Zhu Z, Yan Z. Apoptosis induction in K562 human myelogenous leukaemia cells is connected to the modulation of Wnt/ β -catenin signalling by BHX, a novel pyrazoline derivative. *CELL PROLIFERAT* 2018 Jan 16.
151. Sun Y, Wang C, Meng Q, Liu Z, Huo X, Sun P, Sun H, Ma X, Peng J, Liu K. Targeting P-glycoprotein and SORCIN: Dihydropyridinone strengthens anti-proliferative efficiency of adriamycin via MAPK/ERK and Ca²⁺-mediated apoptosis pathways in MCF-7/ADR and K562/ADR. *J Cell Physiol* 2018 Apr;233(4):3066-3079.
152. Zhou DC, Su YH, Jiang FQ, Xia JB, Wu HY, Chang ZS, Peng WT, Song GH, Park KS, Kim SK, Cai DQ, Zheng L, Qi XF. CpG oligodeoxynucleotide preconditioning improves cardiac function after myocardial infarction via modulation of energy metabolism and angiogenesis. *J Cell Physiol* 2018 May;233(5):4245-4257.
153. Liu L, Luo C, Luo Y, Chen L, Liu Y, Wang Y, Han J, Zhang Y, Wei N, Xie Z, Wu W, Wu G, Feng Y. MRPL33 and its splicing regulator hnRNPK are required for mitochondria function and implicated in tumor progression. *Oncogene* 2018 Jan 4;37(1):86-94.
154. Zeng X, Yang J, Hu O, Huang J, Ran L, Chen M, Zhang Y, Zhou X, Zhu J, Zhang Q, Yi L, Mi M. Dihydropyridinone Ameliorates Nonalcoholic Fatty Liver Disease by Improving Mitochondrial Respiratory Capacity and Redox Homeostasis through Modulation of SIRT3 Signaling. *ANTIOXID REDOX SIGN* 2018 Feb 21.
155. Zhou H, Du W, Li Y, Shi C, Hu N, Ma S, Wang W, Ren J. Effects of melatonin on fatty liver disease: The role of NR4A1/DNA-PKcs/p53 pathway, mitochondrial fission, and mitophagy. *J Pineal Res* 2018 Jan;64(1).
156. Cheng C, Wang T, Song Z, Peng L, Gao M, Hermine O, Rousseaux S, Khochbin S, Mi JQ, Wang J. Induction of autophagy and autophagy-dependent apoptosis in diffuse large B-cell lymphoma by a new antimalarial artemisinin derivative, SM1044. *CANCER MED-US* 2018 Feb;7(2):380-396.
157. Li H, Yang XL, Song HL, Zhang S, Long XZ. Effects of direct current on *Klebsiella* spp. viability and corresponding resistance gene expression in simulative bio-electrochemical reactors. *Chemosphere* 2018 Apr;196:251-259.
158. Wang W, Shi Q, Dou S, Li G, Shi X, Jiang X, Wang Z, Yu D, Chen G, Wang R, Xiao H, Hou C, Feng J, Shen B, Ma Y, Han G. Negative regulation of Nod-like receptor protein 3 inflammasome activation by T cell Ig mucin-3 protects against peritonitis. *Immunology* 2018 Jan;153(1):71-83.
159. Shen C, Li M, Zhang P, Guo Y, Zhang H, Zheng B, Teng H, Zhou T, Guo X, Huo R. A Comparative Proteome Profile of Female Mouse Gonads Suggests a Tight Link between the Electron Transport Chain and Meiosis Initiation. *Mol Cell Proteomics* 2018 Jan;17(1):31-42.
160. Zhang Y, Wei J, Xu J, Leong WS, Liu G, Ji T, Cheng Z, Wang J, Lang J, Zhao Y, You L, Zhao X, Wei T, Anderson GJ, Qi S, Kong J, Nie G, Li S. Suppression of Tumor Energy Supply by Liposomal Nanoparticle-Mediated Inhibition of Aerobic Glycolysis. *ACS APPL MATER INTER* 2018 Jan 24;10(3):2347-2353.
161. Zhu J, Wang Q, Li C, Lu Y, Hu H, Qin B, Xun Y, Zhu Y, Wu Y, Zhang J, Wang S. Inhibiting inflammation and modulating oxidative stress in oxalate-induced nephrolithiasis with the Nrf2 activator dimethyl fumarate. *FREE RADICAL BIO*

- MED 134:9-22. 2019 Apr
162. Nan LP, Wang F, Ran D, Zhou SF, Liu Y, Zhang Z, Huang ZN, Wang ZY, Wang JC, Feng XM, Zhang L Naringin alleviates H₂O₂-induced apoptosis via the PI3K/Akt pathway in rat nucleus pulposus-derived mesenchymal stem cells. *Connect Tissue Res* 2019 Jul 11:1-14
163. Chen Z, Gao Y, Lv B, Sun F, Yao W, Wang Y, Fu X Hypoionic Shock Facilitates Aminoglycoside Killing of Both Nutrient Shift- and Starvation-Induced Bacterial Persister Cells by Rapidly Enhancing Aminoglycoside Uptake. *Front Microbiol* 10:2028. 2019 Sep 6
164. Youbo Zhao, Menghuan Li, Xuemei Yao, Yang Fei, Zhenghong Lin, Zhengguo Li, Kaiyong Cai, Yanli Zhao, Zhong Luo HCAR1/MCT1 Regulates Tumor Ferroptosis through the Lactate-Mediated AMPK-SCD1 Activity and Its Therapeutic Implications *Cell Rep* 2020 Dec 8;33(10):108487.
165. Yanna Zhao, Boyan Lv, Fengqi Sun, Jiafeng Liu, Yan Wang, Yuanyuan Gao, Feng Qi, Zengyi Chang, Xinmiao Fu Rapid Freezing Enables Aminoglycosides To Eradicate Bacterial Persisters via Enhancing Mechanosensitive Channel MscL-Mediated Antibiotic Uptake *mBio* 2020 Feb 11;11(1):e03239-19.
166. Yinlong Yang, Yingying Yue, Nannan Song, Cuiling Li, Zenglin Yuan, Yan Wang, Yue Ma, Hui Li, Fengyu Zhang, Weiwei Wang, Haihong Jia, Peng Li, Xiaobing Li, Qi Wang, Zhe Ding, Hongjie Dong, Lichuan Gu, Bingqing Li The YdiU Domain Modulates Bacterial Stress Signaling through Mn²⁺-Dependent UMPylation *Cell Rep* 2020 Sep 22;32(12):108161.
167. Jie Zhou, Jie Yang, Yu-Meng Wang, Hong Ding, Tu-Shuai Li, Zhi-Hong Liu, Li Chen, Rui-Qing Jiao, Dong-Mei Zhang, Ling-Dong Kong IL-6/STAT3 signaling activation exacerbates high fructose-induced podocyte hypertrophy by ketohexokinase-A-mediated tristetraproline down-regulation *Cell Signal* 2021 Oct;86:110082.
168. Cheng Zeng, Tingting Zou, Junyan Qu, Xu Chen, Suping Zhang, Zhenghong Lin Cycloviobuxine D Induced-Mitophagy through the p65/BNIP3/LC3 Axis Potentiates Its Apoptosis-Inducing Effects in Lung Cancer Cells *Int J Mol Sci* 2021 May 29;22(11):5820.
169. Xiangfei Li, Rumeng Han, Teng Bao, Tolbert Osire, Xian Zhang, Meijuan Xu, Taowei Yang, Zhiming Rao Citrulline deiminase pathway provides ATP and boosts growth of *Clostridium carboxidivorans* P7 *Biotechnol Biofuels* 2021 Oct 16;14(1):204.
170. Meng Duan, Hainan Chen, Linjie Yin, Xiao Zhu, Petr Novák, Yuncheng Lv, Guojun Zhao, Kai Yin Mitochondrial apolipoprotein A-I binding protein alleviates atherosclerosis by regulating mitophagy and macrophage polarization *Cell Commun Signal* 2022 May 7;20(1):60.

Version 2024.03.12