



碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-168-3301或800-8283301  
订货e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

## Triton X-100/曲拉通X-100

产品编号	产品名称	包装
ST795	Triton X-100/曲拉通X-100	100ml

### 产品简介:

- Triton X-100中文名为曲拉通X-100, 分子式为 $t\text{-Oct-C}_6\text{H}_4\text{-(OCH}_2\text{CH}_2)_x\text{OH}$ ,  $x=9-10$ , 平均分子量为647。
- 常用的非离子性去垢剂, 进口分装。

### 包装清单:

产品编号	产品名称	包装
ST795	Triton X-100/曲拉通X-100	100ml
—	说明书	1份

### 保存条件:

室温保存。

### 注意事项:

- 本产品对人体有害, 操作时请小心, 并注意有效防护以避免直接接触人体或吸入体内。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

### 使用本产品的文献:

1. Xin H, Sha X, Jiang X, Chen L, Law K, Gu J, Chen Y, Wang X, Fang X. The brain targeting mechanism of Angiopep-conjugated poly(ethylene glycol)-co-poly( $\epsilon$ -caprolactone) nanoparticles. *Biomaterials*. 2012 Feb;33(5):1673-81.
2. Wang SJ, Liu WJ, Wu CJ, Ma FH, Ahmad S, Liu BR, Han L, Jiang XP, Zhang SJ, Yang LG. Melatonin suppresses apoptosis and stimulates progesterone production by bovine granulosa cells via its receptors (MT1 and MT2). *Theriogenology*. 2012 Oct 15;78(7):1517-26.
3. Zhou X, Tu J, Li Q, Kolosov VP, Perelman JM. Hypoxia induces mucin expression and secretion in human bronchial epithelial cells. *Transl Res*. 2012 Dec;160(6):419-27.
4. Fei X, Jiang S, Zhang S, Li Y, Ge J, He B, Goldstein S, Ruiz G. Isolation, Culture, and Identification of Amniotic Fluid-Derived Mesenchymal Stem Cells. *Cell Biochem Biophys*. 2013 Nov;67(2):689-94.
5. Xiao YD, Huang YY, Wang HX, Wu Y, Leng Y, Liu M, Sun Q, Xia ZY. Thioredoxin-Interacting Protein Mediates NLRP3 Inflammasome Activation Involved in the Susceptibility to Ischemic Acute Kidney Injury in Diabetes. *Oxid Med Cell Longev*. 2016;2016:2386068.
6. Wu Z, Liu H, Ren W, Dai F, Chang J, Li B. VE-cadherin involved in the pulmonary microvascular endothelial cell barrier injury induced by angiotensin II through modulating the cellular apoptosis and skeletal rearrangement. *Am J Transl Res*. 2016 Oct 15;8(10):4310-9.
7. Yuan Z, Li Y, Hu Y, You J, Higashisaka K, Nagano K, Tsutsumi Y, Gao J. Chitosan nanoparticles and their Tween 80 modified counterparts disrupt the developmental profile of zebrafish embryos. *Int J Pharm*. 2016 Dec 30;515(1-2):644-56.
8. Hu J, Li G, Qu L, Li N, Liu W, Xia D, Hongdu B, Lin X, Xu C, Lou Y, He Q, Ma D, Chen Y. TMEM166/EVA1A interacts with ATG16L1 and induces autophagosome formation and cell death. *Cell Death Dis*. 2016 Aug 4;7(8):e2323.

Version 2016.12.08