

右旋糖酐(MW4,000)


产品编号	产品名称	包装
ST3369-1g	右旋糖酐(MW4,000)	1g
ST3369-5g	右旋糖酐(MW4,000)	5g
ST3369-25g	右旋糖酐(MW4,000)	25g

产品简介:

CAS Number	Chemical Formula	Molecular Weight	Purity	Grade
9004-54-0	[C ₆ H ₁₀ O ₅] _n	~4000	-	Reagent grade

➤ 基本信息(General Information):

Name (Chinese)	右旋糖酐(MW4,000)
Name (English)	Dextran (MW4,000)
Specifications	Reagent grade
Chemical Formula	[C ₆ H ₁₀ O ₅] _n
Synonym (Chinese)	葡聚糖; 甘露三糖
Synonym (English)	3',6'-dihydroxy-3-Oxospiro(isobenzofuran-1-(3H), 9'-[9Hxanthen]-5 (or 6)-ylcarbamothiate
Beilstein Registry No.	-
EINECS Number	-
MDL Number	MFCD00130935
UNSPSC Code	-



➤ 产品描述(Description):

Application	advanced drug delivery cell analysis life science and biopharma solubility enhancement
General description	右旋糖酐是一种多糖。葡聚糖液体表现为胶体; 也用作血浆替代品。它可能具有硅化作用, 有助于涂覆原始浆膜表面。

➤ 性质(Properties):

Technique(s)	cell based assay: suitable cryopreservation: suitable
Solubility	~50mg/ml (H ₂ O)

➤ 安全信息(Safety Information):

Hazard Pictogram Codes	-
Signal Word	-
Hazard Statements	-
Precautionary Statements	-
Personal Protective Equipment	Eyeshields, Gloves, type N95 (US)
Hazard Codes (Europe)	-
Risk Codes (Europe)	-
Safety Codes (Europe)	-
RIDADR	-
WGK Germany	2

RTECS	-
Flash Point (F)	-
Flash Point (C)	-

包装清单:

产品编号	产品名称	包装
ST3369-1g	右旋糖酐(MW4,000)	1g
ST3369-5g	右旋糖酐(MW4,000)	5g
ST3369-25g	右旋糖酐(MW4,000)	25g
—	说明书	1份

保存条件:

室温保存, 三年有效。

注意事项:

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

Version 2024.06.27